


Confident Quantitation

Any compound, any matrix, any user.



ThermoFisher
SCIENTIFIC

Tomorrow's Quantitation with Triple Quadrupole MS

The world leader in serving science

Confident Quantitation

1

An Introduction to the TSQ portfolio

2

Features and Benefits

3

Robust Solution

Confident Quantitation

1

An Introduction to the TSQ portfolio

2

Features and Benefits

3

Robust Solution

Environmental and Food Safety
Clinical Research
Pharma QA/QC



TSQ Fortis

- Mass Range m/z 5 – 3000
- Max Resolution **0.4 FWHM**
- Max 30,000 transitions per run
- Polarity Switching < 20 msec
- Dynamic interscan time
- 600 SRM/sec
- TNG software
- Chromeleon support
- **50,000:1 S/N**

Food Safety
Pharma
Clinical Research
Forensic Toxicology



TSQ Quantis

- Mass Range m/z 5 – 3000
- Max Resolution **0.4 FWHM**
- Max 30,000 transitions per run
- Polarity Switching < 20 msec
- Dynamic interscan time
- 600 SRM/sec
- TNG software
- Chromeleon support
- **150,000:1 S/N**

Pharma/Biopharma
Environmental and Food Safety
Omics



TSQ Altis

- Mass Range m/z 5 – 2000
- Max Resolution **0.2 FWHM**
- Max 30,000 transitions per run
- Polarity Switching < 20 msec
- Dynamic interscan time
- 600 SRM/sec
- TNG software
- Chromeleon support
- **500,000:1 S/N**

TSQ Fortis: Affordable Productivity, For Everyone



Active Ion Management Plus (AIM+) - The next step in precision design delivers the ultimate in ion management, inception to detection, from the OptaMax™ ion source housing to the enhanced electron multiplier. Incorporates segmented quadrupoles with hyperbolic surfaces and enhanced RF Electronics to further optimize ion management precision, reliability, speed, and reproducibility.

Enhanced dual-mode electron multiplier detector
ensures excellent linearity and dynamic range

NEW!

Segmented Quadrupoles with hyperbolic surfaces
for enhanced performance with both SRM and H-SRM (0.4 FWHM.)

NEW!

Active collision cell with axial DC field
facilitates more SRMs/sec

Matrix Separator Ion Guide (MSIG)
Ensures Robustness while efficiently transmitting the ion beam

NEW!

Ion beam guide with neutral blocker
Reduces chemical background

NEW!

OptaMax™ NG
APCI ready



TSQ Quantis: Confidence, Day After Day



Active Ion Management Plus (AIM+) - The next step in precision design delivers the ultimate in ion management, inception to detection, from the OptaMax™ ion source housing to the enhanced electron multiplier. Incorporates segmented quadrupoles with hyperbolic surfaces and enhanced RF Electronics to further optimize ion management precision, reliability, speed, and reproducibility.

Enhanced dual-mode electron multiplier detector
ensures excellent linearity and dynamic range

NEW!

Stacked ring ion guide (SRIG)
Increases ion flux

Segmented Quadrupoles
with hyperbolic surfaces for enhanced performance with both SRM and H-SRM (0.4 FWHM)

NEW!

OptaMax™ NG
APCI ready

NEW!

Ion beam guide with neutral blocker
Reduces chemical background

Active collision cell with axial DC field
facilitates more SRMs/sec

TSQ Altis: Sensitivity with Robustness, No Compromises

AIM+
TECHNOLOGY

Active Ion Management Plus (AIM+) - The next step in precision design delivers the ultimate in ion management, inception to detection, from the OptaMax™ ion source housing to the enhanced electron multiplier. Incorporates segmented quadrupoles with hyperbolic surface and enhanced RF Electronics to further optimize ion management precision, reliability, speed, and reproducibility.

Ion beam guide with neutral blocker
Reduces chemical background

High capacity ion transfer tube (HCTT)
Increases ion flux

Electrodynamic ion funnel (EDIF)
Increases ion flux

Segmented Quadrupoles
with hyperbolic surface for enhanced performance with both SRM and H-SRM (0.2 FWHM)

Active collision cell with axial DC field
facilitates more SRMs/sec

OptaMax™ NG
APCI ready

NEW!


NEW!

Enhanced dual-mode electron multiplier detector
Ensures excellent linearity and dynamic range

NEW!

OptaMax NG Source Housing

Benefits: Reliable and consistent performance with improved usability!



• Re-designed APCI discharge assembly

- Built-in to every source (separate APCI sprayer required for APCI mode)
- Re-designed on/off switch (to improve usability)

• Re-designed HESI Sprayer

- Needle adjustment is no longer possible during acquisition (locked position)
- Tool available to help the user to correctly set needle protrusion

• Usability and Consistency

- Vertical adjustment moved to the side for easier access
- New drain insert with improved latching and locating pin to prevent rotation
- Improved sprayer alignment and stability
- New finer threads on HESI and APCI sprayers to make installation easier

Segmented Quadrupoles

Benefits: Increased Sensitivity (more significant at higher mass range)
Flat tuning for consistent and robust performance


- The use of RF only pre-filters (segments) between the entrance lens and the quadrupole minimizes the effects of fringe fields, leading to improved transmission (and therefore sensitivity) at unit and higher resolution.
- With the RF only pre-filter, the tuning of several lenses is flat across mass range allowing the voltage to be set and not tuned. This helps reducing the complexity of the tune and making the systems more consistent.



Detector

Benefits: Increased electron multiplier lifetime. Increased Uptime!


- Increased number of dynodes (21) for extended lifetime.
- Improved electron multiplier calibration routine.
- Excellent linearity and dynamic range across the mass range.
- Reduced number of service visits leading to more uptime.



RF Circuitry

Benefits: More compounds in the same run or longer dwells on existing method

- New main RF/DC electronics
- Analyze more compounds in the same time window or better Quantitation results with better ion statistics (more scans across your chromatographic peak)
- Up to **600 SRM/sec**



Confident Quantitation

1

Introduction to TSQ Fortis - Critical Features


2

Features and Benefits

3

Robust Solution – Content

Features that Enable Every Analytical Laboratory



What Makes the new Triple Quads Robust?



Robustness

The ability to perform at the expected level (LOD, LOQ, MRL) under adverse conditions (complex samples, limited sample preparation) for the desired period of time without maintenance.




OptaMax NG	ITT/Sweep Cone	MSIG	Ion Beam Guide
			
Source housing	Ion Transfer Tube and Sweep Cone	Matrix Separator Ion Guide	Ion beam guide with neutral blocker
			
<i>Reliable and consistent performance with improved usability!</i>	<i>No need to break vacuum for basic maintenance!</i>	<i>Eliminates neutrals while ensuring efficient transmission and sensitivity</i>	<i>Blocks neutrals - keeps the ion path clean and reduces chemical background!</i>

What Makes the new Triple Quads Consistent?



Consistency

The ability to report equivalent results every time independently of the system or the user (day-to-day and system-to-system)


OptaMax NG	Optics	Segmented Quads
		
Re-designed Source Housing	Improved Optics Alignment	Segmented Quads with hyperbolic surfaces
<i>Reliable and consistent performance with improved usability!</i>	<i>More consistent performance from system-to-system!</i>	<i>Flat tuning across the mass range!</i>

What Makes the new Triple Quads Reliable?



Reliability

Increased time between service calls
(mean time between failure- MTBF)
Increased uptime




Factors Contributing to Enhanced Speed



Speed


Critical feature when targeting more analytes in the same method or when increasing throughput by reducing the analysis time.

- **Active Collision Cell with axial DC field**
 - 90° cell design for noise reduction
- **New main RF/DC electronics**
 - Analyze more compounds in the same time window or better Quantitation results with better ion statistics (more scans across your chromatographic peak)
 - Up to **600 SRM/sec**





Pesticides in Leek: Robust, Reliable, Fast Quantitation Workflows

**Increase Productivity - Monitoring more compounds per unit time!
Excellent Quantitative Performance at Lower Dwell Times!
Robust data enables Maximizing Uptime!**



SRM Visualization from Instrument Control Software 3.0 displaying the number of transitions per unit time



**~ 160 Transitions Monitored Simultaneously with
Polarity Switching. Excellent Reproducibility (% RSD 2.3) below the MRL**

Sensitivity



Sensitivity

The ability to accurately and precisely detect and/or quantify an analyte

Different criteria can be used to establish a limit of quantitation.

- **Multiple Product Ions** – typically one quantifier ion and multiple confirming ions
- **Excellent precision** of replicate measurements at reporting level
- **Excellent accuracy** at reporting level
- **Extended Linearity**

Active Ion Management Plus (AIM+)

The next step in precision design delivers the ultimate in ion management, inception to detection, from the OptaMax NG source housing to the enhanced electron multiplier. Incorporates segmented quadrupoles and enhanced RF Electronics to further optimize ion management precision, reliability, speed, and reproducibility.




TSQ Altis: Demonstration of sensitivity for lipidomics and metabolomics


LPA 17:0 @ 0.01 ng/mL

TSQ Quantiva

TSQ Altis




RSD 3.4% (n=3)



RSD 0.9% (n=3)


Glutamic Acid



LPA at 0.01 ng/mL with RSD below 1%! 2 x greater response on TSQ Altis! Confirming Ion detected on TSQ Altis!
A 2.5 fold response increase was observed when transferring the metabolomics assay from TSQ Quantiva to TSQ Altis.

TSQ Altis: Sensitivity, Precision and Dynamic Range


Superior Sensitivity and
Reproducibility - Fluticasone
Propionate in matrix!
25 fg on column with 4.2% CV!



UHPLC : Thermo Scientific™ Vanquish™ Flex Binary UHPL System
Flow rate: 0.75 mL/min
Solvent : 75% MeOH:25% H_2O with 0.01% NH_4OH
MS: TSQ Altis Triple Quadrupole MS
Software: TraceFinder Software 4.1

Experimental Details

Superior Sensitivity and
Reproducibility – Testosterone in
human plasma
5 fg on column with 3% CV!



UHPLC : Thermo Scientific™ Vanquish™ Flex Binary UPLC System
Flow rate: 0.4 mL/min
Solvent A: 0.5 mM NH_4F in Water
Solvent B: MeOH
MS: TSQ Altis Triple Quadrupole MS
Software: TraceFinder Software 4.1


Experimental Details

For Research Use Only

TSQ Altis: Quantitation of Therapeutic Drugs in Plasma

Compound Name	TSQ Quantiva LOQ (pg/mL)	TSQ Altis LOQ (pg/mL)
Desomorphine	5	5
Desmethyldoxepin	10	2.5
Flecainide	2.5	1
Midazolam	5	2.5
Imipramine	10	2.5
Amitriptyline	10	2.5
Fluoxetine	5	5
Diazepam	5	2.5


Quantitation of Desmethyldoxepin in plasma
2.5 pg/mL



~ 2.5 X average sensitivity improvement over TSQ Quantiva

Application Note 64977

TSQ Quantis: Sensitivity Improvement for Pesticide Residue Analysis




3 x average sensitivity improvement across mass range!

Application Note 64971

TSQ Quantis: Sensitivity in Regulatory Methodology (EPA 544)

TSQ Endura



In 2014, half-a-million people were without drinking water in the state of Ohio due to high toxin level!





Algae blooms in Western Lake Erie

EPA 544 released in February 2015



TSQ Quantis




~ 4 x average sensitivity improvement over TSQ Endura

Application Note 64968

How to demonstrate sensitivity? Which system is right for you?


Sensitivity is a lot more than peak area! It should be demonstrated as a precise and accurate measurement!

TSQ Fortis




RSD = 3% (n=5)

TSQ Endura



RSD = 4% (n=5)

TSQ Quantum Access Max




RSD = 24% (n=5)

What to expect when evaluating Limits of Quantitation based on precision and accuracy?

TSQ Quantis > 3x TSQ Endura (2x) TSQ Fortis > 3-5x TSQ Quantum Access Max


How to demonstrate sensitivity? Which system is right for you?

TSQ Quantis



RSD = 5% (n=5)

TSQ Quantiva



RSD = 5% (n=5)

What to expect when evaluating Limits of Quantitation based on precision and accuracy?

TSQ Altis > 2-3x TSQ Quantiva (2x) TSQ Quantis

Specificity



Resolution

The ability to filter the m/z of interest and obtain a better signal to noise ratio in complex matrices.

*New segmented quadrupoles
with hyperbolic surfaces*

TSQ Altis – 0.2 Da FWHM


TSQ Quantis – 0.4 Da FWHM

TSQ Fortis – 0.4 Da FWHM



Superior Sensitivity with H-SRM (0.2 Da FWHM) – GPSVFPLAPSSK

GPSVFPLAPSSK - Peptide from monoclonal antibody hinge region




25 ng GPSVFPLAPSSK and IS



Experimental Details

LC : Ultimate NCS-3500RS
Flow rate: 25 μ L/min
Solvent A: 2% ACN in H₂O w/0.1% FA
Solvent B: 5% H₂O in ACN w/0.1% FA



Ease of Use: Simple Maintenance

Simple user maintenance for additional uptime



5 min operation!
No need to break vacuum!

OptaMax NG



Re-designed
source housing

NEW!

*New tool available to ensure
consistent positioning of
HESI needle*

ITT/Sweep Cone



Ion transfer tube
and sweep cone

*No need to break vacuum for
basic maintenance
5 minute Operation!*

Ease of Use: Method Setup made easy with Method Templates

The screenshot displays the Method Editor software interface, showing various panels for method setup. The main panel displays a table of scan data, a heatmap of dwell time per trajectory, and a line graph showing the number of transitions per cycle. The interface is organized into sections: Global Parameters, Scan Parameters, and Summary. The Scan Parameters section includes a table with columns for Compound, Retention Time (min), RT Window (min), Polarity, Precursor (m/z), Product (m/z), and Collision Energy (V). The table lists various compounds such as Chlorzoxazone, Pregabalin, and others. The Summary section includes MS/MS Properties and other method settings.

Compound	Retention Time (min)	RT Window (min)	Polarity	Precursor (m/z)	Product (m/z)	Collision Energy (V)
016 Chlorzoxazone	9.676	0.5	Positive	445.18	281.034	18.24
017 Pregabalin	9.895	0.5	Positive	421.207	213.111	25.42
018 Pregabalin	9.895	0.5	Positive	421.207	231.111	24.76
019 Pregabalin	9.895	0.5	Positive	421.207	384.183	14.9
022 Abacavir Metabolite	9.952	0.5	Positive	399.526	395.111	24.31
021 Abacavir Metabolite	9.952	0.5	Positive	399.526	307.189	19.58
024 Abacavir Metabolite	9.952	0.5	Positive	399.526	367.262	12.88
023 Acemetacin	10.024	0.5	Positive	339.195	138	15.78
026 Acemetacin	10.024	0.5	Positive	339.195	143.054	24.76
027 Acemetacin	10.024	0.5	Positive	339.195	171.125	14.9
028 Acemetacin	10.024	0.5	Positive	339.195	289.111	14.25
029 Acemetacin H4	10.026	0.5	Positive	339.222	138.01	41
030 Acemetacin H4	10.026	0.5	Positive	339.222	143.06	26
032 Acemetacin H4	10.026	0.5	Positive	339.222	171.07	16
031 Acemetacin H4	10.026	0.5	Positive	339.222	179.014	23.75
033 Acemetacin	10.03	0.5	Positive	329.89	294.111	34.22
034 Acemetacin	10.03	0.5	Positive	329.89	319	21.88
035 Dabofosin H4	10.245	0.5	Positive	394.238	106.982	39.83
036 Dabofosin H4	10.245	0.5	Positive	394.238	177.04	14.25
037 Dabofosin H4	10.245	0.5	Positive	394.238	339.195	16.25
038 Penazapam	10.294	0.5	Positive	307.18	96.899	23.2
039 Penazapam	10.294	0.5	Positive	307.18	147.034	20.47
040 Penazapam	10.294	0.5	Positive	307.18	181.183	17.79

Recommended conditions available to the user:
Optimized conditions specific for the MS model used

Instrument Method Templates available in the latest ICSW – version 3.1

- **Clinical Research**
 - Vitamin D
 - Steroids
 - Immunosuppressants
- **Forensic Toxicology**
 - ETG/ETS
- **Food Safety**
 - Pesticides
 - Veterinary Medicines
- **Environmental**
 - Perchlorate
 - HAAs
 - EPA 544
 - EPA 545
 - PFCs
- **Pharma QA/QC**
 - Phenylephrine



TSQ Fortis

- **Clinical Research**
 - Steroids
 - Immunosuppressants
- **Forensic Toxicology**
 - Drugs of Abuse
 - ETG/ETS
- **Food Safety**
 - Pesticides
 - Veterinary Medicines
- **Environmental**
 - HAAs
 - EPA 544
 - EPA 545
- **Pharma**
 - Phenylephrine
 - Pharma Panel
- **Peptide Analysis**
 - PRTC Heavy
 - PRTC Hela




TSQ Quantis

- **Clinical Research**
 - Steroids
 - Immunosuppressants
- **Food Safety**
 - Pesticides
 - Veterinary Medicines
- **Environmental**
 - PFCs
- **Pharma**
 - Fluticasone
 - Pharma Panel
- **Peptide Analysis**
 - PRTC Heavy
 - PRTC Hela




TSQ Altis

Ease of Use: Conversion Tool



Automated Conversion Tool – conversion of source parameters from our competitors into TSQ parameters

Ease of Use: Compound Optimization



Automated Compound Optimization Tool – Optimization of SRM parameters via infusion or Autosampler!

Flexibility

- **LC Options**

- **High Sample Throughput**

- Multi channel – Transcend and Prelude
- CTC PAL 3 – reduced time between injections, barcode reading

- **High Performance LC**

- UHPLC – Vanquish Horizon and Vanquish Flex

- **MS Ion Sources**

- H-ESI
- APCI
- OptaMax – APCI ready (only need to change the sprayer – 2 min operation)
- OptaMax Duet – Automated HESI/APCI methods (no need to change sprayer)
- APPI
- Easy spray source and flex for nanoLC applications



- **Software Options**

- TraceFinder Software

- Method development to automated report generation – offering comprehensive solution for every target application

- Chromeleon

- Data acquisition and processing

- Skyline

- Compatibility with TraceFinder



Confident Quantitation

1 *Introduction to TSQ Fortis - Critical Features*

2 *Features and Benefits*

3 **Robust Solution**

Markets and Opportunities



Clinical Research

In Clinical Research (RUO), many analytical methods are used. But, MS is gaining popularity owing to speed, selectivity, and sensitivity



Environmental Safety

MS is widely used as a detection tool. But, presence of polar compounds either require IC separation, or a time consuming derivatisation step.



Food Safety

LC-MS is a popular choice to address an ever increasing list of analytes across multiple types of complex matrices. Selectivity, sensitivity, robustness play critical role



Pharma/ Biopharma

From high throughput screening to QA/QC based quantitation, from targeted quan in discovery Bioanalysis to DMPK – robust, reliable MS based detection is critical for confident quantitation



Targeted Proteomics

Ability to develop sensitive and selective quantitation assays for multiple peptides from several protein digests ensure LC-MS/MS as a platform of choice

Environmental Safety – Critical Challenges



Increased Confidence for Environmental Safety

Ultimate robustness with outstanding data quality, comprehensive and reliable workflows that ensures every laboratory achieve their business and scientific goals with ease

- **Be Prepared to Address Challenges, Everyday**
 - Robust, reproducible, reliable workflows for quantitation of hundreds of contaminants in a variety of matrices
 - High quality results enabling every environmental safety laboratory
- **Address Productivity Goals with Workflow Solutions**
 - **Speed** – analyze more compounds/injection or reduce overall run times
 - **Robustness** – less maintenance, increased uptime
 - **Ease of use** – easy method setup, data review and customized reporting
 - **Sensitivity** – detect contaminants with class leading sensitivity

Demonstration of Sensitivity - Environmental

Detection of haloacetic acids according to the EPA guidelines

Ion Chromatography

Workhorse for every analytical lab analyzing polar molecules
Spans high-pressure and address cost/sample challenges

TSQ Fortis MS

Robustness offering unmatched productivity with best-in-class sensitivity
Addresses critical quantitation challenges for every analytical laboratory


TraceFinder Software 4.1

One software for multiple LC-MS platform



Application Note 65196

Compound	MDL (µg/L)
MCAA	0.03
MBAA	0.03
DCAA	0.02
DBAA	0.02
BCAA	0.03
TCAA	0.06
BDCAA	0.05
DBCAA	0.15
TBAA	0.15
Dalapon	0.03
Bromate	0.02



Perchlorate in Water



Best-in-class sensitivity with exceptional robustness

Ion Chromatography

Workhorse for every analytical lab analyzing polar molecules
Spans high-pressure and address cost/sample challenges

TSQ Fortis MS

Robustness offering unmatched productivity with best-in-class sensitivity
Addresses critical quantitation challenges for every analytical laboratory.

TraceFinder Software 4.1

One software for multiple LC-MS platform
One software from method development to report generation



Authors

Reza Wajnski, Gaudy Melny, Jonathan Beck, ThermoFisher Scientific, San Jose, CA
Deborah Stathakopoulou, ThermoFisher Scientific, Houston, TX
Heidi Soboni, Lee Park, Archie Huard, Santa Clara Valley Water District, Santa Clara, CA

Keywords

Perchlorate, LC-MS/MS, Ions, perchlorate MS, TSQ Fortis MS, quantitative analysis solution

Goal

Development and easy implementation of a robust, reliable, and reproducible workflow solution for the analysis and quantitation of perchlorate in water using a triple-stage quadrupole mass spectrometer (MS).

Application benefits

- Development of a robust workflow for the analysis and quantitation of perchlorate in water with ion chromatography (IC) and the Thermo Scientific™ TSQ Fortis™ triple-stage quadrupole mass spectrometer (QqQ)
- Leveraging enhanced performance of a robust QqQ with best-in-class sensitivity to reduce cost per sample for every analytical environment

Introduction

Clean drinking water in today's world is becoming more scarce, and any contamination can result in long-lasting damage to human health. Besides purifying water by means of mechanical processes, disinfection also plays an essential role in ensuring the supply of clean drinking water. Drinking water goes through an extensive disinfection process to ensure high quality. However, contaminants continue to cause some severe health risks. In 1997, low-level perchlorate contamination (500 ng/mL or parts per billion) was discovered in the western U.S. Since that time, it has been found in



ThermoFisher
SCIENTIFIC

Application Note 65201



TSQ Fortis (& Ion Chromatography): Environmental Analysis

Haloacetic Acids by IC-MS/MS

Compound	MDL (µg/L)
MCAA	0.03
MBAA	0.03
DCAA	0.02
DBAA	0.02
BCAA	0.03
TCAA	0.06
BDCAA	0.05
DBCAA	0.15
TBAA	0.15
Dalapon	0.03
Bromate	0.02




Perchlorate by IC-MS/MS




Application Note 65196

Application Note 65201

TSQ Quantis: Separation of Perfluorinated Alkyl Substances, 50 pg/mL



PFASs C4 – C10 at 50 pg/mL (0.15 pg on-column); LODs are ~10 pg/mL (0.03 pg O.C.)



C8-PFAS at 50 pg/mL (0.15 pg on-column); LODs are ~10 pg/mL (0.03 pg O.C.)

Vanquish Flex Binary UHPLC system

Delay Column: 3.0 x 50 mm, 5 μ m BDS Hypersil C8
Analytical Column: 2.1x100 mm, 2.6 μ m Accucore C18
Column Temp: 30 C
Mobile Phase: [A] H₂O + 10 mM Am. Acetate; [B] MeOH
Injection Volume: 3 μ L


PFAS can cause reproductive and developmental, liver and kidney, and immunological issues

Increased cholesterol levels among exposed populations, low infant birth weights, effects on the immune system, cancer, and thyroid hormone disruption

Robust, reproducible, sensitive quantitative assays for detection and quantitation of PFAS were developed using Vanquish Flex binary UHPLC system and TSQ Quantis MS

TSQ Quantis MS

Ionization Mode: HESI, Negative ion mode
MS Acquisition Mode: Selective Reaction Monitoring (SRM) **Cycle time: 0.45 s**
Quad Isolation (Q1,Q3) = Unit (0.7 Da FWHM)



TSQ Altis: Confident Quantitation of challenging analytes in environmental matrices


IC-MS/MS solutions for environmental analysis Quantitation of **Glyphosate** at 5 ng/L, 3% RSD

Experimental Details

ICS 5000+
Flow rate: 0.3 mL/min
Eluent Source: Eluent Generator
Mobile Phase: KOH

MS: TSQ Altis

Software: TraceFinder Software 4.1






Confident Quantitation for Food Safety

Ultimate robustness with selectivity, reproducibility enabling comprehensive workflows to ensure every food safety laboratory

- **Be Prepared to Address Challenges, Everyday**
 - Robust, reproducible, reliable workflows for quantitation of hundreds of contaminants in a variety of matrices
 - Address upcoming analytical demands of increased number of contaminants with methods having shorter run times
- **Address Productivity Goals with Workflow Solutions**
 - **Speed** – analyze more compounds/injection or reduce overall run times
 - **Robustness** – less maintenance, increased uptime
 - **Ease of use** – easy method setup, data review and customized reporting
 - **Sensitivity** – detect contaminants with class leading sensitivity

Pesticide Quantitation with Fortis



Sample Prep:



- QuEChERS based

LC Conditions:

- Thermo Scientific™ Accucore™ aQ column (2.1x100x2.6 m)
- 2 μ L injection
- Run Time 15 min


MS Conditions:

- TSQ Fortis MS (comprehensive database with all optimized SRMs)
- Positive/negative switching




TSQ Platform: Robust, Reliable, Fast Quantitation Workflows

Excellent Quantitative Performance at Lower Dwell Times!



SRM Visualization from Instrument Control Software 3.0 displaying the number of transitions per unit time




~ 160 Transitions Monitored Simultaneously with
Polarity Switching. Excellent Reproducibility (% RSD 2.3) below the MRL

Application Note 64971

TSQ Quantis: Demonstration of Robustness – Food Safety




Atrazine QC monitored in leek for more than 400 injections with 4.5% RSD . Red lines represent $\pm 20\%$ response at $10 \mu\text{g}/\text{Kg}$. Yellow lines show the time the system was placed in standby mode for 12h to demonstrate consistent performance after standby period




Application Note 64971

TSQ Quantis: Fipronil Quantitation- Matrix Matched Calibration for UHPLC-MS/MS

Fipronil



Fipronil sulfone



Compound name	Recovery (%) 0.5 µg/kg spike level	Recovery (%) 1 µg/kg spike level	Recovery (%) 5 µg/kg spike level
Fipronil	104	89	99
Fipronil sulfone	99	95	102

Compound name	LOD [µg/kg]	LOQ [µg/kg]	Repeatability (%) 0.5 µg/kg spike level	Repeatability (%) 5 µg/kg spike level
Fipronil	0.1	0.5	8.5	6.1
Fipronil sulfone	0.1	0.5	7.7	6.4

Ion Ratio Calculations	Fipronil (SRM 249.96)	Fipronil (SRM 398.845)	Fipronil sulfone (SRM 282.00)	Fipronil sulfone (SRM 243.845)
Ion ratio* (%) Standard 100 ng/ml	24.4	13.1	78.1	23.6
Ion ratio* (%) Spike 0.5 µg/kg	26.0	14.8	83.4	29.4
Ion ratio repeatability (RSD %)* 0.5 µg/kg spike level	9.1	16.3	4.4	18.9
Ion ratio repeatability* (RSD %) 5 µg/kg spike level	4.4	4.5	1.8	5.8

Reproducibility and Long-term Stability Test

thermoscientific

APPLICATION BROTHER 72453

Rapid analysis of fipronil and fipronil sulfone in eggs by liquid chromatography and triple quadrupole mass spectrometry

Authors
Sascha von Soos, Steffen Mühling
and Michael Gocula
Thermo Fisher Scientific, Special
Solutions Center, Dreieich, Germany

Keywords
Fipronil, fipronil sulfone, eggs
LC-MS, Accucore HPLC, TSO Quantin,
UltiMate 3000 RS-C

Goal
Develop a quick and simple method
for the determination of fipronil
and fipronil sulfone in eggs using
an in-house modified QuEChERS
acetonitrile extraction protocol and
LC-MS/MS determination.


Introduction
Recently, it was reported that nations of eggs contaminated with the
pesticide fipronil have been distributed to more than 17 countries. On
July 20th 2017, it was made public that in some cases the pesticide fipronil
was mixed with another formulation and sprayed on chickens against ticks,
fleas and lice. As the determined levels were in some cases substantially
higher (up to 1.2 mg/kg) than the EU MRL of 0.006 mg/kg for the sum
of fipronil and fipronil sulfone,¹⁾ there is a demand for quick and efficient
methods for the determination of both substances in egg matrix and
especially in chicken meat.

The first purpose is a quick and simple method for the determination of
fipronil and fipronil sulfone in eggs using an in-house modified QuEChERS
acetonitrile extraction protocol.


Experimental
Sample preparation
Egg samples purchased in a local store were extracted using the procedure
described in Figure 1.

- Sample preparation consumables**
- 50 mL conical sterile polypropylene centrifuge tubes, PVN 352652
 - 15 mL conical sterile polypropylene centrifuge tubes, PVN 339690
 - Thermo Scientific™ HyperSep™ dispersive SPE M444 pouch 3000 mg
magnesium sulfate and 1000 mg NaCl, 50 µm, PVN 9010E-340
 - Magnesium sulfate, 99%, for analysis anhydrous, 1210672

ThermoFisher
SCIENTIFIC



Vet Drugs Analysis and Quantitation



- Multi-class veterinary drugs analysis showing
 - Fit for purpose Thermo Scientific™ Accucore™ RP-MS reversed-phase column for robust analysis, great peak shape for wide range of compound classes
 - Generic QuEChERS extraction applied to bovine, salmon file, and milk; easy-to-use, low cost, with no extraction concentration
 - Robust, sensitive, reproducible results for absolute recovery, precision, and low MDLs for most analytes

Sample Prep:


- QuEChERS based approach

LC Conditions:

- Thermo Scientific™ Accucore™ RP-MS (2.1x100x2.6 um)
- 5 µL injection


MS Conditions:

- TSQ Fortis MS (comprehensive database with all optimized SRMs)
- Positive/negative switching



Triclabendazole Sulfoxide
@ 3.4 µg/Kg in beef extract

RSD = 9% (n=3)






Confident Quantitation for Clinical Research and Forensic Toxicology









Address high throughput capabilities with robust, reliable data addressing sensitivity requirements to enable every laboratory achieve their business and scientific goals

- **Be Prepared to Address Challenges, Everyday**
 - Robust, reproducible, reliable methods for fast quantification critical analytes in biological matrices
- **Address Productivity Goals with Workflow Solutions**
 - **Speed** – analyze more compounds/injection or reduce overall run times
 - **Robustness** – less maintenance, increased uptime
 - **Ease of use** – easy method setup, data review and customized reporting
 - **Sensitivity** – detect contaminants with class leading sensitivity

TSQ Fortis: Immunosuppressants in Whole Blood



~2000+ injections showing %CV 6%



For Research Use Only

ThermoFisher
SCIENTIFIC

TSQ Fortis: Demonstration of Robustness – Clinical Research

Immunosuppressants





+2000 injections over 6 days
(CV% 6%)

Red lines represent $\pm 20\%$ of
calculated amounts (ng/mL)




Yellow line represent the period in
which the ion transfer tube was
cleaned (user basic maintenance
~5 min operation) to demonstrate
consistent performance before and
after user maintenance.






Technical Note 65206




TSQ Fortis: Testosterone in Plasma



TSQ Fortis: Clinical Research – Testosterone in plasma




5 pg/mL in plasma (RSD = 2.3%)



Technical Brief 65207

Example SRM Chromatograms – EtS & EtG



EtS & EtG show measurable amounts in Blank Urine samples

- TSQ Quantis MS provides excellent quantitative performance for LC-MS/MS of EtS, EtG and Barbiturates in diluted urine in under 4 minutes.
 - LLOQs of 15 ng/mL and 75 ng/mL were achieved for EtS and EtG, respectively, in diluted urine.
 - LLOQs of 200 ng/mL (urine, diluted 1:10) were achieved for 5 barbiturates
- All target compounds can be analyzed with the same LC method.



Thermo Scientific™ Vanquish™ Horizon UHPLC system


Column: 3.0 x 100 mm, 3 um Synchronis AQ
 Column Temp: 30 C
 Mobile Phase: [A] H₂O + 0.1% HAc; [B] MeOH
 Flow Rate: 0.5 mL/min
 Gradient: see table
 Injection Volume: 10 uL (EtG/EtS); 5 uL (Barbs)

TSQ Quantis Triple Quadrupole MS

Ionization Mode: HESI, Negative ion mode
 MS Acquisition Mode: Selective Reaction Monitoring (SRM) – see SRM Table next slide
 Cycle time: 0.30 s
 Quad Isolation: Q1 = Unit (0.7 Da FWHM); Q3 = Unit or Wide (1.2 Da; Barbs)


TSQ Quantis: Robust Quantitation of ETG and ETS

Excellent Precision for the most challenging quantitative assays




50 pg on column for Ethyl Sulfate in Urine

Technical Note 64970




EtS -d5, %RSD 1.5




EtG -d5, %RSD 2.9

Experimental Details

UHPLC : Vanquish Flex Binary
Flow rate: 0.4mL/min
Solvent A: 0.1%FA in H₂O
Solvent B: 0.1%FA in MeOH



~100 Drugs of Abuse: Diluted Urine @ 1 mL/min




SRM chromatograms of ~100 drugs of abuse in under 1.3 minutes [THC-COOH elutes at 1.21 min, inset]

Vanquish Horizon UHPLC System

Column: 2.1 x 50 mm, 1.9 μ m Hypersil Gold AQ
 Column Temp: 40 C
 Mobile Phase: [A] H₂O + 0.1% HCOOH; [B] ACN + 0.1% HCOOH
Flow Rate: 1.0 mL/min (no split)
 Gradient: see table
 Injection Volume: 2 μ L

Sensitivity, acquisition speed and ruggedness of the TSQ Quantis MS, and the reproducibility of the Vanquish Horizon UHPLC system, make the measurement of ~100 drugs of abuse in diluted urine in under 1.5 minutes possible.



Opiate isomers at m/z 286 are well separated in under 12 s [typical LC peak = 1.2 s wide]

TSQ Quantis MS

Ionization Mode: HESI, Positive ion mode
 MS Acquisition Mode: Selective Reaction Monitoring (SRM) **Cycle time: 0.15 s**
 Quad Isolation (Q1,Q3) = Unit (0.7 Da FWHM)



Tomorrow's quantitation with LC-MS/MS: fast screening and quantitation of drugs of abuse in urine for forensic toxicology

Key Features:

- Multi-targeted analysis (up to 1000) in a single run
- Easy to implement and maintain address high productivity throughout the LC-MS/MS lifecycle

Introduction:

The TSQ Quantis MS is a powerful, flexible, and easy-to-use LC-MS/MS system for forensic toxicology. It is designed for high-throughput screening and quantitation of drugs of abuse in urine. The system features a robust design, easy-to-use interface, and advanced data analysis software. It is ideal for laboratories looking to improve their workflow and reduce the time to results.

Enabling Clinical Research Applications



Quantification of tricyclic antidepressants in human plasma or serum by liquid chromatography-tandem mass spectrometry for clinical research

Authors: Claudio De Nard, Thermo Fisher Scientific GmbH, Germany; Sergio Rodriguez, Thermo Fisher Scientific, Ltd (UK), France

Application benefits:

- Simple, efficient sample preparation for protein precipitation
- 15 analytes in a single quantitative method

Introduction: An analytical method for clinical research for the quantification of 15 tricyclic antidepressants in human plasma or serum by liquid chromatography-tandem mass spectrometry. Analytes: Amitriptyline, Doxepin, Imipramine, Nortriptyline, Propriety, Imipramine, Propriety, Propriety, Propriety, Propriety, Propriety, Propriety, Propriety, Propriety, Propriety, Propriety, Propriety.

Keywords: Tricyclic antidepressants, liquid chromatography, tandem mass spectrometry, clinical research, quantitative analysis

Goal: Implementation of an analytical method for the quantification of 15 tricyclic antidepressants in human plasma or serum using a Thermo Scientific™ TSO Quantum™ triple quadrupole tandem mass spectrometer.



Quantification of benzodiazepines in human plasma or serum by liquid chromatography-tandem mass spectrometry for clinical research

Authors: Claudio De Nard, Thermo Fisher Scientific GmbH, Germany; Sergio Rodriguez, Thermo Fisher Scientific, Ltd (UK), France

Application benefits:

- Simple, efficient sample preparation for protein precipitation
- 15 analytes in a single quantitative method

Introduction: An analytical method for clinical research for the quantification of 15 benzodiazepines in human plasma or serum by liquid chromatography-tandem mass spectrometry. Analytes: Alprazolam, Clonazepam, Clobazam, Clonazepam, Clonazepam, Clonazepam, Clonazepam, Clonazepam, Clonazepam, Clonazepam, Clonazepam, Clonazepam, Clonazepam, Clonazepam, Clonazepam, Clonazepam.

Keywords: Benzodiazepines, liquid chromatography, tandem mass spectrometry, clinical research, quantitative analysis

Goal: Implementation of an analytical method for the quantification of 15 benzodiazepines in human plasma or serum using a Thermo Scientific™ TSO Quantum™ triple quadrupole tandem mass spectrometer.



Quantification of antiepileptics in human plasma or serum by liquid chromatography-tandem mass spectrometry for clinical research

Authors: Claudio De Nard, Thermo Fisher Scientific GmbH, Germany; Sergio Rodriguez, Thermo Fisher Scientific, Ltd (UK), France

Application benefits:

- Simple, efficient sample preparation for protein precipitation
- 15 analytes in a single quantitative method

Introduction: An analytical method for clinical research for the quantification of 15 antiepileptics in human plasma or serum by liquid chromatography-tandem mass spectrometry. Analytes: Carbamazepine, Phenytoin, Valproic acid, Phenytoin, Phenytoin, Phenytoin, Phenytoin, Phenytoin, Phenytoin, Phenytoin, Phenytoin, Phenytoin, Phenytoin, Phenytoin, Phenytoin, Phenytoin.

Keywords: Antiepileptics, liquid chromatography, tandem mass spectrometry, clinical research, quantitative analysis

Goal: Implementation of an analytical method for the quantification of 15 antiepileptics in human plasma or serum using a Thermo Scientific™ TSO Quantum™ triple quadrupole tandem mass spectrometer.



Quantification of antimycotics in human plasma or serum by liquid chromatography-tandem mass spectrometry for clinical research

Authors: Claudio De Nard, Thermo Fisher Scientific GmbH, Germany; Sergio Rodriguez, Thermo Fisher Scientific, Ltd (UK), France

Application benefits:

- Simple, efficient sample preparation for protein precipitation
- 15 analytes in a single quantitative method

Introduction: An analytical method for clinical research for the quantification of 15 antimycotics in human plasma or serum by liquid chromatography-tandem mass spectrometry. Analytes: Itraconazole, Voriconazole, Isavuconazole, Isavuconazole, Isavuconazole, Isavuconazole, Isavuconazole, Isavuconazole, Isavuconazole, Isavuconazole, Isavuconazole, Isavuconazole, Isavuconazole, Isavuconazole, Isavuconazole, Isavuconazole.

Keywords: Antimycotics, liquid chromatography, tandem mass spectrometry, clinical research, quantitative analysis

Goal: Implementation of an analytical method for the quantification of 15 antimycotics in human plasma or serum using a Thermo Scientific™ TSO Quantum™ triple quadrupole tandem mass spectrometer.



Quantification of antidepressants in human plasma or serum by liquid chromatography-tandem mass spectrometry for clinical research

Authors: Claudio De Nard, Thermo Fisher Scientific GmbH, Germany; Sergio Rodriguez, Thermo Fisher Scientific, Ltd (UK), France


Application benefits:

- Simple, efficient sample preparation for protein precipitation
- 15 analytes in a single quantitative method

Introduction: An analytical method for clinical research for the quantification of 15 antidepressants in human plasma or serum by liquid chromatography-tandem mass spectrometry. Analytes: Amitriptyline, Doxepin, Imipramine, Nortriptyline, Propriety, Imipramine, Propriety, Propriety, Propriety, Propriety, Propriety, Propriety, Propriety, Propriety, Propriety, Propriety.

Keywords: Antidepressants, liquid chromatography, tandem mass spectrometry, clinical research, quantitative analysis

Goal: Implementation of an analytical method for the quantification of 15 antidepressants in human plasma or serum using a Thermo Scientific™ TSO Quantum™ triple quadrupole tandem mass spectrometer.



- Quantitation workflows for a range of analytes
- Antidepressants, antimycotics, benzodiazepines, antiepileptics in human plasma
- Low LLOQs achieved for each analyte in complex matrix condition in an easy-to-implement workflow solution

TSQ Fortis and TSQ Quantis – Enabling Clinical Research

Productivity

Address high productivity demands for every clinical research applications with unmatched speed

Ease-of-Use

Achieve quality data regardless of user expertise with exemplary ease



Sensitivity

Address best-in-class sensitivity that addresses all requirements while improving price-performance ratio of the workflow solution

Robustness

High robustness offered by the triple quadrupole ensures very low variability for thousands of injections

For Research Use Only

Pharmaceutical Applications




Confident Quantitation for Pharmaceutical QA/QC Applications

Achieve robust, reliable, and reproducible assays of drug candidates in water and biological matrices with easy-to-implement quantitative solutions


- **Be Prepared to Address Challenges, Everyday**
 - Robust, reproducible, reliable methods for fast quantification of drug candidates
 - Address both small and large molecule types
- **Address Productivity Goals with Workflow Solutions**
 - **Robustness** – less maintenance, increased uptime
 - **Ease of use** – easy method setup, data review and customized reporting
 - **Speed** – analyze more compounds/injection or reduce overall run times
 - **Sensitivity** – detect contaminants with class leading sensitivity

TSQ Fortis: Demonstration of Sensitivity – Pharma QA/QC




- Demonstration of sensitivity, precision, accuracy and robustness for Pharma QA/QC.
- The same quantitation workflow can also be transferred to quantify Phenylephrine in biological matrices

Application Note 65200



TSQ Quantis: Sensitive Quantitation of Albuterol

Excellent Linearity and Reproducibility across the Dynamic Range






Experimental Details

UHPLC : Thermo Scientific™ Vanquish™ Flex Binary system Flow rate: 0.5 mL/min
Solvent A: 0.1% FA 4 mM ammonium formate in H₂O
Solvent B: 0.1% FA 4 mM ammonium formate in MeOH

Peptide Quantitation with TSQ Quantis Triple Quadrupole MS

Experimental Details


UHPLC : Thermo Scientific™ Easy nLC 1200 system
Flow rate: 300 nL/min
Solvent A: 2%ACN 0.1%FA in H₂O
Solvent B: 10%H₂O 0.1 %FA in ACN



Consistent Ion Ratios
Across Wide Mass
Range for Confident
Peptide Quantitation

TSQ Altis: Fluticasone Propionate


Superior Sensitivity and Reproducibility - Fluticasone Propionate in matrix!
25 fg on column with 4.2% CV!



Experimental Details

UHPLC : Thermo Scientific™ Vanquish™ Flex Binary system
Flow rate: 0.75 mL/min
Solvent : 75% MeOH:25% H_2O with 0.01% NH_4OH
MS: TSQ Altis Triple Quadrupole MS
Software: TraceFinder Software 4.1

TSQ Altis: Protein Quantitation – Benefits of Higher Selectivity and Sensitivity



Experimental Details




LC : Thermo Scientific™ Easy nLC 1200
 Flow rate: 300nL/min
 Solvent A: 2%ACN with 0.1%FA in H₂O
 Solvent B: 10%H₂O with 0.1 %FA in ACN

Quantitation of Mixture of Large Molecules





- Column: 2.1 x 50 mm, 1.5 μ m Accucore Vanquish C18
- Column Temp: 60 C
- Mobile Phase: [A] H₂O + 0.1% Formic Acid; [B] ACN + 0.1% Formic Acid
- Injection Volume: 10 μ L
- Sample Temp: 10 C



- Ionization Mode: HESI, Positive ion modes
- MS Acquisition Mode: Selective Reaction Monitoring (SRM) – see table below
- Cycle Time: 0.8 s
- Quad Isolation (Q1,Q3) = Unit (1.2 Da FWHM)



Optimized LC/MS – C-Peptide, Parathyroid Hormone, and Insulin



Confident Quantitation – Today and Tomorrow

TSQ Altis



TSQ Quantis



TSQ Fortis

