



ThermoFisher
SCIENTIFIC



Raise the Bar

Application of QuEChERS and Solid Phase Extraction in Food Safety

Jenny Huang
Senior Product Manager

SEA and Taiwan


Aug. 30-31, 2018

- Analytical Challenges
- Meeting the Challenge
- QuEChERS applications
 - Multi residues pesticides in agriculture food
 - Multiple Mycotoxins in Grain
- SPE applications
 - Multiple veterinary drugs screening in animal products
 - Hormone residues in animal-derived food



Analytical Challenges

- Diverse and complex food matrices
- Increasingly large number of target compounds
- Low limits of detection
- High throughput (Efficient)
- Cost effective



How do we overcome the Challenges?

- Sample preparation
 - Removal of suppressants
 - Reduced matrix effects
- High chromatographic resolution
 - Accurate quantification
 - Accurate qualification
 - Speed in analysis
- Selectivity of detection
 - Accurate quantification
 - Accurate qualification
 - Speed in analysis
- Reproducible products



Commonly employed techniques:

- Filtration
- **Conventional SPE** – silica, polymeric, Hypercarb
 - Targeted extraction
- **QuEChERS**
 - Multi residue extraction
- Liquid/Liquid Extraction




Why Use SPE for Food Analysis ?

- Sample Enrichment
 - To achieve sub ppb detection limits, usually the sample must be concentrated 10-100 times
 - using only liquid extraction (LE), as much as **100 mL** of solvent must be evaporated
 - using SPE to concentrate the extract, only **a few mL** of solvent must be evaporated
- Sample Cleanup
 - Significant cleanup of food matrix can often be accomplished using SPE
 - Multiple SPE cartridges can often be used for optimum cleanup

SPE - Considerations

The choice of SPE mechanism (RP, NP, IEX, MM) depends upon;

- Physicochemical differences between the analyte and matrix
- Analyte solubility
- Analyte polarity (LogP/LogD)
- Analyte charge state (pKa analyte/pH of solution)



- Strategy 1: Retention, cleanup, elution, most effective cleanup and enrichment for individual compound class analysis
- Strategy 2: Pass-through cleanup, more effective cleanup for multiresidue analysis
- Strategy 3 : Dispersive cleanup, ex. QuEChERS, acceptable cleanup for multiresidue screening
 1. Combine sorbent, sample matrix and solvent into a vessel
 2. Sample is filtered or centrifuged. **Matrix interferences are retained by sorbent**
 3. Filtrate or supernatant is collected for analysis **Analytes are in the filtrate or supernatant**

What is QuEChERS (dispersive SPE)?

- QuEChERS or dispersive SPE is a sample preparation method which is typically employed for the extraction of pesticides from food matrices
- QuEChERS stands for:
 - **Q**uick
 - **E**asy
 - **C**heap
 - **E**ffective
 - **R**obust
 - **S**afe
- Typical formats are:
 - Sorbent pre packed in 50, 15 and 2mL centrifuge tubes
 - Sorbent pre-packed in mylar pouches
- Different sorbent configurations are designed to target specific food matrices



Sample Homogenisation Before Extraction

- Unavoidable steps in the analysis
- Prerequisite to obtain representative sub-samples
- Critically important



How does QuEChERS work?


QuEChERS is typically a 4 step process:

- Stage 1: Select the appropriate product based on the matrix type
- Stage 2: Extraction (liquid/liquid) - 50mL centrifuge tube
- Stage 3: Sample clean-up to remove interferences that can interfere with subsequent analysis - 2 and 15mL centrifuge tube
- Stage 4: Analysis - typically performed via GCMS (LCMS can also be present)




What is QuEChERS? *Extraction & Clean-up in a Tube*

1) Extraction



2) Dispersive SPE



Note: Add sample to the tube, then solvent, then sorbent then mix, to avoid agglomeration

QuEChERS Method Variations

| Method | Description |
|-----------------------------------------------|------------------------------------------------------------------------------------------------------------|
| Original QuEChERS Method – introduced in 2003 | Uses Sodium Chloride to enhance extraction |
| Dispersive AOAC 2007.01 Method | Uses Sodium Acetate as a buffer replacing Sodium Chloride |
| Dual Phase Variation | Uses PSA & GCB to remove high levels of chlorophyll and plant sterols |
| European Version | Similar to AOAC method – uses sodium chloride, sodium citrate dihydrate and disodium citrate sesquihydrate |

Stage 1: Selecting the Method

1. Weigh 15 g of homogenised (hydrated at least 80%) sample in a 50 mL centrifuge tube
2. Add 15 mL acetonitrile (or 1:1 acetone/hexane, ethyl acetate) and IS
3. Shake briefly
4. Add 4 g anhydrous magnesium sulfate and 1g sodium chloride
5. Shake by hand for 1 minute
6. Centrifuge at 5,000 rpm for 5 min
7. Transfer a portion of supernatant to a QuEChERS clean up tube
8. Shake for 30 sec
9. Centrifuge for 1 min at 6,000 rpm
10. Transfer 0.5 mL aliquot for analysis



Stage 1: Selecting the Method

1. Weigh 15 g of homogenised (hydrated at least 80%) sample in a 50 mL centrifuge tube
2. Add 15 mL 1% acetic acid in acetonitrile (or 1:1 acetone/hexane, ethyl acetate) and IS
3. Shake briefly
4. Add 6 g anhydrous magnesium sulfate and 1.5 g anhydrous sodium acetate
5. Shake by hand for 1 min
6. Centrifuge at 3,700 rpm for 5 min
7. Transfer a portion of supernatant to a QuEChERS clean up tube
8. Shake for 30 sec
9. Centrifuge for 1 min at 3,700 rpm
10. Transfer 0.5 mL aliquot for analysis and add TPP and either formic acid or toluene / magnesium sulfate

Stage 1: Selecting the Method

- Using sodium acetate as a buffer protects base sensitive compounds
 - Folpet
 - Dichlofluanid
 - Chlorthanonil
 - Dicofol
 - Captan
 - Tolyfluanid

AOAC Method –Disadvantages

- The presence of acetic acid
 - PSA will absorb acetic acid
 - Less sample clean-up
 - Higher baseline
- Only use this method if looking at specific compounds



Stage 1: Selecting the Method





Using dual layer PSA/GCB cartridge

- Extraction phase as AOAC method
- Pre-rinse the SPE cartridge with 5 mL of toluene
- Add an aliquot of the supernatant to the cartridge
- Start collection
- Elute with 6-12 mL of 3:1 acetone:toluene
- Concentrate for GCMS analysis –or–
- Concentrate to near dryness and reconstitute in mobile phase for LC analysis

Stage 1: Selecting the Method

1. Weigh 15 g of homogenised (hydrated at least 80%) sample in a 50 mL centrifuge tube
2. Add 15 mL acetonitrile (or 1:1 acetone/hexane, ethyl acetate) and IS
3. Shake briefly
4. Add 6 g anhydrous magnesium sulfate, 1.5 g sodium chloride, 1.5 g sodium citrate tribasic dihydrate, 0.75 g sodium citrate dibasic
5. Shake by hand for 1 min
6. Centrifuge at 5,000 rpm for 5 min
7. Transfer a portion of supernatant to a QuEChERS clean up tube
8. Shake for 30 sec
9. Centrifuge for 1 min at 6,000 rpm
10. Transfer 0.5 mL aliquot for analysis

Stage 2: Select the Right Product

| Matrix Type | Examples | Sorbent Requirements for Clean-Up |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|
| General Matrices | <ul style="list-style-type: none"> • Apples • Cucumber • Melon  | MgSO ₄ , PSA Removal of excess water organic acids, fatty acids, sugars |
| Fatty Matrices | <ul style="list-style-type: none"> • Milk • Cereals • Fish  | MgSO ₄ , PSA, C18 Additional removal of lipids & sterols |
| Pigmented Matrices | <ul style="list-style-type: none"> • Lettuce • Carrot • Wine  | MgSO ₄ , PSA, C18, GCB Additional removal of pigments & sterols |
| High Pigmented Matrices | <ul style="list-style-type: none"> • Spinach • Red Peppers  | MgSO ₄ , PSA, C18, GCB, Chlorofiltr™ Additional removal of chlorophyll |

Stage 2: Select the Right Product

| Material | Purpose |
|---------------------------------|-----------------------------------------------|
| Magnesium Sulphate | Removal of excess water |
| PSA (Primary / Secondary Amine) | Removal of organic acids, fatty acids, sugars |
| C18 | Removal of lipids & sterols |
| GCB (Graphitized Carbon Black) | Removal of pigments & sterols |
| Chlorofiltr™ | Removal of chlorophyll |

QuEChERS Dispersive Solid Phase Extraction for the GC-MS Analysis of Pesticides in Cucumber

Anila I Khan, Thermo Fisher Scientific, Runcorn, Cheshire, UK

Application Note 20349

Key Words

QuEChERS, pesticide residue analysis, cucumber, food safety

Abstract

QuEChERS dispersive SPE is a simple, fast and quantitative sample preparation method. This application demonstrates the effectiveness of this technique in the GC/MS analysis of pesticides in cucumber, using a Thermo Scientific TraceGOLD TG-5MS GC column for analysis.

The recoveries for the spiked pesticides in cucumber matrix at 50 ng/g were between 75.2 to 119.6% with relative standard deviations ranging from 2.1 to 8.9% using the QuEChERS method described in EN15662.

Introduction

QuEChERS (Quick, Easy, Cheap, Effective, Rugged and Safe) is a dispersive Solid Phase Extraction (SPE) technique for extracting multi-residue pesticides from fruits and vegetables. The advantages of this methodology are speed, ease of execution, minimal solvent requirement and cost. The QuEChERS methodology was developed by Anastassiades et al¹ and has become widely used in food safety analyses.

The method is:

- **Quick** – high sample throughput, typically 8 samples can be prepared in under 30 min
- **Easy** – it requires less handling of extracts than other techniques i.e. fewer steps are required
- **Cheap** – less sorbent material is needed and less time is required to process samples compared to other techniques
- **Effective** – the simple technique gives high and accurate recovery levels for a range of different compound types.
- **Rugged** – the method can detect a large number of pesticides including charged and polar pesticides
- **Safe** – unlike other techniques, it does not require the use of chlorinated solvents. Extraction is typically carried out using acetonitrile, which is both GC and LC compatible.




The sample preparation approach described in the European EN15662 QuEChERS procedure² was used for extracting pesticides from cucumber. This is a two stage process: sample extraction, followed by dispersive SPE.

In the sample extraction stage, the food sample is homogenized to increase the available surface area of the sample to provide optimal extraction efficiencies. The homogenized sample is placed in the extraction tube containing magnesium sulfate and salts (sodium chloride, sodium citrate tribasic dihydrate, sodium citrate dibasic sesquihydrate). Magnesium sulfate ensures that, upon addition of acetonitrile, a phase separation is induced between water and organic solvent with the pesticides of interest being extracted into the organic phase. When acetonitrile is poured into the extraction tube containing the homogenized sample,


- QuEChERS sample preparation
- Analysis using 5% diphenyl / 95% dimethylpolysiloxane GC phase

Sample Preparation

The methodology described in this application note is for the preparation of calibration standards and sample spike (Figure 1).



Application – Pesticides in Cucumber



| Pesticides | t_r (min) | Linearity | Nominal concentration ng/g | Measured concentration (n=6) ng/g. | Average % Recovery (n=6) | Recovery %RSD (n=6) |
|----------------------------|-------------|-----------|----------------------------|------------------------------------|--------------------------|---------------------|
| 1. Dichlobenil | 4.52 | 0.9988 | 50 | 58.7 | 117.4 | 2.1 |
| 2. Tribromoanisol | 6.60 | 0.9990 | 50 | 54.3 | 108.5 | 6.0 |
| 3. Sulfotep | 6.95 | 0.9984 | 50 | 59.8 | 119.6 | 2.3 |
| 4. Hexachlorobenzene | 7.49 | 0.9983 | 50 | 55.2 | 110.4 | 2.8 |
| 5. Parathion | 10.90 | 0.9979 | 50 | 53.0 | 106.0 | 5.9 |
| 6. Triphenylphosphine (IS) | 13.41 | - | - | - | - | - |
| 7. EPN | 17.90 | 0.9985 | 50 | 46.1 | 92.1 | 6.7 |
| 8. Azinphos methyl | 19.20 | 0.9984 | 50 | 37.6 | 75.2 | 4.9 |
| 9. Permethrin isomer a | 21.38 | 0.9987 | 50 | 49.8 | 99.5 | 8.9 |
| 10. Permethrin isomer b | 21.58 | 0.9985 | 50 | 50.9 | 101.9 | 4.8 |
| 11. Fenvalerate isomer a | 24.60 | 0.9985 | 50 | 47.7 | 95.4 | 8.9 |
| 12. Fenvalerate isomer b | 25.02 | 0.9973 | 50 | 50.6 | 101.2 | 7.2 |
| 13. Deltamethrin | 25.84 | 0.9949 | 50 | 51.7 | 103.3 | 6.8 |

- TraceGOLD TG-5MS GC base line resolution 12 pesticides
- Linear range of 25-1000ng/g (R^2 0.99)
- QuEChERS method:
 - Recoveries >75%
 - Reproducibility <9%RSD.

Application Note:
ANSC
PESTGRAPES 0709

Key Words

- QuEChERS
- Food Safety
- Pesticides Residue Analysis
- TRACE TR-5MS
- TRACEGuard

QuEChERS Dispersive Solid Phase Extraction for the GC/MS Analysis of Pesticides in Grapes

Anilá Khan, Laisa Pereira, Stephen Aspy, Rob Bunn, Ruth Lewis, Thermo Fisher Scientific, Runcorn, UK

Introduction

QuEChERS (Quick, Easy, Cheap, Effective, Rugged and Safe) is a dispersive Solid Phase Extraction (SPE) technique for extracting multi-residues of pesticides from fruits and vegetables. The advantages of this methodology are speed, ease of execution, minimal solvent requirement and cost to perform when compared with conventional solid phase extraction techniques.

The QuEChERS methodology was developed by Anastassiades *et al* and has become widely used in food safety analyses.¹ The method is:

- Quick – high sample throughput, typically 8 samples can be prepared in just under 30 min
- Easy – it requires less handling of extracts than other techniques and no laborious steps are involved
- Cheap – less sorbent material is needed and less time is needed to process samples compared to other techniques
- Effective – the simple technique gives high and accurate recovery levels for a range of different compound types e.g. polar pesticides, pH dependent compounds
- Rugged – the method can detect a large number of pesticides including pH dependent and polar pesticides
- Safe – unlike other techniques, it does not require any chlorinated solvents. Extraction is typically carried out using acetonitrile, which is both GC and LC amenable

The QuEChERS procedure is usually a two stage process: sample extraction, followed by dispersive SPE. In the sample extraction stage, the food sample is homogenized to maximize the available surface area of the sample for better extraction efficiencies. The homogenized sample is placed in the extraction tube containing magnesium sulphate and sodium acetate. Magnesium Sulfate ensures that upon addition of acetonitrile, a phase separation is induced between water and organic solvent with the pesticides of interest being extracted into the organic phase. When acetonitrile is poured into the extraction tube containing the homogenized sample, an exothermic reaction occurs between the magnesium sulphate and water, which can lead to low recoveries of the pesticides. This effect can be reduced by adding the salt and sample to the extraction tube while this is immersed in an ice bath or by weighing the sample into a FEP tube and then adding the solvent and salts slowly. The tube is then capped, shaken vigorously and centrifuged. The second stage of the QuEChERS method uses dispersive SPE, which involves transferring a portion of the acetonitrile extract to a clean-up tube containing a combination of sorbents for removal of unwanted sample components. This may be followed by solvent exchange to improve compatibility of samples to

GC analysis, and additional sample clean-up to reduce matrix effects and therefore improve method robustness.

Internal standards are used to minimize errors that might be introduced in the different steps of the QuEChERS method, as well as compensate for GC injection variability. Furthermore, adding analyte protectants such as sorbitol can be useful for labile pesticides at intermediate pH, which can be prone to decomposition in the GC injector port.

The pesticides analyzed are a mixture of organophosphate, organochlorine, pyrethroid, benzenoid, triazole and dicarboximide compounds. Lehotay reviewed the LC and GC analyses of pesticides in produce and the type of pesticide that is likely to be found in each matrix.² The requirements for pesticide residue analysis in fruit and vegetables are established by organizations such as World Health Organization, Japanese Food Chemical Research Foundation, EEC Directives, and the US-EPA.^{3,4,5,6} These organizations establish which pesticides need to be determined in different produce and the Method Regulatory Limits (MRLs). The pesticides determined in this study are all listed by the four regulatory organizations and all have minimum MRLs of 50 ng/g (ppb). The recoveries of the pesticides in grapes are based on this value.

Goal

To demonstrate QuEChERS dispersive SPE as a simple, fast and quantitative sample preparation method for the GC/MS analysis of pesticides in grapes.

Additionally, demonstrate the suitability of the Thermo Scientific TRACE TR5-MS analytical column combined with the Thermo Scientific TRACEGuard guard column for pesticide analysis.

Experimental

Sample Preparation

Reagents

- Green grapes obtained from the local supermarket
- Acetonitrile HPLC grade (Thermo Fisher P/N A0626V17)
- Hexane GC grade (Thermo Fisher P/N H0355V17)
- Acetone GC grade (Thermo Fisher P/N A0600V17)
- Glacial acetic acid HPLC grade (Thermo Fisher P/N A0406V815)



| Pesticide (peak number, name) | Rt (min) | Linearity (R ²) | Average % recovery (n=3) | Recovery % RSD |
|-------------------------------|----------|-----------------------------|--------------------------|----------------|
| 1. Chlorpyrifos methyl | 14.92 | 0.9961 | 76 | 14.5 |
| 2. Metalaxyl | 15.40 | 0.9908 | 100 | 7.9 |
| 3. Malathion | 15.94 | 0.9989 | 103 | 19.0 |
| 4. Chlorpyrifos | 16.12 | 0.9919 | 77 | 21.5 |
| 5. Penconazole | 17.23 | 0.9978 | 94 | 7.2 |
| 6. Procymidon | 17.45 | 0.9968 | 106 | 4.0 |
| 7. Dicofol | 20.21 | 0.9964 | 110 | 9.9 |
| 8. Permethrin Isomer a | 21.42 | 0.9980 | 102 | 3.9 |
| 9. Permethrin Isomer b | 21.55 | 0.9965 | | |

organophosphate, organochlorine, pyrethroid, benzenoid, triazole recoveries of greater than 76% & and RSD's less than 11%




Figure 2: TIC for the GC/MS analysis of grapes spiked with 1 ng/µL of each pesticide

Thermo Fisher Scientific Products

| Step | Description | Capacity (mL) | Cat. No. | Quantity |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------|----------|
| Original Method | | | | |
| Extraction | 4000mg magnesium sulfate, 1000mg sodium chloride | 50 | 60105-211 | 250 Pack |
| Clean-Up | 150mg magnesium sulfate, 50mg PSA, 50mg C18 | 2 | 60105-204 | 100 Pack |
| | 150mg magnesium sulfate, 50mg PSA, 50mg C18, 50mg GCB | 2 | 60105-223 | 100 Pack |
| | 1200mg magnesium sulfate, 400mg PSA, 400mg C18 | 15 | 60105-225 | 50 Pack |
| AOAC 2007.1 | | | | |
| Extraction | 6000g magnesium sulfate, 400mg PSA, 400mg C18, 400mg GCB | 50 | 60105-210 | 250 Pack |
| Clean-Up | 150mg magnesium sulfate, 50mg PSA | 2 | 60105-203 | 100 Pack |
| | 150mg magnesium sulfate, 50mg PSA, 50mg C18 | 2 | 60105-204 | 100 Pack |
| | 150mg magnesium sulfate, 50mg PSA, 50mg C18, 50mg GCB | 2 | 60105-223 | 100 Pack |
| | 900mg magnesium sulfate, 300mg PSA, 150mg GCB | 15 | 60105-205 | 50 Pack |
| | 900mg magnesium sulfate, 300mg PSA, 150mg C18 | 15 | 60105-206 | 50 Pack |
| | 1200mg magnesium sulfate, 400mg PSA | 15 | 60105-224 | 50 Pack |
| | 1200mg magnesium sulfate, 400mg PSA, 400mg C18 | 15 | 60105-225 | 50 Pack |
| | 1200mg magnesium sulfate, 400mg PSA, 400mg C18, 400mg GCB | 15 | 60105-226 | 50 Pack |
| European EN15662 | | | | |
| Extraction | 6000g magnesium sulfate, 1500mg sodium chloride, 1500mg sodium citrate tribasic dihydrate, 750mg sodium citrate dibasic sesquihydrate | 50 | 60105-212 | 250 Pack |
| | 4000mg magnesium sulfate, 1000mg sodium chloride, 1000mg sodium citrate tribasic dihydrate, 500mg sodium citrate dibasic sesquihydrate | 50 | 60105-216 | 250 Pack |
| Clean-Up | 150mg magnesium sulfate, 25mg PSA | 2 | 60105-219 | 100 Pack |
| | 150mg magnesium sulfate, 25mg PSA, 2.5mg GCB | 2 | 60105-221 | 100 Pack |
| | 150mg magnesium sulfate, 25mg PSA, 7.5mg GCB | 2 | 60105-222 | 100 Pack |
| | 900mg magnesium sulfate, 150mg PSA | 15 | 60105-215 | 50 Pack |
| | 900mg magnesium sulfate, 150mg PSA, 45mg GCB | 15 | 60105-217 | 50 Pack |
| | 900mg magnesium sulfate, 150mg PSA, 15mg GCB | 15 | 60105-218 | 50 Pack |
| | 900mg magnesium sulfate, 150mg PSA, 150mg C18 | 15 | 60105-227 | 50 Pack |
| Dual Phase Method | | | | |
| | 200mg GCB on top, 400mg PSA on bottom, separated by a frit | 6 | 60105-207 | 30 Pack |
| | 250mg GCB on top, 500mg PSA on bottom, separated by a frit | 6 | 60105-208 | 30 Pack |
| | 500mg GCB on top, 500mg PSA on bottom, separated by a frit | 6 | 60105-209 | 30 Pack |

HyperSep Dispersive SPE Multipacks

| Description | Cat. No. | Quantity |
|----------------------------------------------------------------------------------------------------------------------------------------|-----------|----------|
| 4000mg anhydrous magnesium sulfate, 1000mg sodium chloride | 60105-332 | 50 Pack |
| 4000mg anhydrous magnesium sulfate, 1000mg sodium chloride, 500mg sodium citrate dibasic sesquihydrate, 1000mg sodium citrate tribasic | 60105-333 | 50 Pack |
| 4000mg anhydrous magnesium sulfate, 1000mg sodium acetate | 60105-334 | 50 Pack |
| 6000mg anhydrous magnesium sulfate, 1500mg sodium acetate | 60105-335 | 50 Pack |
| 6000mg anhydrous magnesium sulfate, 1500mg sodium chloride | 60105-336 | 50 Pack |
| 6000mg anhydrous magnesium sulfate, 1500mg magnesium sulfate, 1500mg sodium citrate dihydrate, 750mg disodium citrate sesquihydrate | 60105-337 | 50 Pack |
| 8000mg anhydrous magnesium sulfate, 2000mg of sodium chloride | 60105-338 | 50 Pack |
| 8000mg anhydrous magnesium sulfate, 3500mg of sodium chloride | 60105-339 | 50 Pack |



Multiple Mycotoxins in Grain Using QuEChERS

Determination of Multiple Mycotoxins in Grain Using a QuEChERS Sample Preparation Approach and LC-MS/MS Detection

Jon Bardsley, Mike Oliver, Thermo Fisher Scientific, Runcorn, UK

Application Note 21121

17 Mycotoxins in Grain on Accucore aQ

Key Words

Mycotoxins, food, HyperSep, QuEChERS, dispersive SPE, Accucore aQ, TSQ Vantage

Goal


To demonstrate a fast, easy, and cost-effective approach for the determination of 16 mycotoxin residues in grain-based food using QuEChERS sample preparation with a Thermo Scientific™ Accucore™ aQ HPLC column and a Thermo Scientific™ TSQ™ Vantage™ triple quadrupole mass spectrometer for HPLC separation.

Introduction


Mycotoxins are toxic secondary metabolites that are naturally produced by several species of fungi on agricultural products, particularly grain-based products. Mycotoxins are chemically stable and cannot be destroyed during food processing and heat treatment. Therefore, they may occur in the field, in raw materials during storage, and in processed foods. To date, more than 300 mycotoxins, possessing varying degrees of toxicity, have been identified, although only a relatively few of these are widely accepted as presenting a significant food or animal feed safety risk.¹ Studies have shown that




limited, final extracts may still contain large amounts of matrix components that can negatively affect the detection system.




AFB₁




AFB₂




AFM₁



AFG₁



AFG₂



AFB_{2a}

Multiple Mycotoxins in Grain Using QuEChERS

• Sample extraction

1. Place 5 g of sample into a 50 mL centrifuge tube.
2. Add 10 mL water.
3. Vortex briefly and allow to hydrate for at least 15 minutes.
4. Add 250 μ L of 1 μ g/mL internal standard. Add 10 mL acetonitrile containing 2% formic acid.
5. Shake/vortex samples for 15 minutes to extract mycotoxins.
6. Add the contents of a HyperSep Dispersive SPE Mylar Pouch (P/N 60105-340) containing the extraction salts.
7. Immediately shake for 1 minute.
8. Centrifuge for 5 min at $\geq 3000 g$.

• Sample cleanup

1. Transfer 1 mL of supernatant to a dSPE tube (HyperSep Dispersive SPE Clean-up Product **P/N 60105-204**). 150mg anhydrous magnesium sulfate, 50mg PSA, 50mg C18
2. Vortex for 30 seconds.
3. Centrifuge for 5 min at $\geq 3000 g$.
4. Transfer 500 μ L of purified supernatant to a 5 mL test tube.
5. Evaporate the acetonitrile extract to dryness and reconstitute with 500 μ L methanol / water (50:50, v/v) for better chromatographic performance.
6. Filter the extract, using a 0.2 μ m syringe filter (P/N F2513-2), directly into an autosampler vial (P/N 60180-600).

Multiple Mycotoxins in Grain Using QuEChERS

For improved sensitivity at low ppb concentrations, the dSPE step can be scaled-up by following the steps below:

1. Transfer 8 mL of supernatant to a 15 mL HyperSepDispersive SPE Clean-up Product (1200 mg anhydrousMgSO₄, 400 mg PSA, and 400 mg C18, 15 mL centrifuge tube, P/N 60105-204).
2. Vortex for 30 seconds.
3. Centrifuge for 5 minutes at ≥ 3000 g.
4. Transfer 5 mL of supernatant to a 5 mL test tube.
5. Evaporate the sample to dryness at 40–50 ° C under a gentle stream of nitrogen.
6. Reconstitute sample in 1 mL methanol / water (50:50, v/v).

Multiple Mycotoxins in Grain Using QuEChERS

| | |
|-------------------------|------------------------------------------------------------------------------------------------|
| Instrumentation | Thermo Scientific™ Dionex™ UltiMate™ 3000 |
| Column | Accucore aQ, 100 × 2.1 mm, 2.6 μm (P/N 17326-102130) |
| Guard column | Accucore aQ, 10 × 2.1 mm, 2.6 μm (P/N 17326-012105) |
| Run time | 17 min (including 4 min equilibration) |
| Temperature | 45 °C |
| Injection volume | 5 μL |
| Autosampler temperature | 10 °C |
| Wash solvent | Methanol / water (1:1, v/v) |
| Flow rate | 400 μL/min |
| Waste divert | Mobile phase was diverted to waste from 0–1.5 and 13–17 min to reduce ion source contamination |
| Mobile phase A | 10 mM ammonium formate in water |
| Mobile phase B | Methanol |

LC gradient.

| Time (min) | A (%) | B (%) |
|------------|-------|-------|
| 0.0 | 100 | 0 |
| 1.0 | 75 | 25 |
| 4.0 | 75 | 25 |
| 5.0 | 60 | 40 |
| 8.0 | 60 | 40 |
| 8.5 | 40 | 60 |
| 9.5 | 40 | 60 |
| 10.0 | 0 | 100 |
| 13.0 | 0 | 100 |
| 13.2 | 100 | 0 |
| 17.0 | 100 | 0 |

Multiple Mycotoxins in Grain Using QuEChERS



Accuracy and precision data for the 16 mycotoxins fortified at two concentrations.

| | 20 µg/kg | | 100 µg/kg | |
|------------------------|--------------|---------|--------------|---------|
| | Recovery (%) | RSD (%) | Recovery (%) | RSD (%) |
| Nivalenol | 71.4 | 11.2 | 67.2 | 6.5 |
| Deoxynivalenol | 106.7 | 4.1 | 97.0 | 2.8 |
| 3-Acetyldeoxynivalenol | 100.4 | 3.9 | 97.2 | 1.9 |
| Fusarenon X | 96.3 | 3.9 | 96.2 | 3.8 |
| Neosolaniol | 100.5 | 3.3 | 99.4 | 2.0 |
| Diacetoxyscirpenol | 102.6 | 2.8 | 99.0 | 2.3 |
| Altemariol | 94.8 | 4.9 | 85.9 | 5.4 |
| β-zearalanol | 94.5 | 9.2 | 92.7 | 4.6 |
| α-zearalanol | 93.9 | 10.5 | 89.0 | 3.5 |
| Zearalenone | 92.4 | 9.4 | 87.6 | 4.5 |
| Ochratoxin A | 93.8 | 3.0 | 94.7 | 3.8 |
| T-2 toxin | 96.2 | 4.5 | 94.2 | 2.8 |
| Aflatoxin B1 | 97.0 | 2.7 | 91.7 | 5.4 |
| Aflatoxin B2 | 97.4 | 2.9 | 91.4 | 4.8 |
| Aflatoxin G1 | 95.0 | 3.3 | 92.0 | 4.1 |
| Aflatoxin G2 | 95.5 | 3.1 | 93.9 | 2.7 |

Matrix-matched calibration curves of the 16 mycotoxins were prepared at concentrations of 5, 12.5, 25, 50, and 100 ng/mL and were found to give linear responses over the entire concentration range with correlation coefficients (R²) typically greater than 0.995. The signal-to-noise ratio (S/N) at the lowest calibration level (5 ng/mL; 10 ng/g) was found to be >10 for all 16 compounds. Therefore, the LOQ was estimated to be ≤10 ng/g in this study.

Animal-Derived Food multiple vet. drug residue analysis


Pork,
Chicken,
liver, duck



HyperSep Retain PEP pass through application for multiple vet. drug residue analysis

Accurately weigh 2.5g sample into a 50mL centrifuge tube, first add 1mL water to mix, then add 9mL containing 0.5% formic acid-acetonitrile, vortex to disperse and mix; Hold homogenizer for 1~2min

If homogenizer is not available, it is recommended to shake for 20 min instead. Ultrasonic extraction at 40 C for 15 min, centrifugation at 4 C for 12000rpm/min for 5 min. The supernatant is to be purified by SPE.



SPE pass through application for multiple vet. drug residue analysis

- Column: **Accucore RP-MS (2.1 x 100 mm, 2.6 μm)**
- Positive mode :
 - mobile phase A: water (0.1% formic acid);
 - mobile phase B: acetonitrile (0.1% A) Acid);
- Negative mode :
 - mobile phase A: water (0.03% ammonia);
 - mobile phase B: acetonitrile (0.03% ammonia)
- Flow rate: 0.3 mL/min;
- Injection volume: 5 μL;
- Column temperature: 30 ° C.
- The mobile phase elution gradient is shown in Table


| 时间 /min | A | B | 流速: μL/min |
|---------|----|----|------------|
| 0.00 | 95 | 5 | 300 |
| 15.0 | 5 | 95 | 300 |
| 17.0 | 5 | 95 | 300 |
| 17.1 | 95 | 5 | 300 |
| 20.0 | 95 | 5 | 300 |

| 参数名称 | 设置 |
|----------------------|--------------------|
| 喷雾电压 | 3500/3000v (+/-) |
| 雾化温度 | 400°C |
| 离子传输管温度 | 350°C |
| 鞘气 | 40 |
| 辅助气 | 5 |
| 反吹气 | 1 |
| Full Scan Resolution | 70000 |
| Full Scan Mass Range | 100-1000m/z |
| ddMS2 Resolution | 17500 |
| MS Isolation | 2.0 m/z |

- Mass spectrometry conditions: ESI source, unlike triple quadrupole, Orbitrap high-resolution mass spectrometry does not require optimization of mass spectrometry conditions and collision energy for each compound, so the time spent on method development is greatly reduced. The ion source parameters and mass spectrometry scan parameters are shown in Table

SPE pass through application for multiple vet. drug residue analysis

Comparison of the clean up effects of the three sample prep methods



The data indicating that the sample background and lipid interferences purified by the HyperSep Retain-PEP SPE column are obvious.

HyperSep Retain-PEP can effectively remove most of the fat and phospholipids, and the use of pass through strategy sample purification method can greatly improve sample processing.

SPE pass through application for multiple vet. drug residue analysis

对猪肉、猪肝、鸡肉、鸭肉等基质中添加了194种常见兽药，加标浓度为1.0 μ g/Kg和5.0 μ g/Kg的回收率做了统计，

回收率基本在50-120%之间，以下列出了部分化合物的回收率数据结果。

| 英文名称 | 中文名称 | Add 1ppb | | Add 5ppb | |
|-----------------------------------------|----------|----------|-------------|----------|-------------|
| | | 测定值 | Recovery(%) | 测定值 | Recovery(%) |
| 17-Methyltestosterone | 甲基睾酮 | 0.68 | 68% | 3 | 60% |
| 2-Aminoflubendazole | 2-氨基氟苯咪唑 | 0.74 | 74% | 4.2 | 84% |
| 4-EPICHLORTETRACYCLINE HYDROCHLORIDE | 差向金霉素 | 0.62 | 62% | 3.3 | 66% |

第 7 页

The results showed that the addition concentration was 1.0 μ g/Kg, and the detection rate of all compounds in the matrix was >84%; for addition 5.0 μ g/Kg and the detection rate of all compounds in the matrix was >92%. The compound recovery is between 50-120%.

Sample prep for hormone residues in Animal-Derived food

Weigh 5 g sample, add 3 mL of 10 % sodium bicarbonate and 10 mL ethyl acetate, Homogenized 30 s. Shake 10 min (4°C), centrifuge @6000 r/min 10 min, transfer upper organic layer to ashigatabin. Repeat 10 mL Ethyl Acetate extraction again. combine upper organic layer and evaporate to dryness under 40°C, and reconstituted with 30%MeOH and dilute to 5mL

Thermo HyperSep C₁₈
500mg,3mL,60108-304

Condition and Equilibrium
3 mL MeOH · 3 mL Water

Load sample
1 mL/min

Wash
2 mL 30% MeOH

Elute
5 mL MeOH

Evaporate and reconstitution
1 mL 50% Water/ACN

Thermo Hypersil Gold C₁₈
2.1x150mm,5µm, 25005-152130

Sample prep for hormone residues in animal-derived food

- LC Conditions
- Column: Thermo Hypersil Gold (150 × 2.1 mm, 5 μm)
- Gradient see Table 1, A: 0.1%FA B:MeOH
- Flow Rate: 250 μL/min
- Injection Volume: 10 μL

表1 流动相梯度洗脱条件

| Time (min) | A (%) | B (%) |
|------------|-------|-------|
| 0 | 90 | 10 |
| 5 | 10 | 90 |
| 12 | 10 | 90 |
| 12.1 | 90 | 10 |
| 14 | 90 | 10 |

- **MS Conditions:**
- ESI: Positive mode
- SRM: Scan mode,
- Spray volt: 4500V;
- Ion transfer temp : 350 °C;

表2 11种合成类固醇类激素的SRM质谱条件

| | 药物名称 | 母离子 | 碎片离子 (碰撞能量 V) |
|----|-------------------------------------|-------|-----------------------------------|
| 1 | 睾酮 (Testosterone) | 289.2 | 97.2 (24), 109.2 (26), 253.2 (16) |
| 2 | 甲基睾酮 (Methyl testosterone) | 303.2 | 97.2 (30), 109.1 (28), 285.1 (17) |
| 3 | 黄体酮 (Progesterone) | 315.2 | 97.2 (22), 109.1 (27) |
| 4 | 群勃龙 (Trenbolone) | 271 | 199 (26), 253 (20) |
| 5 | 勃地龙 (Boldenone) | 287.2 | 121 (25), 135 (15), 269 (12) |
| 6 | 诺龙 (Nandrolone) | 275.2 | 109 (28), 257 (15) |
| 7 | 美雄酮 (Methandienone) | 301 | 121 (27), 149 (15), 283 (12) |
| 8 | 司坦唑醇 (Stanozolol) | 329.2 | 81 (42), 121 (36), 95 (38) |
| 9 | 丙酸诺龙 (Nandrolone propionate) | 331.1 | 239 (17), 257 (14) |
| 10 | 丙酸睾酮 (Testosterone propionate) | 345.2 | 97.1 (29), 109.1 (28) |
| 11 | 苯丙酸诺龙 (Nandrolone phenylpropionate) | 407.2 | 105.2 (29), 257.1 (17) |

Sample prep for hormone residues in Animal-Derived food

LC/MS/MS 色谱图




图2 11种合成类固醇类激素LC/MS/MS色谱图 (1.0 ng/mL)

This method can be applied to the extraction of synthetic steroid hormone in pork, beef, mutton, chicken, egg, milk etc. animal-derived foods . The recovery is in the range of 50 - 105 %.

Thermo Scientific Food Communities: Resources

- View application notes, on-demand webinars, product information, and many more resources on our Pesticides and Food Communities Libraries:

www.thermoscientific.com/pesticides

www.thermoscientific.com/foodandbeverage

The screenshot displays the Thermo Scientific Food and Beverage Community website. On the left, a navigation menu lists various categories under 'Food and Beverage', including Chemical and Elemental Analysis, Authenticity and Adulteration, Environmental Pollutants, Natural Toxins and Biotoxins, Nutritional and Food Label Testing, Packaging and Food Contact Materials, Pesticide Residues, Processing Contaminants, Veterinary Drugs, Dietary Supplements and Nutraceuticals Analysis, Fermentation, Microbiology Testing, Packaged Product Inspection, and Process Analysis and Control. The main content area features a green header with the text 'Welcome to your Food and Beverage Community'. Below this, there are tabs for 'Application notes' and 'Sample Preparation'. A list of application notes is visible, including 'AN 1097: Determination of P... Extraction and Gas Chromat...', 'AB 152: Extraction of Organic Solvent Extraction', 'Food Safety Applications No...', 'Accelerated Solvent Extracti... Chlorinated Compounds, PO...', 'Determination of Pesticides i... Sample Preparation LC-MS/I...', 'QuEChERS Dispersive Solid Cucumber', 'Multi-Residue Pesticide Anal... a Triple Quadrupole GC-MS/I...', 'Analysis of Pesticide Residu... and Single Quadrupole GC/M...', and 'Multi-residue Pesticide Analysis in Green Tea by a Modified QuEChERS Extraction and Ion...'. To the right, there is a section for 'Pesticide Residues' with a green background and a basket of apples. Text reads 'Diverse pesticide compounds demand an array of testing technologies'. Below this is a 'Featured webinars' section with a video thumbnail and the text 'Top 9 Pesticide Residues Webinars Learn strategies and solutions from food safety testing experts for the successful determination of pesticides in food.' There is also a 'Latest application notes' section with a video thumbnail and the text 'Three-fold Increase in Productivity for Pesticide Residue Analysis in Baby Food Using Fast Triple Quadrupole GC-MS/MS Learn how a fast, easy, and robust workflow was used to analyze pesticide residues in baby food.' At the bottom right, there is a 'Pesticide Analysis Center of Excellence' section with a globe icon and the text 'The Pesticide Analysis Center of Excellence is a consolidation of our expertise and technologies into a resource to guide through even the most challenging of your pesticides analyses.' A 'Visit Now' button is present. Below that is a 'Tips for Pesticide Analysis' section with the text 'Improve pesticide analysis'.