

The Powerful Technology for Genomic Analysis

27th March 2018, PSU

- INTRODUCTION
- Sbeadex : DNA extraction kit
- KASP : method for genotyping
- MASS ARRAY : Multi-gene Analysis

- **INTRODUCTION**
- Sbeadex : DNA extraction kit
- KASP : method for genotyping
- MASS ARRAY : Multi-gene Analysis

About LGC



LGC is an international leader in the extended life sciences sector, including human healthcare, agri-food & the environment.

We provide a comprehensive range of reference materials, proficiency testing schemes, genomics reagents and instrumentation, as well as research and measurement services. Our scientific tools and solutions enable organisations to advance research, develop new products and form an essential part of their quality and compliance procedures.

Our 2,300 employees include internationally-recognised scientists who are experts in their field. LGC is headquartered in London and operates out of 19 countries worldwide. We are extensively accredited to quality standards such as GMP, GLP, ISO 13485, ISO 17034, ISO 17043, ISO/IEC 17025 and ISO 9001.

LGC has been home to the UK Government Chemist for more than 100 years and is the UK National Measurement Laboratory and Designated Institute for chemical and bio measurement.

LGC has been privately-owned since 1996 and has diversified through internal investment and acquisition to be an international leader in its chosen niche markets. LGC is now owned by funds affiliated with KKR.

We help customers conform to international statutory and regulatory standards. Science is at the heart of all we do - for a safer world.



LGC has been home to the UK Government Chemist for more than 100 years and is the UK National Measurement Laboratory and Designated Institute for chemical and bio measurement... Privately-Owned since 1996.

AGENA BIOSCIENCE : COMPANY AT A GLANCE



SEQUENOM[®]

(BEFORE 2014)



HQs : San Diego, CA


- INTRODUCTION
- **Sbeadex** : DNA extraction kit
- KASP : method for genotyping
- MASS ARRAY : Multi-gene Analysis

Sbeadex : DNA extraction kit for whom seeking both high yield and purified DNA.

Science
for a safer world




- INTRO TO DNA
- Workflow overview
- Summary



To get pure DNA: an underestimated challenge

DNA is part of a complex,
multicomponent cell.



High molecular weight compounds

- Polysaccharides
- Polypeptides
- Polyphenols
- Lipids
- Others

Low molecular weight compounds

- Dye stuff
- Secondary metabolites
(alcaloides, terpenes)
- Bivalent metal ions
- Preservatives
- Others

→ Storage of samples:

frozen, fresh & storage buffers (ie. Oragene)

→ Anti-Coagulants (blood):

EDTA; Citrate; Heparine

→ Environmental: young, old, illness, etc.


→ Time & throughput


→ Sample weight & volume

→ Human impact



- INTRO TO DNA
- Workflow overview
- Summary







Columns




Magnetic beads

Common workflow:







Steal balls:
3, 4, 5 mm




+/- 65°C
>/= 10 min



10 min.;
2500 x g



Kleargene spin plates: Columns




Kleargene feature	Specification
Extraction technology	Spin columns
Binding material & technology	Silica → Polarity binding
Final extraction wash buffer	(70%) Ethanol
Extraction format	96 & 384
Starting lysate volume	96 format : 300 µl 384 format: 200 µl
Automation	Semi-automated; i.e. Genespin
Sample material	Plant samples & Rodent tails
Application	Genotyping




sbeadex: Backbone technology




sbeadex magnetic
particle




Lysis




Bind 1





Wash 1



Bind 2



Wash 2
(water)




Elution







Lysis



Add beads
to lysis






Bind DNA to
beads



Collect DNA on
beads

Microplate with 96 samples of DNA bound to sbeadex

Set of magnetic rods



Manual extractions possible.

>96 Samples?

Automation.

sbeadex feature	Specification
Extraction technology	Proprietary magnetic beads
Binding material & technology	2-step Binding mechanism; Silica & DNA adaptor
Final extraction wash buffer	Pure water
Extraction format	96
Starting lysate volume	sbeadex mini: 20 -30 mg sbeadex maxi: 80 -100 mg
Automation	Fully-automated: oKtopure; KingFisher; open liquid handlers
Sample material	Plant, Livestock, Forensics, Plasmids
Application	All: KASP; Sequencing; arrays; PCR...

Product	Lysis	Bind 1	Wash 1	Bind 2	Wash 2	Wash 3	Dry	Elution
sbeadex	→	→	→	→	Water			
Kleargene	→	→	→		EtOH	→	→	
Qiagen	→	→	→		EtOH		→	

→ Unique 2-step binding mechanism

→ Final water wash

→ No Drying step




Plant samples

- Grains
- Vegetables
- Oilplants
- Fruits & trees
- Sample types:
 - Leaves
 - Seeds
 - Flowers
 - Roots, etc.



Livestock samples

- Blood & tissue
- Hair
- Semen & others
- Species:
 - Cattle, pig, sheep, goat, horses
 - Horses, insects etc
 - Chicken & other birds
 - Fish



Human samples

- Low volume blood (up to 2 ml)
- High volume blood (up to 10 ml)
- Tissue
- Special applications
 - Forensics
 - Hair
 - Stool



85%

AgBio-


-Plants- & livestock-



Forensic

5%

- Police and special applications-



5%

Human

- Blood, cohorts & diagnostics-



Others


5%

- (Soil) bacteria, plasmids, etc. -

sbeadex: Validated protocols for

Plant species	Leaves	Seeds	Plant species	Leaves	Seeds
Apricot (<i>Prunus armeniaca</i>)	✓		Onion (<i>Allium cepa</i>)	✓	✓
Barley (<i>Hordeum vulgare</i>)	✓	✓	Parsley (<i>Petroselinum crispum</i>)	✓	✓
Beet, sugar (<i>Beta vulgaris</i>)	✓		Peach (<i>Prunus persica</i>)	✓	
Canola / Oilseed (<i>Brassica napus</i>)	✓	✓	Pepper (<i>Capsicum annuum</i>)	✓	✓
Chicory (<i>Cichorium intybus</i>)	✓		Potato (<i>Solanum tuberosum</i>)	✓	
Corn (<i>Zea mays</i>)	✓	✓	Rice, Asian (<i>Oryza sativa</i>)	✓	✓
Cotton (<i>Gossypium</i>)	✓	✓	Rubber (<i>Hevea brasiliensis</i>)	✓	✓
Cucumber (<i>Cucumis sativus</i>)	✓	✓	Soybean (<i>Aphis glycines</i>)	✓	✓
Flax (<i>Linum usitatissimum</i>)	✓		Sunflower (<i>Helianthus annuus</i>)	✓	✓
Grape (<i>Vitis vinifera</i>)	✓	✓	Tobacco leaves (<i>Nicotiana tabacum</i>)	✓	✓
Lettuce (<i>Lactuca sativa</i>)	✓		Tomato (<i>Solanum lycopersicum</i>)	✓	✓
Muskmelon (<i>Cucumis melo</i>)	✓	✓	Wheat (<i>Triticum L.</i>)	✓	✓






One protocol – all livestock samples




- INTRO TO DNA
- Workflow overview
- Summary






Sbeadex

- **Targets:** plant & livestock
- **High DNA quality**
- **All down-stream applications**
- **Flexibility:** allows protocol optimizations
- **oKtopure**



Kleargene

- **Targets:** plant & rodent tails
- **Highest throughputs:** up to 20 k extr./ day
- **96/ 384 format**
- **Genotyping/ PCR**
- **Genespin**



Instru- mentation

- **oKtopure & Genespin**
- **Workhorses** for DNA extraction
- **Full and semi-automated** extractions
- **Low costs for plastics** (tips and plates)

- INTRODUCTION
- Sbeadex : DNA extraction kit
- **KASP** : method for genotyping
- MASS ARRAY : Multi-gene Analysis

KASP : A new era of KASP method for genotyping


Science
for a safer world



- Introduction to KASP
- Components for KASP
- KASP assay design
- Summary



Over 360,000 unique KASP assays developed in our laboratories



- >25,000 Validated Human assays.
- Multiple assays to key genetic diseases such as Cystic Fibrosis and Breast & Ovarian Cancer
- Extensive Panels for key field crops (~13,600 assays)

Title : QTL Mapping and Molecular Breeding for Developing Stress Resilient Maize for Sub-Saharan Africa (Oct 2014)

- KASP good at fingerprinting due to simplicity of data and low drop out rate
- CIMMYT routinely uses KASP, generating in excess of **a million data points** annually for different purposes
- The **KASP** average genotyping **error rate** in positive control DNA samples varied between 0.7 to 1.6%, which is **lower than** that observed in **the Illumina GoldenGate** (2.0-2.4 %)
- Because of the continuous change in SNP genotyping technology and high upfront costs, large-scale projects could benefit from exploiting genotyping service providers and avoid significant investment in setting up in-house SNP genotyping platforms, such as the KASP assay

For ALSPAC, the entire cohort (10,145 participants, including 38 carriers of the rare A allele) was genotyped using KASP with a **genotyping accuracy of 100%**

A rare variant in *APOC3* is associated with plasma triglyceride and VLDL levels in Europeans

One-step, codominant detection of imidazolinone resistance mutations in weedy rice (*Oryza sativa* L.)

KASP has **low cost**, **high throughput**, and **high specificity** and sensitivity as has been demonstrated in massive SNP genotyping studies in biomedical research [23], genome-wide SNP platforms for rice genotyping [24] and molecular markers-assisted wheat breeding

Thus, the method (**KASP**) validated here for timely and accurate detection of RWR is **a valuable and cost-effective tool** for decision making in Clearfield rice management and regional surveillance of RWR in the framework of a sustainable use of this production system.


Molecular Marker Applications in Oat (*Avena Sativa* L.) Breeding and Germplasm Diagnostics

GBS will likely provide a good source for future **KASP™** marker design

Therefore, we **recommend KASP™** for rapid cultivar identification and GBS for more conclusive identification and/or for purity assessment"


The present work can also be considered as a resource for breeders. Firstly, **KASP™** and GBS have been developed for rapid cultivar identification and purity assessments

- Kompetitive Allele Specific PCR
- Endpoint PCR-based genotyping chemistry
- Suitable for single nucleotide polymorphisms (SNPs) and Insertions / Deletions (InDels).
- Fluorescent signal is generated during the PCR - read at the end of the thermal cycle.
- Signal plotted on a cluster plot.




Taqman[®] vs KASP

Taqman[®] components:



KASP Components



Taqman[®]

KASP

Detect single nucleotide polymorphisms (SNPs)

Fluor-Quencher Probe

KASP Master MIX

1 Forward Primer

2 Forward Primer

2 dual labelled probes

No probes

complex regions/Gene expression

complex regions

large insertion/deletion




- Introduction to KASP
- Components for KASP
- KASP assay design
- Summary





KASP: components of the reaction

- The three key components for the KASP reaction are discussed in detail in the first video in the KASP series:




- The primers are contained within KASP Assay mix, and run with KASP Master mix
- Please note: KASP Assay mix works exclusively with KASP Master mix, and cannot be used with any other mix.

Assay components:



How does KASP work?


- Now that you have an understanding of the components in a KASP reaction, it is time to delve deeper into how KASP chemistry actually works
- The second video in the KASP series explains the details of this:



The video thumbnail shows the components of a KASP reaction: a test tube labeled 'KASP Assay mix', a brown bottle labeled 'KASP Master mix', and a circular icon labeled 'DNA Template' with 'AC' inside. Below, three microcentrifuge tubes are shown with labels 'T/M A:A', 'T/M C:C', and 'T/M A:C'. A video player interface at the bottom right shows a duration of '9:09'.

LGC's KASP Genotyping Chemistry Explained:
LGC Group video
6,472 views

Schematic drawing of the KASP method (Kompetitive Allele Specific PCR)



- Introduction to KASP
- Components for KASP
- **KASP assay design**
- Summary



KASP Assay design

- The main types of polymorphisms that can be genotyped using KASP are:

Different types of polymorphisms

○ SNP differentiation _____^A _____ or _____^G _____ ?

○ Single base insertion or deletion _____^A _____ or _____ - _____ ?

○ Small insertions or deletions

(2-26 bases) _____^{AA} _____ or _____^{AGCTGGTCA} _____ ?

○ Large insertions or deletions


(>26 bases) _____^{AA} _____ to _____^A - - - - -^A - - - - - _____ ?
n = 26 – 200,000 bases



KASP Assay design




Standard KASP assay

- Atypical KASP assay designed to a SNP
- The SNP is an A/G
 - One allele-specific primer is designed to the G base (detected with FAM)
 - One allele-specific primer is designed to the A base (detected with HEX)
 - A common reverse primer is also designed that will work with both allele-specific primers.




Note: despite the example the first letter inside the bracket is generally labelled with FAM and the one after the / symbol is labelled with Hex.

Please specify if the design requires specific labelling of each allele.

KEY	
	Allele-specific primer 1 (FAM)
	Allele-specific primer 2 (HEX)
	Common reverse primer



KASP Assay design




Typical results from a KASP SNP assay




Single base insertion / deletion

- KASP assays can be designed to detect single base insertions or deletions
 - One allele-specific primer is designed to incorporate the indel sequence (detected with FAM)
 - One allele-specific primer is designed to the standard sequence (detected with HEX)
 - A common reverse primer is also designed that will work with both allele-specific primers.



Insertion / deletion of 2-26 bases

- Standard KASP assays (3 primers) can be designed to detect insertions or deletions of up to 26 bases
 - One allele-specific primer is designed to incorporate the InDel sequence (detected with FAM)
 - One allele-specific primer is designed to the standard sequence (detected with HEX)
 - A common reverse primer is also designed that will work with both allele-specific primers.



KASP Assay design

Large insertion or deletion: overview

Insertion



Deletion





- Assays for large insertions / deletions require four primers rather than the typical three:
 - The first allele-specific primer is designed to incorporate part of the insertion sequence (detected with FAM), with a reverse primer located within the inserted sequence
 - The second allele-specific primer is designed across the region where the insertion would be (detected with HEX), with a reverse primer downstream of the insertion point.

KASP Assay design

Large insertion or **deletion**: assay design

- There is no upper size limit for the InDel aside from the requirement for the sequence to all be present on the same strand of DNA.



KEY

- Allele-specific primer 1 (FAM)
- Allele-specific primer 2 (HEX)
- Reverse primer 1
- Reverse primer 2
- Represents large insertion



KASP on qPCR instruments

- We know that KASP works well on most qPCR machines
- These include:



7900

ABI



Via7




7500



Agilent

Mx3000



BIO-RAD

CFX



Roche

LC480



- Introduction to KASP
- Components for KASP
- KASP assay design
- Summary



Summary

- KASP real-time mix is suitable for detect : agriculture animal and human targets
- assays are correctly designed : highly specific and sensitive detection
- analysis for detection : gene expression and copy number variation
- Use of KASP real-time for SNP genotyping and pathogen/species detection




- INTRODUCTION
- Sbeadex : DNA extraction kit
- KASP : method for genotyping
- **MASS ARRAY : Multi-gene Analysis**




MassARRAY: The Powerful Technology for Multi-gene Analysis


- INTRO
- MASS ARRAY DEBRIEF
- APPLICATIONS




Cost-Effective, Targeted Genetic Analysis Robust, Flexible, and High Throughput



SKYROCKET PUBLICATION OF MASSARRAY



- INTRO
- MASS ARRAY DEBRIEF
- APPLICATIONS





Chip prep module
96 well plate



384-well plate platform

Proven Versatility of the MassARRAY[®] System

- Proprietary mass spectrometry-based detection system with robust chemistries and advanced data analysis software for sensitive, accurate analyses



MULTI-PURPOSES



Genotyping

Somatic Mutation


Methylation


Copy Number Variation

Ultrasensitive Detection



MassARRAY™ Applications





The figure above illustrates an example of a 40-plex genotyping assay. The MassARRAY Assay Design software typically multiplexes up to 40 targets per reaction.

- Highly scalable throughput
 - 96 or 384-pad SpectroCHIP® Arrays*
 - From 1 to 10 arrays per day
 - 10s – 100s – 1000s samples per day


Marker x Sample flexibility

- Multiplex 5 - 40 targets per well
- Run 1 to 24 reactions per sample
- Screen 10s to 100s mutations per sample

FLEXIBLE DNA SAMPLING




- Just 5-10 ng DNA per well (reaction)
- FFPE, fresh and frozen tissue samples
- Liquid Biopsy
- Blood plasma and serum
- FTA cards
- Whole genome amplified DNA
- Buccal cells
- Ear punches
- Semen
- Hair pulls
- Micro-dissected cells



DNA FORMS OF NUCLEOBASES "A, C, G, T"





DIFFER IN MASS

dAMP = 313.2 Da

dTMP = 304.2 Da

dCMP = 289.2 Da


dGMP = 329.2 Da

MALDI-TOF


“Mass as Signature”



Matrix Assisted Laser Desorption/Ionization (MALDI)



Processes




- Molecular mass is a unique intrinsic property.
- The mass of a nucleic acid string is determined by its length and composition (ATCG).

GENOTYPING


- Example of an [A/G] SNP

PCR Reaction (Amplification) → **SAP Reaction** → **Single Base Extension**




Process Workflow : 7-8 hours

Process can be run overnight and interrupted at any point, post-PCR, post-SAP or post-Extend




MassARRAY[®] Process




Amplification




PCR Product




iPLEX[®]
Extension
Reaction




MALDI-TOF Mass
Spectrometry
Analysis






- The Power of Multiplexing Multiple Markers in a Single Well

ADVANCE SOFTWARE PROCESSING BASED ON "MASS"




Data Analyzer

Software application suite for Genotyping
& Somatic Mutation Analysis



- The Power of Multiplexing Multiple Markers in a Single Well



Generating Report

The screenshot displays the TyperAnalyzer@seqcom035 software interface, which is used for analyzing sequencing data. The interface is divided into several main sections:

- Project Explorer:** Shows a hierarchical tree of data files, including HapMap assays and samples, and various recall files.
- Traffic Light:** A grid of colored circles (green, yellow, red) representing data points for different assays across samples. A yellow box labeled "CHIP LIST" is overlaid on this section.
- Assay Summary:** A detailed view for the selected assay (rs2278725), showing call statistics such as Conservative Calls (378), Moderate Calls (1), Aggressive Calls (0), and User Calls (0). A yellow box labeled "Assay LIST" is overlaid on this section.
- Mass Spectrum:** A plot showing Intensity versus Mass for the assay rs2278725. The x-axis ranges from 6500 to 7100, and the y-axis ranges from 0 to 14. A red vertical line indicates the expected peak at mass 6600.
- Assay Table:** A table listing various assays and their call descriptions, such as rs2274739 (T, A.Conservative) and rs2278725 (CT, A.Conservative).
- Call Table:** A table showing expected peaks, probe names, and sequences, such as rs2278725 (ATCCTATCC) and rs2345450 (ATCCTATCC).



Data Analysis

TypewriterAnalyzer SEQUENOM®

File View Options Tools Help

Project Explorer Traffic Light Details

Date Assay Customer

Customer/Project/Plate/Experiment/C

1 2 3 4 5 6
 A ◆ ● ● ● ● ●
 B ● ● ● ● ● ●
 C ● ● ● ● ● ●
 D ● ● ● ● ● ●

Chip Experiment
 L0017691_(1) 1

Chip Summary

Chip L0017691_(1)
 Conservative Calls 444
 Moderate Calls 56
 Aggressive Calls 4
 User Calls 0
 Calls 504
 No Calls 24
 Total Possible Calls 528
 Call Rate 95.5
 Negative Controls 0

Assay

- Cum_Chr1_1143447
- Cum_Chr1_17305528
- Cum_Chr1_28662279
- Cum_Chr1_28667990
- Cum_Chr1_28999369
- Cum_Chr2_19978357
- Cum_Chr2_8104400
- Cum_Chr3_18775
- Cum_Chr3_35005861
- Cum_Chr3_37762964
- Cum_Chr4_17847401

Intensity

Cum_Chr1_28662279

Mass

Well: A02 Sample: KYTCU-1

Assay	Call	Description
Cum_Chr1_11...	G	A.Conservative
Cum_Chr1_17...	T	A.Conservative
Cum_Chr1_28...	T	A.Conservative
Cum_Chr1_28...	A	A.Conservative
Cum_Chr1_28...	C	B.Moderate
Cum_Chr2_19...	A	A.Conservative
Cum_Chr2_81...	T	A.Conservative

Expected Peaks	Name	Sequence	Mass
Probe	UEP.Cum_Ch...	CACACTTAAGCTTTGTTCATC	6331.10
Analyte	C	CACACTTAAGCTTTGTTCATCC	6578.30
Analyte	G	CACACTTAAGCTTTGTTCATCG	6618.30
Analyte	T	CACACTTAAGCTTTGTTCATCT	6658.20

Automate... Assay Su... Chip Sum... Post Processing Clusters Call Cluster Plot Histogram D|Details

Done SCRL



Data Analysis

TypewriterAnalyzer

File View Options Tools Help

Project Explorer Traffic Light Details

Date Assay Customer

Customer/Project/Plate/Experiment

	1	2	3	4	5	6
A	●	●	●	●	●	●
B	●	●	●	●	●	●
C	●	●	●	●	●	●
D	●	●	●	●	●	●

Chip L0017691_(1) (1) Experiment 1

Chip Summary

Chip L0017691_(1)

Conservative Calls 444

Moderate Calls 56

Aggressive Calls 4

User Calls 0

Calls 504

No Calls 24

Total Possible Calls 528

Call Rate 95.5

Negative Controls 0

Negative Control Calls 0

Automated... Assay Sum... Chip Sum...

6449.174, -0.726

Cum_Chr1_28662279

Intensity

Mass

Well: A03 Sample: KYTCU-2

Assay	Call	Description
Cum_Chr1_11...	G	A.Conservative
Cum_Chr1_17...	T	A.Conservative
Cum_Chr1_28...	GT	A.Conservative
Cum_Chr1_28...	A	A.Conservative

Expected Peaks	Name	Sequence	Mass
Probe	UEP.Cum_Ch...	CACACTTAAGCTTTGTTTCATC	6331.10
Analyte	C	CACACTTAAGCTTTGTTTCATCC	6578.30
Analyte	G	CACACTTAAGCTTTGTTTCATCG	6618.30

Post Processing Clusters Call Cluster Plot Histogram Details


Ready

SEQUENOM

SCRL



Data Analysis



Data Analysis

B [Compatibility Mode] - Excel

FILE HOME INSERT PAGE LAYOUT FORMULAS DATA REVIEW VIEW

Clipboard Font Alignment Number Styles Cells Editing

MS Sans Serif 10

General

Clipboard Font Alignment Number Styles Cells Editing

E2

	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	U	V	W	X	Y
1	rs1228662	rs1228667	rs1228999	rs2219978	rs2281044	rs3218775	rs3235005	rs3237762	rs4217847	rs4220081	rs4257284	rs5217301	rs5224960	rs5257581	rs6212587	rs6225093	rs6230700	rs7219148	rs7243002	rs728138860	
2	T	A	C	A	T	G	G	G	G	A	T	G	C	C	G	C	A	A	T	A	
3	GT	A	C	A	A	G	G	G	G	A	CT	G	CT	T	A	T	A	G	G	A	
4	T	A	C	A	TA	G	G	G	G	A	CT	G	CT	TC	AG	CT	A	GA	GT	A	
5	T	A	C	A	TA	G	G	G	G	A	T	G	C	TC	AG	CT	A	GA	GT	A	
6	GT	A	C	A	TA	G	G	G	G	A	CT	G	CT	TC	AG	CT	A	GA	GT	A	
7	T	A	C	G	T	G	G	A	G	T	T	A	C	T	A	C	A	G	G	A	
8	T	A	C	A	T	G	G	G	G	A	T	G	C	C	G	C	A	A	T	A	
9	T	A	C	A	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
10	T	A	C	AG	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
11	T	A	C	AG	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
12	T	A	C	AG	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
13	T	A	C	AG	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
14	T	A	C	A	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
15	T	A	C	AG	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
16	T	A	C	A	T	G	G	GA	G	T	GA	C	TC	AG	C	A	GA	GT	A		
17	T	A	C	A	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
18	T	A	C	AG	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
19	T	A	C	AG	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
20	T	A	C	AG	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
21	T	A	C	A	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
22	T	A	C	A	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
23	T	A	C	A	T	G	G	GA	G	T	GA	C	TC	AG	C	A	GA	GT	A		
24	T	A	C	A	T	G	G	GA	G	TA	T	GA	C	TC	AG	C	A	GA	GT	A	
25																					
26																					
27																					
28																					


Genotypes

READY 100%



Assay Design

- Fully customized options



- Agena or a CSP can design a custom panel for you



- Patient-specific panels for ultrasensitive detection
- Disease-specific panels for tumor profiling
- Gene-specific panels for therapy selection
- Mutation-specific panels for validation



<https://www.agenacx.com/Home>

HOME

AgenaCx

Returning Users
If you are a returning user from the old AgenaCx Site, and this is your first visit to our new site, please [update your password here](#) before logging in.

EMAIL

PASSWORD

REMEMBER ME

LOG IN

[Register](#) | [Lost your password?](#)

[Back to AgenaCx Home](#)

LOGIN OR REGISTER

AgenaCx provides you with the ability to:

- Access Agena Bioscience's online resources
- View and download the latest AgenaCx data files
- Exchange data files securely
- Contact your local Agena Bioscience representative
- View tutorials and training videos
- Log a support request

Our goal is to ensure your success with AgenaCx.

Welcome to the new AgenaCx site!
[Click here to view a short video!](#)



- PRODUCT SUPPORT
- ONLINE TOOLS**
- MEDIA LIBRARY
- MY FILES
- SUPPORT

ONLINE TOOLS

ASSAY DESIGN SUITE V2.0




LAUNCH

Assay Design Suite is a comprehensive and powerful tool for designing genotyping, somatic mutation, and ultra sensitive assays. The user-friendly interface integrates the design steps of importing RS numbers or sequences, retrieving and formatting sequences, finding proximal SNPs, identifying optimal areas for PCR primers, designing, and verifying multiplexed assays into a single software package. Assay Design Suite also enables users to design control SNPs into the panel, specify high priority SNPs, and use multiple design iterations for increased multiplexing efficiency. The enhanced exporting functionality generates output files for easy oligo ordering.

ATTENTION! - EpiDesigner Access

We are experiencing difficulties with access to EpiDesigner from this page; and are working to repair the problem. Temporarily, you may access EpiDesigner at <http://www.epidesigner.com>. We apologize for any inconvenience.

EPIDESIGNER



LAUNCH

EpiDesigner is a tool for designing quantitative methylation assays for



The screenshot displays the Agena Bioscience web interface. On the left, a 'Create New Design' sidebar contains three buttons: 'New Genotyping Design' (highlighted with a red box), 'New Somatic Mutation Design', and 'New UltraSEEK Design'. The main content area shows the 'Genotyping' design configuration form with the following fields:

- Design Type: Genotyping
- Design Name: Untitled
- Version: [dropdown]
- Current input: None (View button)
- rs or FASTA: [Edit Text Input button]
- Files: [File Upload button]
- Redesign options: [dropdown]
- Preset: [dropdown]
- Organism: Human
- Database: Dec 2013(GRCh38/hg38)
- Chemistry: iPLEX
- Multiplex Level: 36
- [Advanced Settings button]

To the right of the form is a progress table with five steps:

1	Retrieve & Format Sequences	Not Started
2	Find Proximal SNPs	Not Started
3	Identify Optimal Primer Areas	Not Started
4	Design Assays	Not Started
5	Validate	Not Started

Below the table are 'Export All' and 'Begin Run' buttons. At the bottom, a 'Design Assays' section shows 'Validate' with a checkmark and 'Rejects Completed'. A 'Design History' button is located at the bottom center.




Distinct advantages Of MassARRAY for Nucleic Acid Analysis

1. **Don't use fluorescence**
 - Mass of the actual bioanalyte is detected - 4 decimal place accuracy
 - No non-specific background issues – background is a different mass
2. **System is quantitative**
 - Many biological phenomena need to be accurately quantified
 - Allele ratios, gene copy number, methylation
3. **Multiplexed assays**
 - Provides high throughput
 - Cost effectiveness
4. **Very sensitive and flexible**
 - High accuracy – published 99.7%
 - High genotyping call rates (+98%)
 - Numbers of samples and markers are easily scaled
 - Simple and flexible assay design with little optimization required
 - Comprehensive Genetic Analysis >> Somatic mutation, Rare mutation, SNP, Epigenetic



- INTRO
- MASS ARRAY DEBRIEF
- APPLICATIONS



Sensitivity Across the Spectrum

Standard iPLEX

LOD: 5 – 10%

Germline
Inherited mutations

iPLEX HS

LOD: 1 – 10%

Somatic
Standard biopsy

UltraSeek™





LOD: 0.1 – 2%

Somatic
Liquid biopsy



Application of MassARRAY[®] in Human Genetics

Clinical Research Solutions : Targeted Biomarker Panels

	Application	Panel	Content
 ONCOLOGY	Tumor Profiling	OncoCarta™ v.1 OncoFOCUS™ v.3	19 oncogenes with 238 hot spot mutations EGFR, KRAS, NRAS, BRAF and KIT with 300+ mutations
	Lung Adenocarcinoma	LungCarta™ LungFusion™	26-gene with 200+ mutations ROS1, RET, ALK translocations
	Melanoma	MelaCarta™	20-gene panel with 70 mutations
	Colon	ColoCarta™	6-gene panel with 32 mutations
	Gynecological	GyneCarta™	13-gene panel with 168 mutations
	ccfDNA/CTC	UltraSEEK™	12-gene panel with 26 mutations
 GENETICS	Carrier Screening	Cystic Fibrosis	CFTR with 23+ mutations
	Blood Group Typing	Hemoid™	16 blood group systems, 23 platelet and neutrophil antigens
 PHARMACO GENOMICS	Drug Metabolism	ADME PGx ADME CYP2D6 ADME CYP2C19 ADME CYP2C9/VKORC1	36-gene panel with 191 SNPs CYP2D6 panel with 35 SNPs CYP2C19 with 31 SNPs CYP2C9 (36 SNPs), VKORC1 (9 SNPs)
 BIOMARKER VALIDATION		iPLEX® Pro Sample ID	44 SNP panel with gender ID, DNA QC
		iPLEX® Pro Exome ID	44-SNPs spanning intron-exon junctions with gender ID, DNA QC



Application	Panel	Content
Colon	IPLEX HS COLON	KRAS, NRAS, BRAF, EGFR, and PIK3CA with 86 mutations
Lung	OncoFOCUS™ v.3	EGFR, KRAS, BRAF, NRAS and KIT with 200+ mutations
	LungFUSION™	ALK, ROS1 and RET translocations
	IPLEX HS LUNG™	5-genesEGFR with 70 mutations
Melanoma	MelaCarta™	20-gene panel with 70 mutations
Gynecological	GyneCarta™	13-gene panel with 168 mutations
ccfDNA/CTC	UltraSEEK™	12-gene panel with 26 mutations



- Molecular Sequence Detection & Typing
 - Ability to detect and type microbes at the molecular level
 - Epidemiology & transmission studies
 - Bacterial/viral strain and sequence typing
 - Viral detection and typing
 - Outbreak monitoring and tracking
 - Ability to detect, type and quantify pathogens with the one system
 - High sensitivity & specificity
 - Multiplexing – create panels






Illustration: Chris Budd/Science. Reprinted with permission from Science Vol. 312, page 390 (21 April 2004) © 2004 by AAAS

- Human papillomavirus



- More than 120 different HPV types known
- DNA virus
- High risk types known to cause cancer

SCIENTIFIC REPORTS

OPEN **Rapid Sputum Multiplex Detection of the *M. tuberculosis* Complex (MTBC) and Resistance Mutations for Eight Antibiotics by Nucleotide MALDI-TOF MS**

Received: 22 July 2016
Accepted: 21 December 2016
Published: 30 January 2017

Kang-Yi Su^{1,2}, Bo-Shiun Yan³, Hao-Chieh Chiu^{1,2}, Chong-Jen Yu⁴, So-Yi Chang³, Ruwen Jou⁵, Jia-Long Liu², Po-Ren Hsueh^{2,4,*} & Sung-Liang Yu^{1,2,6,7,8,*}

The increasing incidence of multidrug-resistant (MDR) and extensively drug-resistant (XDR) *Mycobacterium tuberculosis* (MTB) adds further urgency for rapid and multiplex molecular testing to identify the MTB complex and drug susceptibility directly from sputum for disease control. A nucleotide matrix-assisted-laser-desorption-ionization time-of-flight mass spectrometry (MALDI-TOF MS)-based assay was developed to identify MTB (MTBID panel) and 45 chromosomal mutations for resistance to eight antibiotics (MTBDR panel). We conducted a 300 case trial from outpatients to evaluate this platform. An MTBID panel specifically identified MTB with as few as 10 chromosome DNA copies. The panel was 100% consistent with an acid-fast stain and culture for MTB, nontuberculous mycobacteria, and non-mycobacteria bacteria. The MTBDR panel was validated using 20 known MDR-MTB isolates. In a 64-case double-blind clinical isolates test, the sensitivity and specificity were 83% and 100%, respectively. In a 300-case raw sputum trial, the MTB identification sensitivity in smear-negative cases using MALDI-TOF MS was better than the COBAS assay (61.9% vs. 46.6%). Importantly, the failure rate of MALDI-TOF MS was better than COBAS (11.3% vs. 26.3%). To the best of our knowledge, the test described herein is the only multiplex test that predicts resistance for up to eight antibiotics with both sensitivity and flexibility.



Detection of HPV subtypes by mass spectrometry in FFPE tissue specimens: a reliable tool for routine diagnostics

Mark Kriegsmann,¹ Petra Wandernoth,² Katharina Lisenko,³ Rita Casadonte,⁴
Rémi Longuespée,⁴ Norbert Arens,² Jörg Kriegsmann^{2,5}

OPEN ACCESS Freely available online

PLOS ONE

Simultaneous Detection and Identification of Enteric Viruses by PCR-Mass Assay

Jingzi Piao^{1,2*}, Jun Jiang^{3*}, Bianli Xu^{4*}, Xiaohong Wang², Yanfang Guan⁵, Weili Wu³, Licheng Liu⁵,
Yuan Zhang³, Xueyong Huang⁴, Pengzhi Wang⁵, Jinyin Zhao³, Xiaoping Kang², Hua Jiang²,
Yuanyin Cao¹, Yuling Zheng², Yongqiang Jiang^{2*}, Yan Li^{6*}, Yinhui Yang^{2*}, Weijun Chen^{2,3,5*}

¹College of Plant Protection, Shenyang Agricultural University, Shenyang, China, ²State Key Laboratory of Pathogen and Biosecurity, Institute of Microbiology and Epidemiology, Academy of Military Medical Sciences, Beijing, China, ³Key Laboratory of Genome Sciences and Information, Beijing Institute of Genomics, Chinese Academy of Sciences, Beijing, China, ⁴Center for Disease Control and Prevention of Henan Province, Zhengzhou, China, ⁵Beijing Genomics Institute in Wuhan, Wuhan, China, ⁶Affiliated Hospital of Academy of Military Medical Sciences, Beijing, China

Syrmis et al. *BMC Infectious Diseases* 2014, **14**:307
<http://www.biomedcentral.com/1471-2334/14/307>



RESEARCH ARTICLE

Open Access

A comparison of two informative SNP-based strategies for typing *Pseudomonas aeruginosa* isolates from patients with cystic fibrosis

Melanie W Syrmis^{1,2}, Timothy J Kidd^{1,2}, Ralf J Moser³, Kay A Ramsay^{1,2}, Kristen M Gibson¹, Snehal Anuj¹, Scott C Bell^{1,4},
Claire E Wainwright^{1,5}, Keith Grimwood^{1,2}, Michael Nissen^{1,2,6}, Theo P Sloots^{1,2,6} and David M Whitley^{1,2*}



ELSEVIER

Contents lists available at ScienceDirect

Journal of Clinical Virology

journal homepage: www.elsevier.com/locate/jcv



Mass spectrometry-based comparative sequencing to detect ganciclovir resistance in the UL97 gene of human cytomegalovirus

Clara C. Posthuma^{a,*}, Martha T. van der Beek^{a,1}, Caroline S. van der Blij-de Brouwer^a,
Pim L.J. van der Heiden^b, Erik W.A. Marijt^b, Willy J.M. Spaan^a, Eric C.J. Claas^a,
Christa Nederstigt^a, Ann C.T.M. Vossen^a, Eric J. Snijder^a, Aloys C.M. Kroes^a

^aDepartment of Medical Microbiology, Center of Infectious Diseases, Leiden University Medical Center, P.O. Box 9600, 2300 RC Leiden, The Netherlands
^bDepartment of Hematology Leiden University Medical Center, P.O. Box 9600, 2300 RC Leiden, The Netherlands





Genotyping

Somatic Mutation

Methylation

Copy Number Variation

Ultrasensitive Detection



Multi-Purpose Utility of Circulating Plasma DNA Testing in Patients with Advanced Cancers

Geraldine Perkins¹, Timothy A. Yap^{1,2}, Lorna Pope¹, Amy M. Cassidy¹, Juliet P. Dukes¹, Ruth Riisnaes¹, Christophe Massard^{1,2}, Philippe A. Cassier², Susana Miranda¹, Jeremy Clark¹, Katie A. Denholm², Khin Thway¹, David Gonzalez De Castro¹, Gerhardt Attard^{1,2}, L. Rhoda Molife², Stan B. Kaye^{1,2}, Udai Banerji^{1,2}, Johann S. de Bono^{1,2*}

1 Division of Clinical Studies, The Institute of Cancer Research, Sutton, Surrey, United Kingdom, **2** Drug Development Unit, Royal Marsden NHS Foundation Trust, Sutton, Surrey, United Kingdom

Abstract



Tumor genomic instability and selective treatment pressures result in clonal disease evolution; molecular stratification for molecularly targeted drug administration requires repeated access to tumor DNA. We hypothesized that circulating plasma DNA (cpDNA) in advanced cancer patients is largely derived from tumor, has prognostic utility, and can be utilized for multiplex tumor mutation sequencing when repeat biopsy is not feasible. We utilized the Sequenom MassArray System and OncoCarta panel for somatic mutation profiling. Matched samples, acquired from the same patient but at different time points were evaluated; these comprised formalin-fixed paraffin-embedded (FFPE) archival tumor tissue (primary and/or metastatic) and cpDNA. The feasibility, sensitivity, and specificity of this high-throughput, multiplex mutation detection approach was tested utilizing specimens acquired from 105 patients with solid tumors referred for participation in Phase I trials of molecularly targeted drugs. The median cpDNA concentration was 17 ng/ml (range: 0.5–1600); this was 3-fold higher than in healthy volunteers. Moreover, higher cpDNA concentrations associated with worse overall survival; there was an overall survival (OS) hazard ratio of 2.4 (95% CI 1.4, 4.2) for each 10-fold increase in cpDNA concentration and in multivariate analyses, cpDNA concentration, albumin, and performance status remained independent predictors of OS. These data suggest that plasma DNA in these cancer patients is largely derived from tumor. We also observed high detection concordance for critical ‘hot-spot’ mutations (*KRAS*, *BRAF*, *PIK3CA*) in matched cpDNA and archival tumor tissue, and important differences between archival tumor and cpDNA. This multiplex sequencing assay can be utilized to detect somatic mutations from plasma in advanced cancer patients, when safe repeat tumor biopsy is not feasible and genomic analysis of archival tumor is deemed insufficient. Overall, circulating nucleic acid biomarker studies have clinically important multi-purpose utility in advanced cancer patients and further studies to pursue their incorporation into the standard of care are warranted.

Citation: Perkins G, Yap TA, Pope L, Cassidy AM, Dukes JP, et al. (2012) Multi-Purpose Utility of Circulating Plasma DNA Testing in Patients with Advanced Cancers. PLoS ONE 7(11): e47020. doi:10.1371/journal.pone.0047020


Editor: Jose Luis Perez-Gracia, University Clinic of Navarra, Spain

Received: May 15, 2012; **Accepted:** September 7, 2012; **Published:** November 7, 2012

MASSARRAY EARLY BREAST CANCER SCHEMA



CURATIVE SURGERY + ADJUVANT




2

Cancer DNA profiling

1

CTC PROGNOSIS

FOLLOW UP



3

cancer DNA detected



Application of MassARRAY® in Agricultural Genomics

- Parentage verification for plant breeding
- Crop strain validation
- Marker assisted breeding
- Candidate genetic marker evaluation
- Genetic Mapping and QTL analysis
- Phenotype selection
- SNP validation
- Etc.

Ting *et al.* *BMC Genomics* (2016) 17:289
DOI 10.1186/s12864-016-2607-4

BMC Genomics

RESEARCH ARTICLE

Open Access



Fine-mapping and cross-validation of QTLs linked to fatty acid composition in multiple independent interspecific crosses of oil palm

[Ting NC](#)^{1,2}, [Yaakub Z](#)¹, [Kamaruddin K](#)¹, [Mayes S](#)³, [Massawe F](#)², [Sambanthamurthi R](#)¹, [Jansen J](#)⁴, [Low LE](#)¹, [Ithnin M](#)¹, [Kushairi A](#)¹, [Arulandoo X](#)⁵, [Rosli R](#)¹, [Chan KL](#)¹, [Aminuddin N](#)¹, [Sritharan K](#)⁵, [Lim CC](#)⁵, [Nookiah R](#)¹, [Amiruddin MD](#)¹, [Singh R](#)⁶.

Candidate SNP markers (designated SNPE) flanking various genes associated with FA and oil biosynthesis were mined from the P5 genome build. The oil palm SNP assay design and genotyping were performed by a service provider, Agena Bioscience, Inc. (San Deigo, California) using the iPLEX[®] biochemistry on MassArray[®] system [23]. A custom two-multiplexed genotyping assay was designed and optimized for a panel of 40 SNPs using the Assay Design Suite 1.0 software (Agena Bioscience, Inc. San Deigo, California)



HHS Public Access

Author manuscript

Nature. Author manuscript; available in PMC 2014 October 26.


Published in final edited form as:

Nature. 2013 August 15; 500(7462): 340–344. doi:10.1038/nature12356.

The oil palm *Shell* gene controls oil yield and encodes a homologue of SEEDSTICK

Rajinder Singh¹, Eng-Ti Leslie Low¹, Leslie Cheng-Li Ooi¹, Meilina Ong-Abdullah¹, Ting Ngoot Chin¹, Jayanthi Nagappan¹, Rajanaidu Nookiah¹, Mohd Din Amiruddin¹, Rozana Rosli¹, Mohamad Arif Abdul Manaf¹, Kuang-Lim Chan¹, Mohd Amin Halim¹, Norazah Azizi¹, Nathan Lakey², Steven W Smith², Muhammad A Budiman², Michael Hogan², Blaire Bacher², Andrew Van Brunt², Chunyan Wang², Jared M Ordway², Ravigadevi Sambanthamurthi^{1,4}, and Robert A Martienssen^{3,4}

These 80 SNP markers (designated as SNPE) were genotyped in the T128 selfed population using the Sequenom MassArray® iPlex platform.



Imai, I., McCouch, S. R., and McClung, A. M. (2011) Plant and Animal Genomics, San Diego, California.


- Collaborative project among IRRI, JIRCAS, U. of Calcutta, U. of Aberdeen, Arkansas, Duke
- 18,760 SNPs discovered on Affymetrix 44K rice chip
- 84 SNPs for higher resolutions of key traits for MassARRAY panel
- Currently designing MassARRAY panels arranged by chromosome

8 SNP Marker Analysis for Quality and Agronomic Traits (rice)

Loci	Functions	Mutation and traits
<i>sd-1</i>	gibberellin 20-oxidase	C/T mutation reduces plant height and increases yield
<i>sd-del</i>		380-383bp deletion increases yield
<i>Pi-ta</i>	928-aa polypeptide	G: blast resistant; T: susceptible
<i>waxyIN1</i>	granule-bound starch synthase	Starch quality. A: low amylose; G: high amylose
<i>waxyEX6</i>		Starch quality. A: low amylose; T: high amylose
<i>alk3</i>	Starch synthase	Cooking and eating quality. alk3 'G'+alk4"GC": high gelatinization temperature and low alkali spreading
<i>alk4</i>		
<i>fgr</i>	BADH2	8bp deletion results in fragrant rice

Conclusion: MassARRAY is a powerful tool for SNP detection and confirmation in rice and can be applied to segregating populations in breeding programs for traits selection.

MassARRAY iPLEX Gold assays are capable to detect large or small deletions.



Haplotypes of the *BADH1* Gene in 92 Diverse Rice Varieties

Variety	S-1	S-2	S-3	S-4	S-5	S-6	S-8	S-9	S-11	S-14	S-15	S-16	S-17	S-18	S-19	SNP Haplotype	Protein Haplotype	Frequency
Jaya	G	C	G	T	T	A	A	C	C	T	T	T	T	A	T	SH1	PH1	38
ADT43	A	C	G	T	T	A	A	G	C	T	T	C	T	A	T	SH2	PH1	19
Basmati 370	G	C	A	A	C	T	G	G	T	T	C	C	C	C	T	SH3	PH2	17
Taraori Basmati	G	A	A	A	T	T	G	G	T	T	C	C	C	C	T	SH4	PH2	6
Kalanamak 3119	G	C	A	T	T	A	A	G	C	A	T	C	C	A	C	SH5	PH3	2
Taipai 309	G	C	A	A	T	T	G	G	T	T	C	C	C	C	T	SH6	PH2	1
Jyothi	G	C	G	T	C	A	G	C	C	T	T	T	T	A	T	SH7	PH1	1
Pusa 44	G	C	G	T	T	A	A	G	C	T	T	T	T	A	T	SH8	PH1	1
SKR 126	G	C	G	T	T	A	A	G	C	T	T	C	T	A	T	SH9	PH1	1
CSR 10	G	C	G	T	T	T	G	G	T	T	T	C	T	A	T	SH10	PH4	1
IR 64	G	C	G	T	T	A	A	C	C	T	T	T	C	A	T	SH11	PH1	1
Pusa 1266	G	C	A	A	C	T	G	C	T	T	C	C	C	C	T	SH12	PH2	1
Kasturi	G	C	A	T	C	T	G	G	T	T	C	C	C	C	T	SH13	PH2	1
Pusa 1121	A	C	G	T	T	A	A	C	C	T	T	C	T	A	T	SH14	PH1	1
Pant Dhan 4	G	C	G	T	T	A	A	C	C	T	T	T	T	C	T	SH15	PH5	1

Conclusion: There is no association between salt tolerance and the *BADH1* haplotypes

A high-throughput assay for rapid and simultaneous analysis of perfect markers for important quality and agronomic traits in rice using multiplexed MALDI-TOF mass spectrometry

Ardashir K. Masouleh¹, Daniel L. E. Waters¹, Russel F. Reinke² and Robert J. Henry^{1,*}

¹Centre for Plant Conservation Genetics, Southern Cross University, Lismore, NSW 2480, Australia

²Yanco Agricultural Institute, Yanco, NSW 2703, Australia

Received 19 December 2008;

revised 3 February 2009;

accepted 4 February 2009

*Correspondence (fax: +61 266222080;


e-mail: robert.henry@scu.edu.au)

Summary

The application of single nucleotide polymorphisms (SNPs) in plant breeding involves the analysis of a large number of samples, and therefore requires rapid, inexpensive and highly automated multiplex methods to genotype the sequence variants. We have optimized a high-throughput multiplexed SNP assay for eight polymorphisms which explain two agronomic and three grain quality traits in rice. Gene fragments coding for the agronomic traits plant height (semi-dwarf, *sd-1*) and blast disease resistance (*Pi-ta*) and the quality traits amylose content (*waxy*), gelatinization temperature (*alk*) and fragrance (*fgr*) were amplified in a multiplex polymerase chain reaction. A single base extension reaction carried out at the polymorphism responsible for each of these phenotypes within these genes generated extension products which were quantified by a matrix-assisted laser desorption ionization-time of flight system. The assay detects both SNPs and indels and is co-dominant, simultaneously detecting both homozygous and heterozygous samples in a multiplex system. This assay analyses eight functional polymorphisms in one 5 µL reaction, demonstrating the high-throughput and cost-effective capability of this system. At this conservative level of multiplexing, 3072 assays can be performed in a single 384-well microtitre plate, allowing the rapid production of valuable information for selection in rice breeding.

Keywords: matrix-assisted laser desorption ionization-time of flight (MALDI-TOF), multiplex perfect markers, rice, single nucleotide polymorphism.

Plant genotyping II: SNP technology



Description

This book aims to describe some of the important recent developments in plant genotyping. It is based upon a second workshop held recently to review progress in this area. Recent developments focus on high-throughput methods and generally target single nucleotide polymorphism (SNP) discovery and analysis. The topics covered include: SNP discovery in plants; SNPs and their use in maize; rare SNP discovery with endonucleases; sequence polymorphisms in the flanking regions of microsatellite markers; SNP discovery by ecotilling using capillary electrophoresis; genotyping by allele-specific PCR; the MassARRAY system for plant genomics; mutation screening; nanotechnology (the future of cost-effective plant genotyping); functionally associated molecular genetic markers for temperate pasture plant improvement; genotyping for rice eating qualities; towards universal loci for plant genotyping; DNA banks as a resource for SNP genotyping; DNA extraction from plant tissue; future prospects for plant genotyping.

Available in Print →

Book details

Editor(s)

Henry, R. J.

Author Affiliation

Centre for Plant Conservation Genetics
Southern Cross University Lismore, New
South Wales, Australia.

Year of Publication

2008

ISBN

9781845933821

DOI

10.1079/9781845933821.0000

Descriptor(s)

alleles

crop quality

enzyme activity

enzymes

genetic markers

genetic polymorphism

genetic variation

genotypes

Book Chapters

- Chapter: 1 (Page no: 1) **SNP discovery in plants.**
Author(s): Edward, K. J. Poole, R. L. Barker, G. L.
- Chapter: 2 (Page no: 30) **SNPs and their use in maize.**
Author(s): Rafalski, A. Tingey, S.
- Chapter: 3 (Page no: 44) **Rare SNP discovery with endonucleases.**
Author(s): Cross, M. J.
- Chapter: 4 (Page no: 68) **Sequence polymorphisms in the flanking regions of microsatellite markers.**
Author(s): Ablett, G. Henry, R. J.
- Chapter: 5 (Page no: 78) **Snp discovery by ecotilling using capillary electrophoresis.**
Author(s): Elliott, F. Cordeiro, G. Bundock, P. C. Henry, R. J.
- Chapter: 6 (Page no: 88) **Genotyping by allele-specific PCR.**
Author(s): Waters, D. L. E. Bundock, P. C. Henry, R. J.
- Chapter: 7 (Page no: 98) **The MassARRAY system for plant genomics.**
Author(s): Irwin, D.
- Chapter: 8 (Page no: 114) **Mutation screening.**
Author(s): Izquierdo, L.
- Chapter: 9 (Page no: 133) **Nanotechnology: the future of cost-effective plant genotyping.**
Author(s): Pattemore, J. A. Trau, M. Henry, R. J.
- Chapter: 10 (Page no: 154) **Functionally associated molecular genetic markers for temperate pasture plant improvement.**
Author(s): Forster, J. W. Cogan, N. O. I. Dobrowolski, M. P. Francki, M. G. Spangenberg, G. C. Smith, K. F.
- Chapter: 11 (Page no: 187) **Genotyping for rice eating qualities.**
Author(s): Bradbury, L. M. T. Waters, D. L. E. Henry, R. J.
- Chapter: 12 (Page no: 195) **Towards universal loci for plant genotyping.**
Author(s): Pacey-Miller, T.

[polymerase chain reaction](#)

[reviews](#)

[rice](#)

[single nucleotide polymorphism](#)

Subject Code(s)

[FF005 - Field Crops. \(New March 2000\)](#)

[FF020 - Plant Breeding and Genetics](#)

[ZZ360 - Molecular Biology and](#)

[Molecular Genetics. \(Discontinued](#)

[March 2000. Reinstated and Revised](#)

[June 2002\)](#)

Record Number

20083134963

Livestock Genotyping

- U.S. Department of Agriculture (USDA) Markers selected from HapMap project
- Genome wide SNPs with high MAF >0.41 selected for discrimination of inbred dairy herds
- 2 well panels for parentage verification with 100-121 SNPs per International Society for Animal Genetics (ISAG) recommendations
- Additional SNPs/wells for disease susceptibility and trait identification



Parental exclusion probabilities for SNP and microsatellite marker panels were similar, but genotyping was more sensitive for individual identification.

Rohrer, G. A. et al. Animal Genetics. 2007. 38 (3): 253–258

Comparison of 60 SNPs & 10 microsatellite markers

- Designed panel with 80 SNPs
 - 155 boars that represent four US purebred populations tested representative of US purebred Duroc, Hampshire, Landrace and Yorkshire populations
 - 60 SNPs w/ MAF >0.15 selected for the final panel of 60 markers
- Parentage Exclusion Probability
 - One parent: 0.9974 (all data)
 - Ranged from 0.9594 to 0.9963 within breeds.
- Sire Exclusion Probability
 - When dam's genotype is known: 0.99998 (all data)
 - Ranged from 0.99868 to 0.99997 within breeds





Relationships among *calpastatin* single nucleotide polymorphisms, *calpastatin* expression and tenderness in pork longissimus¹

A. K. Lindholm-Perry*, G. A. Rohrer*, J. W. Holl[†], S. D. Shackelford*, T. L. Wheeler*,
M. Koohmaraie[‡] and D. Nonneman*

*USDA, ARS, U.S. Meat Animal Research Center, Clay Center, NE 68933-0166, USA. [†]Smithfield Premium Genetics Group, Rose Hill, NC 28458, USA. [‡]IEH Laboratories and Consulting Group, Lake Forest Park, WA 98155, USA

Summary

Genome scans in the pig have identified a region on chromosome 2 (SSC2) associated with tenderness. *Calpastatin* is a likely positional candidate gene in this region because of its inhibitory role in the calpain system that is involved in postmortem tenderization. Novel single nucleotide polymorphisms (SNP) in *calpastatin* were identified and used to genotype a population ($n = 1042$) of Duroc–Landrace–Yorkshire swine for association with longissimus lumborum slice shear force (SSF) measured at days 7 and 14 postmortem. Three genetic markers residing in the *calpastatin* gene were significantly associated with SSF ($P < 0.0005$). Haplotypes constructed from markers in the *calpastatin* gene were significantly associated with SSF (F -ratio = 3.93; P -value = 0.002). The levels of normalized mRNA expression of *calpastatin* in the longissimus lumborum of 162 animals also were evaluated by real-time RT-PCR and were associated with the genotype of the most significant marker for SSF ($P < 0.02$). This evidence suggests that the causative variation alters expression of *calpastatin*, thus affecting tenderness. In summary, these data provide evidence of several significant, publicly available SNP markers associated with SSF that may be useful to the swine industry for marker assisted selection of animals that have more tender meat.

Keywords *calpastatin*, gene expression, meat quality, pig, single nucleotide polymorphism.

SNP markers identified by the association study should be predictive of pork tenderness in industry populations

Nonneman, D. et al. J. Anim. Sci. 2011: 2010-3556

- 194 SNPs in calpastatin locus on pig chr 2 were identified by re-sequencing
- SNP Genotyping on MassARRAY using 40 pigs
- 37 SNPs significantly associated with slice shear force in USMARC population
- 4 SNPs significantly associated with tenderness, juiciness, chewiness, and other flavor scores



in all 4 populations

Questions



PRECISION & MINIMIZE
HUMAN ERROR



“SPECTRO-CHIP”

A close-up photograph of a laboratory instrument. A hand wearing a blue nitrile glove is holding a small, square, dark-colored chip with a grid of tiny dots. The chip is being inserted into a slot in a metal tray. In the background, there are two larger, white, rectangular trays with a grid of circular wells, resting on a black surface. The overall scene is brightly lit, typical of a laboratory environment.

PATENTED "SPECTROCHIP"

MASSARRAY: KEY HIGHLIGHT




“MULTIPLEXING”

40 Target Specific DNA
Fragments per a Single
Reaction...

9 Laser shots per sample

Agena
BIOSCIENCE

COMPACT



Height	105.1 cm
Depth	67.0 cm
Width	47.2 cm