





ORBITRAP Mass Spectrometer


An Ultimate Qual and Quan Machine

Pongsagon Pothavorn
Scispec Co., Ltd.

Information Rich Data



Accurate Mass in Life Science



LC-MS solutions for all analytical challenges


- Best LC-MS Portfolio





Induced by ion packets moving inside the trap



- Ions trapped in an electrostatic field
- Central electrode kept on high voltage
- Outer electrode is split and able to pick up an image current induced by ion packets moving inside the trap



$$U(r, z) = \frac{k}{2} \times \left\{ z^2 - r^2 / 2 + R_m^2 \times \ln(r / R_m) \right\}$$


Ion Injection and Formation of Ion Rings

- An ion packet of a selected m/z enters the field
- Increasing voltage squeezes ions
- Voltage stabilises and ion trajectories are also stabilized
- Angular spreading forms a ROTATING RING




Fourier Transform-based

- The moving ion rings induce an image current on outer electrodes
- The frequency of harmonic oscillations is proportional to ions' m/z




Orbitrap and Nuclear Magnetic Resonance (NMR)

- Free Induction Decay (FID)



Time Domain $\rightarrow$ Fourier Transform $\rightarrow$ Spectrum (Frequency Domain)

Strategies for Analysis




Rapid and sensitive screening methods able to assign positive hits undoubtedly to particular organic compounds

Typical Mass Accuracy


Type of MS	Mass accuracy	Utility for
Quadrupole	0.1 μ	Identify
Traps	0.1 μ	Identify
TOF	0.0001 μ	Empirical formula/ composition
Sector	0.0001 μ	Empirical formula/ composition
FT-MS	0.0001 μ	Empirical formula/ composition

Isobaric Pesticides


Thiamethoxam: $[M+H]^+ = C_8H_{11}ClN_5O_3S$ (292.02656)



Parathion: $[M+H]^+ = C_{10}H_{15}NO_5PS$ (292.04031)




Isobaric Pesticides 3:1 Mix



Resolution – *Why Is It Important?*

- Enables accurate mass
- Increases confidence of identification
- Improves quantitative accuracy
- Gives access to qualitatively different information



Periodic Table of Elements

Average Mass

Average mass																																																													
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18																																												
1 H Hydrogen 1.00794																	2 He Helium 4.002602																																												
3 Li Lithium 6.941	4 Be Beryllium 9.012182	<table><tr><td>C Solid</td><td colspan="5">Metals</td><td colspan="5">Nonmetals</td></tr><tr><td>Hg Liquid</td><td>Alkali metals</td><td>Alkaline earth metals</td><td>Lanthanoids</td><td>Transition metals</td><td>Poor metals</td><td>Other nonmetals</td><td>Noble gases</td><td colspan="3" rowspan="2"></td></tr><tr><td>H Gas</td><td colspan="2"></td><td>Actinoids</td><td colspan="7" rowspan="2"></td></tr><tr><td>Rf Unknown</td><td colspan="10"></td></tr></table>										C Solid	Metals					Nonmetals					Hg Liquid	Alkali metals	Alkaline earth metals	Lanthanoids	Transition metals	Poor metals	Other nonmetals	Noble gases				H Gas			Actinoids								Rf Unknown											5 B Boron 10.811	6 C Carbon 12.0107	7 N Nitrogen 14.0067	8 O Oxygen 15.9994	9 F Fluorine 18.9984032	10 Ne Neon 20.1797
C Solid	Metals					Nonmetals																																																							
Hg Liquid	Alkali metals	Alkaline earth metals	Lanthanoids	Transition metals	Poor metals	Other nonmetals	Noble gases																																																						
H Gas			Actinoids																																																										
Rf Unknown																																																													
11 Na Sodium 22.98976928	12 Mg Magnesium 24.3050											13 Al Aluminium 26.9815386	14 Si Silicon 28.0855	15 P Phosphorus 30.973762	16 S Sulfur 32.065	17 Cl Chlorine 35.453	18 Ar Argon 39.948																																												
19 K Potassium 39.0983	20 Ca Calcium 40.078	21 Sc Scandium 44.955912	22 Ti Titanium 47.887	23 V Vanadium 50.9415	24 Cr Chromium 51.9961	25 Mn Manganese 54.938045	26 Fe Iron 55.845	27 Co Cobalt 58.933195	28 Ni Nickel 58.6934	29 Cu Copper 63.546	30 Zn Zinc 65.38	31 Ga Gallium 69.723	32 Ge Germanium 72.64	33 As Arsenic 74.92160	34 Se Selenium 78.96	35 Br Bromine 79.904	36 Kr Krypton 83.798																																												
37 Rb Rubidium 85.4678	38 Sr Strontium 87.62	39 Y Yttrium 88.90585	40 Zr Zirconium 91.224	41 Nb Niobium 92.90638	42 Mo Molybdenum 95.96	43 Tc Technetium (97.9072)	44 Ru Ruthenium 101.07	45 Rh Rhodium 102.90550	46 Pd Palladium 106.42	47 Ag Silver 107.8682	48 Cd Cadmium 112.411	49 In Indium 114.818	50 Sn Tin 118.710	51 Sb Antimony 121.760	52 Te Tellurium 127.60	53 I Iodine 126.90447	54 Xe Xenon 131.293																																												
55 Cs Caesium 132.9054519	56 Ba Barium 137.327	57–71										72 Hf Hafnium 178.49	73 Ta Tantalum 180.94788	74 W Tungsten 183.84	75 Re Rhenium 186.207	76 Os Osmium 190.23	77 Ir Iridium 192.217	78 Pt Platinum 195.084	79 Au Gold 196.966569	80 Hg Mercury 200.59	81 Tl Thallium 204.3833	82 Pb Lead 207.2	83 Bi Bismuth 208.98040	84 Po Polonium (209.9824)	85 At Astatine (209.9871)	86 Rn Radon (222.0178)																																			
87 Fr Francium (223)	88 Ra Radium (226)	89–103										104 Rf Rutherfordium (261)	105 Db Dubnium (262)	106 Sg Seaborgium (266)	107 Bh Bohrium (264)	108 Hs Hassium (277)	109 Mt Meitnerium (268)	110 Ds Darmstadtium (271)	111 Rg Roentgenium (272)	112 Uub Ununbium (285)	113 Uut Ununtrium (284)	114 Uuq Ununquadium (289)	115 Uup Ununpentium (288)	116 Uuh Ununhexium (292)	117 Uus Ununseptium	118 Uuo Ununoctium (294)																																			

For elements with no stable isotopes, the mass number of the isotope with the longest half-life is in parentheses.

Design and Interface Copyright © 1997 Michael Dayah (michael@dayah.com). <http://www.ptable.com/>

57 La Lanthanum 138.90547	58 Ce Cerium 140.116	59 Pr Praseodymium 140.90765	60 Nd Neodymium 144.242	61 Pm Promethium (145)	62 Sm Samarium 150.36	63 Eu Europium 151.964	64 Gd Gadolinium 157.25	65 Tb Terbium 158.92535	66 Dy Dysprosium 162.500	67 Ho Holmium 164.93032	68 Er Erbium 167.259	69 Tm Thulium 168.93421	70 Yb Ytterbium 173.054	71 Lu Lutetium 174.9668
89 Ac Actinium (227)	90 Th Thorium 232.03806	91 Pa Protactinium 231.03588	92 U Uranium 238.02891	93 Np Neptunium (237)	94 Pu Plutonium (244)	95 Am Americium (243)	96 Cm Curium (247)	97 Bk Berkelium (247)	98 Cf Californium (251)	99 Es Einsteinium (252)	100 Fm Fermium (257)	101 Md Mendelevium (258)	102 No Nobelium (259)	103 Lr Lawrencium (262)

How's About Mass Accuracy

- Average Mass = summing the [average atomic masses](#) of the constituent elements, H_2O ; $1.00794 + 1.00794 + 15.9994 = 18.01528$.
- Exact Mass = summing the masses of the individual isotopes of the molecule, H_2O ; $1.0078 + 1.0078 + 15.9994 = 18.0106$.

The Others Stories;

- Isotopomer (Isotopic Isomer) = same type of isotope but difference in position, $\text{CH}_3\text{CHDCH}_3$ vs $\text{CH}_3\text{CH}_2\text{CH}_2\text{D}$
- Isotopologues = difference in isotope in the molecules, H_2O HOD
- Monoisotopic = sum of masses in molecule. Using of most abundance or stable isotope.

Mass Accuracy – What for?

C = 12.0000

H = 1.0078


N = 14.0031

O = 15.9949

S = 31.9721

Mass measured	Tolerance [Da]	Suggestions	Calc Mass
32.0	+/- 0.2	O ₂ CH ₃ OH N ₂ H ₄ S	31.9898 32.0261 32.0374 31.9721
32.02	+/- 0.02	CH ₃ OH N ₂ H ₄	32.0261 32.0374
32.0257	+/- 0.002	CH₃OH	32.0261


Determine Fine Isotopic Pattern




Mass Accuracy across the Elution Profile

- 21 scans per elution peak
- External calibration

RI:72 - 1.96




Mass Accuracy [ppm]




Average Isotope Ratio Variation

pg on column	Standards	Apple juice	Baby food	Yogurt	Formula
A + 1					
Q-Exactive, Overall: 1.69 ± 2.30					
10	1.95 ± 2.26	3.17 ± 3.27	3.67 ± 3.33	3.21 ± 2.83	2.18 ± 1.69
100	2.61 ± 4.81	1.95 ± 1.98	1.91 ± 2.19	1.95 ± 1.87	2.10 ± 2.08
500	0.86 ± 0.96	1.07 ± 1.05	1.07 ± 1.18	1.26 ± 1.47	1.18 ± 1.36
2000	1.02 ± 1.79	0.75 ± 0.96	0.89 ± 1.34	0.74 ± 0.97	0.66 ± 0.89
MaXis, Overall: 5.01 ± 7.53					
10	9.20 ± 7.07	13.47 ± 9.06	15.30 ± 11.03	11.78 ± 7.62	11.49 ± 9.44
100	4.85 ± 6.66	7.78 ± 13.99	6.79 ± 7.02	6.94 ± 7.91	5.99 ± 6.25
500	3.05 ± 6.45	5.22 ± 9.58	3.30 ± 3.85	3.23 ± 3.79	3.33 ± 4.34
2000	1.77 ± 2.36	2.79 ± 6.28	2.13 ± 3.13	1.88 ± 2.56	2.03 ± 2.62
A + 2					
Q-Exactive, Overall: 1.59 ± 4.33					
10	5.31 ± 18.09	3.36 ± 5.42	4.38 ± 9.08	5.15 ± 6.56	6.44 ± 5.03
100	1.75 ± 3.01	1.93 ± 2.91	2.24 ± 4.60	1.70 ± 2.37	1.57 ± 1.86
500	1.03 ± 1.26	0.91 ± 0.62	0.86 ± 0.59	1.05 ± 0.81	1.22 ± 1.94
2000	0.81 ± 1.05	0.86 ± 1.20	0.73 ± 0.56	0.82 ± 0.57	0.74 ± 0.53
MaXis, Overall: 3.67 ± 6.47					
10	10.96 ± 9.71	12.89 ± 6.70	19.43 ± 38.22	11.21 ± 5.68	14.92 ± 7.62
100	3.55 ± 4.75	6.09 ± 6.85	6.73 ± 7.02	4.67 ± 4.46	5.22 ± 5.24
500	2.13 ± 3.14	4.02 ± 7.02	3.02 ± 3.17	3.01 ± 4.27	2.78 ± 3.38
2000	1.24 ± 2.06	2.23 ± 4.56	1.69 ± 2.36	1.68 ± 2.57	1.94 ± 3.21



Mass Accuracy




Resolving Power and Mass Accuracy



Res setting @ m/z 200	Transient length, ms	Max. scan speed, Hz
140,000	512	1.5
70,000	256	3
35,000	128	7
17,500	64	12



Long-term mass accuracy with external calibration



- Realistic conditions of an average lab
- Air cooling only!
- Experiment was carried out with temperature variations up to 3°C peak-to-peak, up to $1^{\circ}\text{C}/\text{hour}$

Advantage


- Easy method development for multi-residue analysis especially in complex matrices
- Easy troubleshooting with detection of all adducts, degradation and contaminants
- Higher detection specification
- Simultaneous Qual and Quan analysis

Comparison


	Orbitrap	QQQ	Q-TOF
Sensitivity	Green	Green	Yellow
Resolution	Green	Red	Yellow
Identification	Green	Yellow	Yellow
Unknowns	Green	Red	Yellow
Selectivity	Green	Yellow	Red
Quantitation	Green	Green	Yellow
Retrospective data mining	Green	Red	Yellow
Ease of troubleshooting	Green	Yellow	Yellow
Cost	Yellow	Green	Yellow

- High isolation power for higher discrimination
- High precision for accurate mass identification
- High resolution for more identification
- High mass stability for a long lasting mass calibration
- MSⁿ
- Library availability for easy interpretations

Orbitrap Analyzer - the 'Heart' of a Mass Spectrometer




Standard Orbitrap




High-field Orbitrap

Resolution VS m/z





Resolving Power




- 85% of the peptide parent ions (precursors) are below m/z 800
- Most of the chemical interferences are below m/z 600
- The highest resolution is needed below m/z 800 where Orbitrap technology has it and TOF technology doesn't!

Orbitrap VS QToF





Labelling Techniques



- Isobaric labeling techniques (iTRAQ, TMT) need high resolution at low masses
- Chemical interferences are common when using collision cells!


Intact Protein Analysis

- Complete charge state envelope of IgG 'Humira'
- Major glycosylation forms are baseline separated



Intact Protein Analysis

- Mass measurement accuracy
- Average error for 34 measurements 6.9 ppm
- Standard deviation 6.4 ppm




		ppm mass measurement errors				
Q	Exactive	G0+G0F	G0F+G0F	G0F+G1F	G0F+G2F	G1F+G2F
	1	-10.5	0.7	-10.5	-13.8	-18.0
	1	-3.2	-4.3	-6.9	3.2	N/A
	1	-11.6	-1.1	-8.8	-11.2	-12.0
	1	5.1	-5.0	-2.6	5.1	5.6
	2	-14.3	3.0	-6.9	-5.4	-5.9
	2	-8.6	-2.2	-12.2	-12.5	-12.9
	2	-14.3	-6.6	-12.3	-14.8	-10.1

Confirmation of protein primary structure

Sequence Confirmation of mAB



- ETD fragmentation of an intact IgG 'Humira'

- Resolution settings 240,000 for fragment detection
- Increased sequence coverage
- Localization of modifications (deamidation)



What do we gain by selected ion monitoring?


- Signal visibility is dependent, whether a signal is visible above the spectrum noise
- Spectrum noise is dependent on the ratio of compound within a certain ion population



Sensitivity gain 5 – 10 x with SIM mode


Full Scan Spectrum of Atenolol

AZ_1000ng_ml_100k_1e6_HypersilGoldPFP #246 RT: 3.46 AV: 1 SB: 1 3.25 NL: 1.36E6
T: FTMS + p ESI Full ms [140.00-1800.00]



Full Scan Spectrum of Pyridoxine


AZ_1000ng_ml_100k_1e6_HypersilGoldPFP #92 RT: 1.27 AV: 1 SB: 1 1.04 NL: 1.86E6
T: FTMS + p ESI Full ms [140.00-1800.00]



Alprazolam, Full Scan Experiment

Alprazolam


$Y = 6366.31 + 514.015 \cdot X$ $R^2 = 0.9967$ W: 1/X



Alprazolam SIM Experiment



Alprazolam

$Y = -3135.8 + 552.216 \cdot X$ $R^2 = 0.9982$ $W: 1/X$




- High isolation power for higher discrimination
- High precision for accurate mass identification
- High resolution for more identification
- High mass stability for a long lasting mass calibration
- MS^n
- Library availability for easy interpretations

Orbitrap VS QToF



Analysis of Protein Complexes

- Extending the mass range
- Protein assemblies up to 1 million Da





IgG antibody 150 kDa

HK97 bacteriophage capsomers 253 kDa

Yeast proteasome 730 kDa

***E. coli* GroEI 801 kDa**



Ligand Binding Stoichiometry



***E. coli* GroEI**
801 kDa


Data Dependent Decision Tree

- Decision tree–driven tandem mass spectrometry for shotgun proteomics




Product Dependent Trigger

- ZIC HILIC separation of a glycoprotein digest




Product Dependent Trigger

- HCD fragmentation spectrum of m/z 645.6194
- Oxonium ions observed among top 20 peaks



Extended Top-down Capability




Sequence coverage

Orbitrap Elite using ETD: 81%
LTQ Orbitrap XL using ETD: 43%

Product Dependent Trigger: HCD PD ETD

- ETD fragmentation triggered
 - Peptide sequence information
 - Glycosylation site localization




Elevator Speech



nri00

The orbitrap provides reproducible high resolution accurate mass with superior U-HPLC compatibility at resolution unattainable by QTOFs without compromising the sensitivity and dynamic range in MS or MS-MS data. With orbitrap, you will have fewer false positives, higher quality, better accuracy and more confidence in your quan/qual measurements.

From Discovery to Quantification - do it all with a Q Exactive



Range of Experiments

Quanfirmation = No Compromise!

Research
Low throughput
Discovery

Traditional
Proteomics,
Metabolomics,
Metabolism,
Biomarker Research

Development
Medium throughput
Verification

Translational Research,
Biopharma,
Metabolomics,
Drug Discovery,
Various Biomarker,
EFS Research

Routine
High throughput
Optimized assays

Clinical,
Pharma &
Biopharma Quantitation
EFS

Orbitraps

Exactives
& Ion Traps

Triples

Qualitative

Quanfirmation

Quantitative


Linearity and Precision

Milk Samples	Non Fat				Low Fat (2%)				Whole Fat			
Fortification Levels	50	100	250	500	50	100	250	500	50	100	250	500
Sulphamethazine	0.9964				0.9908				0.9970			
(RSD %)	2.4	7.2	4.9	2.4	6.6	14.5	5.2	5.6	8.9	1.2	5.1	n/a
Oxytetracycline	0.9906				0.9931				0.9909			
(RSD %)	2.7	11.4	11.4	5.5	11.6	11.8	5.9	4.2	9.2	12.9	4.9	2.6
Tetracycline	0.9923				0.9948				0.9903			
(RSD %)	6.2	6.4	5.4	3.7	7.3	4.8	5.9	4.5	4.1	9.5	6.5	5.5
Enrofloxacin	0.9969				0.9969				0.9973			
(RSD %)	9.2	7.9	2.3	0.8	11.1	1.4	1.4	3.1	6.9	8.3	3.0	n/a
Difloxacin	0.9958				0.9907				0.9968			
(RSD %)	12.2	4.3	6.0	2.4	10.8	4.6	2.7	5.5	2.6	6.1	5.1	n/a
Spiramycin	0.9920				0.9740				0.9951			
(RSD %)	11.1	11.8	8.4	4.1	10.9	4.0	10.0	9.4	13.3	6.5	5.2	0.2
Albendazole	0.9984				0.9967				0.9928			
(RSD %)	1.6	1.7	1.7	2.7	6.3	3.2	3.6	4.2	2.6	6.2	1.2	2.9
Phenylbutazone	0.9947				0.9922				0.9663			
(RSD %)	2.9	3.3	0.8	2.6	8.1	4.1	4.9	3.1	0.6	0.9	0.7	0.3
Salinomycine Na	0.9993				0.9966				0.9984			
(RSD %)	1.2	0.7	1.2	1.5	1.5	0.8	3.1	0.4	2.8	3.1	1.3	1.2

Stolker, A.A.M. et al; Anal. and Bioanal. Chem. 2010 accepted for publication

Drug identification using ToxID™ 2.1.0


- Fully automated analysis and reporting
- Drug identification based on
- Molecular weight
- MS2 spectra
- Chromatographic retention time
- Built-in library of about 300 drugs
- Library spectra acquired under real world conditions for robust and accurate ID
- The software uses proven NIST search engine
- Feature to easily create and expand library
- Excellent results review and reporting
- Summary report
- Data review report
- Excel spreadsheet



ToxID Summary Report


Your Company Summary Report

Raw File Name: C:\Documents and Settings\benedict.duretz\My Documents\My Data\Clinical Toxicology\Forensic Toxicology
 Config File Name: C:\Documents and Settings\benedict.duretz\My Documents\My Data\Clinical Toxicology\Forensic Toxicology
 Sample Name: Laboratory:
 Acquisition Start Time: 21/11/2008 11:43:54
 Screening Conditions: Based on Full MS2 scan. m/z window (amu): 0.50, RT window (min): 1.00, MS2 Search (library): Tox_Library_LXQ, Use full MS scan to confirm.



Peak Number	Compound Name	Code	SI	RSI	m/z	Expected RT	Actual RT	Concentration ng/ml	Library Name
1	Furosemide	p	999	999	329.00	8.70	9.40	0.03	Tox_Library_LXQ
2	Meprobamate	p	814	814	219.10	7.20	7.48	259.21	Tox_Library_LXQ
3	Venlafaxine	p	831	835	278.20	9.40	10.27	1.79	Tox_Library_LXQ
4	Diazepam	p	842	843	285.10	10.40	10.76	1.46	Tox_Library_LXQ
5	Chlorpromazine-D3	i	814	814	322.20	17.30	17.49	0.12	Tox_Library_LXQ
6	Prasopam-D5	i	911	916	330.20	11.70	12.59	1.69	Tox_Library_LXQ
7	Haloperidol-D4	i	806	826	380.40	14.90	15.30	1.00	Tox_Library_LXQ

Peak Number	Compound Name	Code	SI	RSI	m/z	Expected RT	Actual RT	Concentration ng/ml	Library Name
1	Furosemide	p	999	999	329.00	8.70	9.40	0.03	Tox_Library_LXQ
2	Meprobamate	p	814	814	219.10	7.20	7.48	259.21	Tox_Library_LXQ
3	Venlafaxine	p	831	835	278.20	9.40	10.27	1.79	Tox_Library_LXQ
4	Diazepam	p	842	843	285.10	10.40	10.76	1.46	Tox_Library_LXQ
5	Chlorpromazine-D3	i	814	814	322.20	17.30	17.49	0.12	Tox_Library_LXQ
6	Prasopam-D5	i	911	916	330.20	11.70	12.59	1.69	Tox_Library_LXQ
7	Haloperidol-D4	i	806	826	380.40	14.90	15.30	1.00	Tox_Library_LXQ



What is Mass Frontier?

- Software for small molecule structural elucidation via mass spectral interpretation
 - Predict fragmentation given a compound structure
 - Annotate spectra with fragment structures
 - Store MSⁿ spectra along with structures, peak annotations, ID numbers, pathway information, etc
 - Match unknown spectra against library entries
 - And **MUCH** more...

Tag Line:


Mass Frontier helps you to go from SPECTRA



Who should get Mass Frontier?


- Anyone who is doing small molecule structural elucidation / confirmation *via* mass spectrometry
- Examples:
 - Metabolite Identification in Drug Metabolism
 - Impurity and Degrading analysis in QC/QA
 - Endogenous Metabolite Identification in Metabolomics
 - Forensic Analyses in Federal and State Agencies
 - Doping Control in Horse Racing
 - Chemistry/Biochemistry/Pharmacy Departments in Universities doing small molecule research
 - Service labs for synthetic chemists

General Unknown Screening using Mass Frontier



Sheldon et al. Determination of Ion Structures in Structurally Related Compounds Using Precursor Ion Fingerprinting. JASMS, 2009, 20, 370-376

MSⁿ Spectral Trees—the **ONLY** Route to Unambiguous Structural Elucidation!



Accurate mass information is powerful
– provides a potential formula

However MSⁿ information still
necessary to distinguish between
structural isomers

Trees can automatically be generated
by Data Dependant LC-MS/MS runs on
our instruments

Component Detection from Mass
Frontier can automatically deconvolute
MSⁿ spectral trees!

*This information collectively, uniquely defines the structure of
the molecule*




Mass Frontier: Toolbox for Structural Elucidation

Chromatogram Processor

Easy Structural Editor


Fragmentation Pathways

Component Detection



Fragmentation Prediction: Three Knowledge Bases


1. General fragmentation rules



2. Mass Frontier Fragmentation Library™

Total number of	Mass Frontier 6.0
Fragmentation Schemes	30.936
Individual Reactions	129.229
Chemical Structures	151.762
Decoded Mechanisms	120.029

3. User Libraries



Your Scientific Specialist

Fragmentation Library™ in 6.0 now covers >99% published literature

	Source	Volume	Year
1.	JASMS (Journal of the American Society for Mass Spectrometry)	1-17	1990-2006
2.	IJMSIP (International Journal of Mass Spectrometry and Ion Physics)	1-53	1968-1983
	IJMSIP (International Journal of Mass Spectrometry and Ion Processes)	54-175	1983-1998
	IJMS (International Journal of Mass Spectrometry)	176-255	1998-2006
3.	RCM (Rapid Communications in Mass Spectrometry)	1-20	1987-2006
4.	JMS (Journal of Mass Spectrometry)	30-41	1995-2006
5.	OMS (Organic Mass Spectrometry)	1-29	1968-1994
6.	JMSSJ (Journal of the Mass Spectrometry Society of Japan)	11-27 29-30 37-48 50-53	1964-1979 1981-1982 1989-2000 2002-2005
7.	MSR (Mass Spectrometry Reviews)	1-25	1981-2006
8.	EJMSBMER (European Journal of Mass Spectrometry in Biochemical, Medicine, and Environmental. Research)	1-2	1980-1982
9.	BMS (Biomedical Mass Spectrometry) BEMS (Biomedical and Environmental Mass Spectrometry) BMS (Biological Mass Spectrometry)	1-12 14-19 20-23	1974-1985 1987-1990 1991-1994
10.	JC (Journal of Chromatography)	181-536	1980-1991
11.	EJMS (European Journal of Mass Spectrometry)	4	1998



Predictive Fragmentation

HighChem Mass Frontier 6.0

File Edit View Tools Search Library Options Help Microsoft Office

Database Manager: 2

Fragments & Mechanisms: 1

1 2 3 4 5 6 7 8 9 10 11 12 13 14 Select possible fragments with m/z 316.1114

m/z 346.1220 m/z 346.1220 m/z 316.1114

m/z 316.1114

118.0525
119.0604
120.0808
121.0886
122.0600
132.0310
133.0396
134.0475
135.0552
136.0757
144.0240
146.0475
147.0553
148.0757
149.0709
150.0913
151.0992
166.0321
167.0389
168.0478
172.0553
180.0478
181.0556
184.0427
195.0223
197.0379
198.0583
199.0662
200.0740
208.0301
209.0379
221.0379
225.0692
227.0849
233.0379
242.0350

Database Manager: 2

Reaction Info

m/z 150.0675 m/z 133.0441 m/z 107.0491 m/z 79.0542 m/z 77.0386

Results for

Work Angel Kevin buluzakal HighChem ESI Neg 2008 [R] HighChem ESI Pos 20

ID	Ac Mol. Mass	Formula	Title
2186	208.1094	C ₁₂ H ₁₀ O ₂ ¹⁺	Seiji Tobita, Susumu Tajima, Yasuko Ishihara
2187	192.1145	C ₁₁ H ₁₀ O ₂ ¹⁺	Seiji Tobita, Susumu Tajima, Yasuko Ishihara
2188	192.1145	C ₁₂ H ₁₀ O ₂ ¹⁺	Seiji Tobita, Susumu Tajima, Yasuko Ishihara
2189	150.0675	C ₉ H ₁₀ O ₂ ¹⁺	Osamu Sekiguchi, Tsutomu Noguchi, Kazuo
2190	150.0675	C ₉ H ₁₀ O ₂ ¹⁺	Osamu Sekiguchi, Tsutomu Noguchi, Kazuo

1. General Rules
2. Literature Library
3. User Libraries

How Do I Annotate Spectral Trees? ...Automatically

The screenshot shows the HighChem Mass Frontier 6.0 interface. The main window displays a mass spectrum with peaks at m/z 147.0543, 177.0559, 178.0420, 201.9983, 248.1120, 276.3029, 289.2111, 343.2263, 452.1914, and 469.3213. A chemical structure tree is visible on the left, and a fragmentation tree is on the right. A context menu is open over the spectrum, showing options like 'Manual Annotation', 'Auto Annotation', 'Auto Annotation Layout', 'Scale All Structure Annotations', 'Delete', 'Force Tolerance', 'Filter', 'Zoom', and 'Show Complete Chromatogram'. The 'Auto Annotation' option is selected, and a sub-menu is open showing 'Generated Fragments', 'Fragmentation Tab', and 'Generate Formulas'. In the foreground, the 'Auto Fragment Annotation of Peaks' dialog box is open, showing the following options:

- Options
- Threshold: 5.00 % of highest peak
- ☒ Clear Existed Fragments
- ☒ Check Precursor
- ☒ Auto Annotation Layout
- ☒ Apply to All Nodes
- Buttons: Restore Defaults, OK, Cancel

Below the spectrum, a table lists the identified peaks:

ID	Ac Mol. Mass	Formula	Title
1	469.3213	<chem>C32H40N2O3</chem> $7H^+$	
2	485.3163	<chem>C32H40N2O2</chem> $7H^+$	
3	455.3057	<chem>C31H38N2O3</chem> $7H^+$	
4	485.3163	<chem>C32H40N2O2</chem> $7H^+$	
5	501.3112	<chem>C32H40N2O3</chem> $7H^+$	

Database Manager: Integrated Knowledge Management

All records of installed libraries are shown in Database Manager

All records are accessible without querying

Spectral and Fragmentation libraries are unified in Database Manager

Searches are universal, independent of data type (structures, m/z values, names, CAS number, biological activity, etc)


One Record: Spectral tree with corresponding fragmentation mechanisms & more!

The left screenshot displays a mass spectrum plot with peaks at m/z 121.22, 154.04, 202.09, 232.09, and 261.10. The right screenshot displays a fragmentation tree diagram for the precursor ion at m/z 306.0674, leading to fragments at m/z 335.0662, 261.0703, and 251.0703.


Both screenshots include a table at the bottom with search results:

Match	ID Num.	Mol. Mass	Formula	Title	
>26	1	306.0674	C ₁₄ H ₁₄ N ₂ O ₄ S	ACEDIASULFONE	
27	94.6	439	307.0836	C ₁₀ H ₁₇ N ₃ O ₆ S	L-Glutathione reduced
28	79.2	63	304.9104	C ₇ H ₅ ClNO	CLOQUINOL
29	63.7	345	152.0473	C ₆ H ₈ O ₃	
30	63.2	520	152.0473	C ₆ H ₈ O ₃	

Compound Discoverer




Flexible Workflow




- Use common, predefined workflows or create your own
- Integrate your own nodes
- New in version 2.0: Nodes for Unknown Detection, Identification, Statistics


Identifying Unknown



Predicted Composition



Conclusion



- The Orbitrap Mass Analyzer is a new type of mass analyzers with its own unique combination of analytical parameters
- Orbitraps are still evolving...
 - Higher speed
 - Higher resolving power and mass accuracy
 - Higher sensitivity
 - More routine applications
- Exciting new applications continue to emerge

2



Summary

- **High resolution** is a key characteristics of MS data enabling
 - Mass accuracy
 - Confident identification
 - Reliable quantitation
- **Data dependent acquisition** offers an elegant simplicity and has proven highly useful for discovery-driven proteomics
- **Mass spectrometry technology** enables comprehensive analysis of proteomics samples
 - Multiple fragmentation techniques
 - MSⁿ capability
- **Quan&Qual** experiments done on a single platform