


ORBITRAP Mass Spectrometer

An Ultimate Qual and Quan Machine

Pongsagon Pothavorn
Scispec Co., Ltd.

Information Rich Data




Accurate Mass in Life Science

FT ICR

Orbitrap MS

TOF MS



- Motion in pre-dominantly magnetic field
- Low energy injection
- MSⁿ possibilities
- Broad-band excitation

$$T \propto m/z$$

- Excitation by injection

- Ions are trapped
- Image current detection
- Fourier transform and data processing
- Significant kinetic energy during detection
- Very long mean free path (many km)

$$T \propto \sqrt{m/z}$$

- Motion in electrostatic fields (m/z-independent well)
- High energy injection
- High energy spread upon fragmentation

$$TOF \propto \sqrt{m/z}$$

- Detection by secondary electron multiplier
- Very high kinetic energy for detection (many kV)
- Significant mean free path (tens of m)



Resolution

Thermo Scientific LTQ FT Ultra

Hybrid Mass Spectrometer
Unprecedented Analytical Power

tions



- Concentration independent ppb mass accuracy
- Widest dynamic range
- Parallel detection
- MSⁿ
- ECD and IRMPD
- Ultra high resolution
- Intelligent Data Dependent acquisition

Unprecedented Analytical Power

The Thermo Scientific LTQ FT Ultra hybrid mass spectrometer delivers unprecedented analytical power for the most demanding applications.

The unmatched mass accuracy eliminates false positive identifications in bottom-up and middle-down proteomics and enables the unambiguous identification of unknown analytes with on-line LC-MSⁿ at any concentration.

Ultra high resolution is essential for the analysis of complex samples such as crude oil, Dissolved Organic Matter (DOM) or intact proteins. This is available with a single mouse click.



solarix MRMS

Contact

Home - Products - Mass Spectrometry and Separations - MRMS - solarix - Overview


- Overview
- Technical Details
- Applications
- Learn more

solarix MRMS - routine IFS for unambiguous results

The solarix MRMS (Magnetic Resonance Mass Spectrometry) instruments are available for different magnetic field strengths of 7T, 12T and 15T.

	7T solarix	12T solarix	15T solarix
maximum resolving power	>10M	>10M	>10M
Mass accuracy (internal)	600 ppb	300 ppb	250 ppb
ESI	yes	yes	yes
MALDI	optional	optional	optional
ETD	optional	optional	optional
detector	ParaCell XR or 2xR	ParaCell XR	ParaCell XR
magnetic field	7T	12T	15T
annual cryogen fill	yes	yes	yes
quench line required	yes	yes	yes

Key features of solarix series instruments



Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

APPLICATIONS

PRODUCTS

RESOURCES

EDUCATION & EVENTS

SERVICES & SUPPORT

ABOUT WATERS



[Home](#) > [Products](#) > [Instruments & Systems](#) > [Mass Spectrometry](#) > [Ion Mobility Mass Spectrometry](#) >


[Vion IMS QTof Ion Mobility Quadrupole Time-of-flight Mass Spectrometry](#)

Vion IMS QTof Ion Mobility Quadrupole Time-of-flight Mass Spectrometry

Beyond Resolution

Does high resolution, accurate mass data always lead to the correct answer first time? Avoid additional confirmatory experiments by utilizing the selectivity of ion mobility. Featuring new geometry and including XS Ion Optics and the QuanTof2 detection system, the Vion IMS QTof has the sensitivity and dynamic range to make ion mobility quantitative and routinely usable, enabling:

- easier data interpretation as ion mobility cleans up and simplifies every spectrum
- confident identification and quantification of analytes
- faster method development and higher sample throughput
- greater selectivity than conventional MS/MS with HDMS^E




LC-MS solutions for all analytical challenges

- Best LC-MS Portfolio



Induced by ion packets moving inside the trap



- Ions trapped in an electrostatic field
- Central electrode kept on high voltage
- Outer electrode is split and able to pick up an image current induced by ion packets moving inside the trap



$$U(r, z) = \frac{k}{2} \cdot \left\{ z^2 - r^2 / 2 + R_m^2 \cdot \ln(r / R_m) \right\}$$


Ion Injection and Formation of Ion Rings

- An ion packet of a selected m/z enters the field
- Increasing voltage squeezes ions
- Voltage stabilises and ion trajectories are also stabilized
- Angular spreading forms a ROTATING RING




Fourier Transform-based

- The moving ion rings induce an image current on outer electrodes
- The frequency of harmonic oscillations is proportional to ions' m/z




Orbitrap and Nuclear Magnetic Resonance (NMR)

- Free Induction Decay (FID)



Time Domain \rightarrow Fourier Transform \rightarrow Spectrum (Frequency Domain)

Strategies for Analysis



Rapid and sensitive screening methods able to assign positive hits undoubtedly to particular organic compounds

Typical Mass Accuracy

Type of MS	Mass accuracy	Utility for
Quadrupole	0.1 μ	Identify
Traps	0.1 μ	Identify
TOF	0.0001 μ	Empirical formula/ composition
Sector	0.0001 μ	Empirical formula/ composition
FT-MS	0.0001 μ	Empirical formula/ composition

Mass Accuracy

- The precision of which the mass is measured by the mass spectrometer.
- Typical way of reporting mass error in ppm (relative measure) or mDa (absolute measure)

Good


$$\text{error} = \left(\frac{786.6003 - 786.60124}{786.60124} \right) \times 10^6$$

786.6003 (-1.19 ppm)

786.70 (+124 ppm)

Not so Good

$$\text{error} = \left(\frac{786.7000 - 786.60124}{786.60124} \right) \times 10^6$$



$$\text{Mass error} = \left(\frac{\text{Measured} - \text{Exact Mass}}{\text{Exact Mass}} \right) \times 10^6$$

C = 12.0000

O = 15.9949


H = 1.0078

S = 31.9721

N = 14.0031

Mass Resolution



- Ability of a mass spectrometer to distinguish between ions of nearly equal m/z ratios (isobars).

$$R = \frac{m}{\Delta m}$$



A diagram showing a single red peak. A horizontal double-headed arrow below the peak is labeled Δm (FWHM). The Δm in the label is circled in green.

m - measured mass

Δm - peak width measured at 50% peak intensity (Full Width Half Maximum)




Commercial High Resolution MS Technology Race




Isobaric Pesticides


Thiamethoxam: $[M+H]^+ = C_8H_{11}ClN_5O_3S$ (292.02656)



Parathion: $[M+H]^+ = C_{10}H_{15}NO_5PS$ (292.04031)




Isobaric Pesticides 3:1 Mix



Resolution – *Why Is It Important?*

- Enables accurate mass
- Increases confidence of identification
- Improves quantitative accuracy
- Gives access to qualitatively different information



Periodic Table of Elements

Average Mass

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
1 H Hydrogen 1.00794																	2 He Helium 4.002602
2 Li Lithium 6.941	4 Be Beryllium 9.012182											5 B Boron 10.811	6 C Carbon 12.0107	7 N Nitrogen 14.0067	8 O Oxygen 15.9994	9 F Fluorine 18.9984032	10 Ne Neon 20.1797
3 Na Sodium 22.98976928	12 Mg Magnesium 24.3050											13 Al Aluminium 26.9815386	14 Si Silicon 28.0855	15 P Phosphorus 30.973762	16 S Sulfur 32.065	17 Cl Chlorine 35.453	18 Ar Argon 39.948
4 K Potassium 39.0983	20 Ca Calcium 40.078	21 Sc Scandium 44.955912	22 Ti Titanium 47.887	23 V Vanadium 50.9415	24 Cr Chromium 51.9961	25 Mn Manganese 54.938045	26 Fe Iron 55.845	27 Co Cobalt 58.933195	28 Ni Nickel 58.6934	29 Cu Copper 63.546	30 Zn Zinc 65.38	31 Ga Gallium 69.723	32 Ge Germanium 72.64	33 As Arsenic 74.92160	34 Se Selenium 78.96	35 Br Bromine 79.904	36 Kr Krypton 83.798
5 Rb Rubidium 85.4678	38 Sr Strontium 87.62	39 Y Yttrium 88.90585	40 Zr Zirconium 91.224	41 Nb Niobium 92.90638	42 Mo Molybdenum 95.96	43 Tc Technetium (97.9072)	44 Ru Ruthenium 101.07	45 Rh Rhodium 102.90550	46 Pd Palladium 106.42	47 Ag Silver 107.8682	48 Cd Cadmium 112.411	49 In Indium 114.818	50 Sn Tin 118.710	51 Sb Antimony 121.760	52 Te Tellurium 127.60	53 I Iodine 126.90447	54 Xe Xenon 131.293
6 Cs Caesium 132.9054519	56 Ba Barium 137.327	57-71	72 Hf Hafnium 178.49	73 Ta Tantalum 180.94788	74 W Tungsten 183.84	75 Re Rhenium 186.207	76 Os Osmium 190.23	77 Ir Iridium 192.217	78 Pt Platinum 195.084	79 Au Gold 196.966569	80 Hg Mercury 200.59	81 Tl Thallium 204.3833	82 Pb Lead 207.2	83 Bi Bismuth 208.98040	84 Po Polonium (208.9824)	85 At Astatine (209.9871)	86 Rn Radon (222.0176)
7 Fr Francium (223)	88 Ra Radium (226)	89-103	104 Rf Rutherfordium (261)	105 Db Dubnium (262)	106 Sg Seaborgium (266)	107 Bh Bohrium (264)	108 Hs Hassium (277)	109 Mt Meitnerium (268)	110 Ds Darmstadtium (271)	111 Rg Roentgenium (272)	112 Uub Ununbium (285)	113 Uut Ununtrium (284)	114 Uuq Ununquadium (289)	115 Uup Ununpentium (288)	116 Uuh Ununhexium (292)	117 Uus Ununseptium	118 Uuo Ununoctium (294)

- C** Solid
- Hg** Liquid
- H** Gas
- Rf** Unknown

Metals						Nonmetals	
Alkali metals	Alkaline earth metals	Lanthanoids	Transition metals	Poor metals	Other nonmetals	Noble gases	
		Actinoids					

For elements with no stable isotopes, the mass number of the isotope with the longest half-life is in parentheses.

Design and Interface Copyright © 1997 Michael Dayah (michael@dayah.com). <http://www.ptable.com/>

57 La Lanthanum 138.90547	58 Ce Cerium 140.116	59 Pr Praseodymium 140.90765	60 Nd Neodymium 144.242	61 Pm Promethium (145)	62 Sm Samarium 150.36	63 Eu Europium 151.964	64 Gd Gadolinium 157.25	65 Tb Terbium 158.92535	66 Dy Dysprosium 162.500	67 Ho Holmium 164.93032	68 Er Erbium 167.259	69 Tm Thulium 168.93421	70 Yb Ytterbium 173.054	71 Lu Lutetium 174.968
89 Ac Actinium (227)	90 Th Thorium 232.03806	91 Pa Protactinium 231.03588	92 U Uranium 238.02891	93 Np Neptunium (237)	94 Pu Plutonium (244)	95 Am Americium (243)	96 Cm Curium (247)	97 Bk Berkelium (247)	98 Cf Californium (251)	99 Es Einsteinium (252)	100 Fm Fermium (257)	101 Md Mendelevium (258)	102 No Nobelium (259)	103 Lr Lawrencium (260)

How's About Mass Accuracy

- Average Mass = summing the [average atomic masses](#) of the constituent elements, H₂O; $1.00794 + 1.00794 + 15.9994 = 18.01528$.
- Exact Mass = summing the masses of the individual isotopes of the molecule, H₂O; $1.0078 + 1.0078 + 15.9994 = 18.0106$.

The Others Stories;

- Isotopomer (Isotopic Isomer) = same type of isotope but difference in position, CH₃CHDCH₃ vs CH₃CH₂CH₂D
- Isotopologues = difference in isotope in the molecules, H₂O HOD
- Monoisotopic = sum of masses in molecule. Using of most abundance or stable isotope.

Mass Accuracy – What for?

C = 12.0000

H = 1.0078


N = 14.0031

O = 15.9949

S = 31.9721

Mass measured	Tolerance [Da]	Suggestions	Calc Mass
32.0	+/- 0.2	O ₂ CH ₃ OH N ₂ H ₄ S	31.9898 32.0261 32.0374 31.9721
32.02	+/- 0.02	CH ₃ OH N ₂ H ₄	32.0261 32.0374
32.0257	+/- 0.002	CH₃OH	32.0261


Determine Fine Isotopic Pattern




Mass Accuracy across the Elution Profile

- 21 scans per elution peak
- External calibration

RT: 7.2 - 1.96




Mass Accuracy [ppm]



Average Isotope Ratio Variation




pg on column	Standards	Apple juice	Baby food	Yogurt	Formula
A + 1					
Q-Exactive, Overall: 1.69 ± 2.30					
10	1.95 ± 2.26	3.17 ± 3.27	3.67 ± 3.33	3.21 ± 2.83	2.18 ± 1.69
100	2.61 ± 4.81	1.95 ± 1.98	1.91 ± 2.19	1.95 ± 1.87	2.10 ± 2.08
500	0.86 ± 0.96	1.07 ± 1.05	1.07 ± 1.18	1.26 ± 1.47	1.18 ± 1.36
2000	1.02 ± 1.79	0.75 ± 0.96	0.89 ± 1.34	0.74 ± 0.97	0.66 ± 0.89
MaXis, Overall: 5.01 ± 7.53					
10	9.20 ± 7.07	13.47 ± 9.06	15.30 ± 11.03	11.78 ± 7.62	11.49 ± 9.44
100	4.85 ± 6.66	7.78 ± 13.99	6.79 ± 7.02	6.94 ± 7.91	5.99 ± 6.25
500	3.05 ± 6.45	5.22 ± 9.58	3.30 ± 3.85	3.23 ± 3.79	3.33 ± 4.34
2000	1.77 ± 2.36	2.79 ± 6.28	2.13 ± 3.13	1.88 ± 2.56	2.03 ± 2.62
A + 2					
Q-Exactive, Overall: 1.59 ± 4.33					
10	5.31 ± 18.09	3.36 ± 5.42	4.38 ± 9.08	5.15 ± 6.56	6.44 ± 5.03
100	1.75 ± 3.01	1.93 ± 2.91	2.24 ± 4.60	1.70 ± 2.37	1.57 ± 1.86
500	1.03 ± 1.26	0.91 ± 0.62	0.86 ± 0.59	1.05 ± 0.81	1.22 ± 1.94
2000	0.81 ± 1.05	0.86 ± 1.20	0.73 ± 0.56	0.82 ± 0.57	0.74 ± 0.53
MaXis, Overall: 3.67 ± 6.47					
10	10.96 ± 9.71	12.89 ± 6.70	19.43 ± 38.22	11.21 ± 5.68	14.92 ± 7.62
100	3.55 ± 4.75	6.09 ± 6.85	6.73 ± 7.02	4.67 ± 4.46	5.22 ± 5.24
500	2.13 ± 3.14	4.02 ± 7.02	3.02 ± 3.17	3.01 ± 4.27	2.78 ± 3.38
2000	1.24 ± 2.06	2.23 ± 4.56	1.69 ± 2.36	1.68 ± 2.57	1.94 ± 3.21

Mass Accuracy





Resolving Power and Mass Accuracy

Res setting @ m/z 200	Transient length, ms	Max. scan speed, Hz
140,000	512	1.5
70,000	256	3
35,000	128	7
17,500	64	12



Long-term mass accuracy with external calibration



- Realistic conditions of an average lab
- Air cooling only!
- Experiment was carried out with temperature variations up to 3°C peak-to-peak, up to 1°C/hour

Advantage


- Easy method development for multi-residue analysis especially in complex matrices
- Easy troubleshooting with detection of all adducts, degradation and contaminants
- Higher detection specification
- Simultaneous Qual and Quan analysis

Comparison


	Orbitrap	QQQ	Q-TOF
Sensitivity	Green	Green	Yellow
Resolution	Green	Red	Yellow
Identification	Green	Yellow	Yellow
Unknowns	Green	Red	Yellow
Selectivity	Green	Yellow	Red
Quantitation	Green	Green	Yellow
Retrospective data mining	Green	Red	Yellow
Ease of troubleshooting	Green	Yellow	Yellow
Cost	Yellow	Green	Yellow

- High isolation power for higher discrimination
- High precision for accurate mass identification
- High resolution for more identification
- High mass stability for a long lasting mass calibration
- MSⁿ
- Library availability for easy interpretations

Orbitrap Analyzer - the 'Heart' of a Mass Spectrometer



Standard Orbitrap




High-field Orbitrap

Resolution VS m/z





Resolving Power

Peptide ID Distribution of Precursors from LC/MS of *E. Coli* digest




- **85% of the peptide parent ions (precursors) are below m/z 800**
- **Most of the chemical interferences are below m/z 600**
- **The highest resolution is needed below m/z 800 where Orbitrap technology has it and TOF technology doesn't!**

Orbitrap VS QToF



Labelling Techniques





**16,000 FWHM is
the minimum required!**


- **Isobaric labeling techniques (iTRAQ, TMT) need high resolution at low masses**
- **Chemical interferences are common when using collision cells!**

Intact Protein Analysis

- Complete charge state envelope of IgG 'Humira'
- Major glycosylation forms are baseline separated




- Relative intensity reproducibility within a few percent



Q Exactive	Relative abundances				
	G0+G0F	G0F+G0F	G0F+G1F	G0F+G2F	G1F+G2F
1	12.9	74.1	100.0	67.0	23.4
1	12.3	76.0	100.0	71.4	29.8
1	12.0	72.8	100.0	66.2	22.0
1	12.2	75.0	100.0	67.0	23.6
2	12.7	75.7	100.0	63.6	21.6
2	13.2	75.4	100.0	64.8	21.0
2	12.9	76.6	100.0	64.7	21.6

Intact Protein Analysis

- Mass measurement accuracy
- Average error for 34 measurements 6.9 ppm
- Standard deviation 6.4 ppm




	ppm mass measurement errors				
Q Exactive	G0+G0F	G0F+G0F	G0F+G1F	G0F+G2F	G1F+G2F
1	-10.5	0.7	-10.5	-13.8	-18.0
1	-3.2	-4.3	-6.9	3.2	N/A
1	-11.6	-1.1	-8.8	-11.2	-12.0
1	5.1	-5.0	-2.6	5.1	5.6
2	-14.3	3.0	-6.9	-5.4	-5.9
2	-8.6	-2.2	-12.2	-12.5	-12.9
2	-14.3	-6.6	-12.3	-14.8	-10.1


Confirmation of protein primary structure

Sequence Confirmation of mAB

- ETD fragmentation of an intact IgG 'Humira'





- Resolution settings 240,000 for fragment detection
- Increased sequence coverage
- Localization of modifications (deamidation)



What do we gain by selected ion monitoring?

- Signal visibility is dependent, whether a signal is visible above the spectrum noise
- Spectrum noise is dependent on the ratio of compound within a certain ion population




Sensitivity gain 5 – 10 x with SIM mode

◆ S/N (FMS) ■ S/N (SIM10)


Full Scan Spectrum of Atenolol

AZ_1000ng_ml_100k_1e6_HypersilGoldPFP #246 RT: 3.46 AV: 1 SB: 1 3.25 NL: 1.36E6
T: FTMS + p ESI Full ms [140.00-1800.00]



Full Scan Spectrum of Pyridoxine


AZ_100ng_ml_100k_1e6_HypersilGoldPFP #92 RT: 1.27 AV: 1 SB: 1 1.04 NL: 1.86E6
T: FTMS + p ESI Full ms [140.00-1800.00]



Alprazolam, Full Scan Experiment

Alprazolam


$$Y = 6366.31 + 514.015 * X \quad R^2 = 0.9967 \quad W: 1/X$$



Alprazolam SIM Experiment

Alprazolam

$$Y = -3135.8 + 552.216 * X \quad R^2 = 0.9982 \quad W: 1/X$$



A Switch is on from QqQ to Orbitrap



Q Exactive Focus quadrupole-Orbitrap MS



TSQ Endura triple quadrupole MS

TABLE 6. Accuracy (% recovery, NIST standard concentration ng/mL) using PRM.

	25OH D2	25OH D3	epi-25OH D3
Level 1	BQL (0.54)	108 (28.8)	103 (1.84)
Level 2	BQL (0.81)	110 (18.1)	90.7 (1.29)
Level 3	97 (13.3)	104 (19.8)	114 (1.18)
Level 4	BQL (0.55)	106 (29.4)	90.9 (26.4)

TABLE 7. Accuracy (% recovery) using SRM.

	25OH D2	25OH D3	epi-25OH D3
Level 1	BQL	93.4	BQL
Level 2	BQL	102	BQL
Level 3	92.7	105	BQL
Level 4	BQL	98	100

TABLE 1. Linearity range (ng/mL)

	25OH D2	25OH D3	epi-25OH D3		
PRM	1-100	1-100	1-100		
SRM	2-100	2-100	2-100		
	Amo	Butalbital	Pento	Pheno	Seco
PRM	5-2000	5-2000	5-2000	25-2000	5-2000
SRM	10-2000	5-2000	10-2000	25-2000	5-2000

TABLE 2. Inter-assay precision (% RSD, concentration ng/mL) using PRM (n=15). BQL, below quantitation limit.

	25OH D2	25OH D3	epi-25OH D3
QC0	BQL	6.4 (8.4)	BQL
QC1	5.3 (6.0)	4.1 (14.4)	4.4 (6.6)
QC2	3.1 (15.0)	3.2 (23.4)	2.9 (15.6)
QC3	4.2 (50.0)	3.5 (58.4)	2.5 (50.6)

TABLE 3. Inter-assay precision (% RSD, concentration is the same as above) using SRM (n=15).

	25OH D2	25OH D3	epi-25OH D3
QC0	BQL	7.6	BQL
QC1	7.6	8.6	6.9
QC2	6.4	4.8	9.4
QC3	5.2	3.4	6.9

Biocrates – 188 Metabolites

Key Features

Quantification of up to 188 Metabolites from Key Metabolite Classes

Metabolite Coverage

- Amino acids (21)
- Biogenic amines (21)
- Hexose (1)
- Acylcarnitines (40)
- Lysophosphatidylcholines (14)
- Phosphatidylcholines (76)
- Sphingolipids (15)

Kit Components

Reagents & Consumables

- Patented 96-well filter plate
- Internal & calibration standards
- Quality controls
- Testmix

Methods & Protocols

- Sample preparation protocol
- Instrument-specific acquisition & quantification methods
- System suitability test

Workflow Manager Met/DQ™

- Process guidance
- Automated technical validation
- Basic statistics (optional)

Instrument Platforms*

- Waters Xevo TQ-XS UHPLC
- Waters Xevo TQ-S UHPLC
- Waters Xevo TQ-S micro UHPLC
- SCIEX 6500+ (U)HPLC
- SCIEX 6500 (U)HPLC
- SCIEX 5500 (U)HPLC
- SCIEX 4500 (U)HPLC
- SCIEX 4000 HPLC
- Thermo TSQ Vantage (U)HPLC



Biocrates – ONLY Orbitrap 408 Metabolites



The AbsoluteIDQ[®] p400 HR Kit is a complete solution for broad lipid and metabolic profiling on high resolution, accurate mass (HRAM) Q Exactive[™] mass spectrometers. It provides quantification of a large range of analytes, with high inter-laboratory, inter-instrument and longitudinal reproducibility.

The kit covers up to 408 metabolites from 11 metabolite classes, which are known to be relevant in a multitude of pathophysiological processes:

Small Molecules

Amino Acids

Biogenic Amines

Hexoses

Polar Lipids

Phosphatidylcholines

Lysophosphatidylcholines

Sphingomyelins

Ceramides

Neutral Lipids

Acylcarnitines



Diglycerides

Triglycerides

Cholesteryl Esters


- High isolation power for higher discrimination
- High precision for accurate mass identification
- High resolution for more identification
- High mass stability for a long lasting mass calibration
- MSⁿ
- Library availability for easy interpretations

Orbitrap VS QToF



Analysis of Protein Complexes

- Extending the mass range
- Protein assemblies up to 1 million Da





IgG antibody 150 kDa

HK97 bacteriophage capsomers 253 kDa

Yeast proteasome 730 kDa

***E. coli* GroEl 801 kDa**



Ligand Binding Stoichiometry



***E. coli* GroEI**
801 kDa


Data Dependent Decision Tree

- Decision tree–driven tandem mass spectrometry for shotgun proteomics




Product Dependent Trigger

- ZIC HILIC separation of a glycoprotein digest






Product Dependent Trigger

- HCD fragmentation spectrum of m/z 645.6194
- Oxonium ions observed among top 20 peaks



Extended Top-down Capability




Sequence coverage

Orbitrap Elite using ETD: 81%

LTQ Orbitrap XL using ETD: 43%

Product Dependent Trigger: HCD PD ETD

- ETD fragmentation triggered
 - Peptide sequence information
 - Glycosylation site localization



Structural Elucidation

Discovery and Characterization of Components




Compound Discoverer 2.1 SP1

Configuring catalog for MEF


© Copyright 2014-2017 Thermo Fisher Scientific Inc.
All rights reserved. This program is protected by copyright
law and international treaties as described in Help About.

thermo
scientific


MS and MS² Spectrum



Extracted Ion Chromatogram and Isotopic Pattern



MS² Spectrum Interpretation




Elevator Speech



nri00

The orbitrap provides reproducible high resolution accurate mass with superior U-HPLC compatibility at resolution unattainable by QTOFs without compromising the sensitivity and dynamic range in MS or MS-MS data. With orbitrap, you will have fewer false positives, higher quality, better accuracy and more confidence in your quan/qual measurements.

From Discovery to Quantification - do it all with a Q Exactive



Range of Experiments

Quanfirmation = No Compromise!

Research
Low throughput
Discovery

Development
Medium throughput
Verification

Routine
High throughput
Optimized assays

Traditional
Proteomics,
Metabolomics,
Metabolism,
Biomarker Research

Translational
Research,
Biopharma,
Metabolomics,
Drug Discovery,
Various Biomarker,
EFS Research

Clinical,
Pharma &
Biopharma
Quantitation
EFS

Orbitraps

Exactives
& Ion Traps

Triples

Qualitative

Quanfirmation

Quantitative


Linearity and Precision

Milk Samples	Non Fat				Low Fat (2%)				Whole Fat			
Fortification Levels	50	100	250	500	50	100	250	500	50	100	250	500
Sulphamethazine	0.9964				0.9908				0.9970			
(RSD %)	2.4	7.2	4.9	2.4	6.6	14.5	5.2	5.6	8.9	1.2	5.1	n/a
Oxytetracycline	0.9906				0.9931				0.9909			
(RSD %)	2.7	11.4	11.4	5.5	11.6	11.8	5.9	4.2	9.2	12.9	4.9	2.6
Tetracycline	0.9923				0.9948				0.9903			
(RSD %)	6.2	6.4	5.4	3.7	7.3	4.8	5.9	4.5	4.1	9.5	6.5	5.5
Enrofloxacin	0.9969				0.9969				0.9973			
(RSD %)	9.2	7.9	2.3	0.8	11.1	1.4	1.4	3.1	6.9	8.3	3.0	n/a
Difloxacin	0.9958				0.9907				0.9968			
(RSD %)	12.2	4.3	6.0	2.4	10.8	4.6	2.7	5.5	2.6	6.1	5.1	n/a
Spiramycin	0.9920				0.9740				0.9951			
(RSD %)	11.1	11.8	8.4	4.1	10.9	4.0	10.0	9.4	13.3	6.5	5.2	0.2
Albendazole	0.9984				0.9967				0.9928			
(RSD %)	1.6	1.7	1.7	2.7	6.3	3.2	3.6	4.2	2.6	6.2	1.2	2.9
Phenylbutazone	0.9947				0.9922				0.9663			
(RSD %)	2.9	3.3	0.8	2.6	8.1	4.1	4.9	3.1	0.6	0.9	0.7	0.3
Salinomycine Na	0.9993				0.9966				0.9984			
(RSD %)	1.2	0.7	1.2	1.5	1.5	0.8	3.1	0.4	2.8	3.1	1.3	1.2

Stolker, A.A.M. et al; Anal. and Bioanal. Chem. 2010 accepted for publication

Drug identification using ToxID™2.1.0


- Fully automated analysis and reporting
- Drug identification based on
- Molecular weight
- MS2 spectra
- Chromatographic retention time
- Built-in library of about 300 drugs
- Library spectra acquired under real world conditions for robust and accurate ID
- The software uses proven NIST search engine
- Feature to easily create and expand library
- Excellent results review and reporting
- Summary report
- Data review report
- Excel spreadsheet



ToxID Summary Report


Your Company Summary Report

Raw File Name: C:\Documents and Settings\benedicta.duret\My Documents\My Data\Clinical Toxicology\Forensic Toxicology
 Config File Name: C:\Documents and Settings\benedicta.duret\My Documents\My Data\Clinical Toxicology\Forensic Toxicology
 Sample Name: Laboratory:
 Acquisition Start Time: 21/11/2008 11:43:54
 Screening Conditions: Based on Full MS2 scan, m/z window (amu): 0.50, RT window (min): 1.00, MS2 Search Strategy: Tox_Library_LXQ, Use full MS scan to confirm.



Peak Number	Compound Name	Code	SI	RSI	m/z	Expected RT	Actual RT	Concentration ng/ml	Library Name
1	Furosemide	p	999	999	329.00	8.70	9.40	0.03	Tox_Library_LXQ
2	Meprobamate	p	814	814	219.10	7.20	7.48	259.21	Tox_Library_LXQ
3	Venlafaxine	p	831	835	278.20	9.40	10.27	1.79	Tox_Library_LXQ
4	Diazepam	p	842	843	285.10	10.40	10.76	1.46	Tox_Library_LXQ
5	Chlorpromazine-D3	i	814	814	322.20	17.30	17.49	0.12	Tox_Library_LXQ
6	Prazepam-D5	i	911	916	330.20	11.70	12.59	1.69	Tox_Library_LXQ
7	Haloperidol-D4	i	806	826	380.40	14.90	15.30	1.00	Tox_Library_LXQ

Peak Number	Compound Name	Code	SI	RSI	m/z	Expected RT	Actual RT	Concentration ng/ml	Library Name
1	Furosemide	p	999	999	329.00	8.70	9.40	0.03	Tox_Library_LXQ
2	Meprobamate	p	814	814	219.10	7.20	7.48	259.21	Tox_Library_LXQ
3	Venlafaxine	p	831	835	278.20	9.40	10.27	1.79	Tox_Library_LXQ
4	Diazepam	p	842	843	285.10	10.40	10.76	1.46	Tox_Library_LXQ
5	Chlorpromazine-D3	i	814	814	322.20	17.30	17.49	0.12	Tox_Library_LXQ
6	Prazepam-D5	i	911	916	330.20	11.70	12.59	1.69	Tox_Library_LXQ
7	Haloperidol-D4	i	806	826	380.40	14.90	15.30	1.00	Tox_Library_LXQ







ToxID Review Report


Company Name ToxID Long Report

Raw File Name: C:\Documents and Settings\maria.kozak\Desktop\Desktop\Application_Notes\ToxID\2\RAW
 Config File Name: C:\Xcaibur\examples\ToxID\ToxID_config_13min.csv
 Sample Name:
 Laboratory: ChemLab
 Acquisition Start Time: 2/13/2007 1:04:54 AM


Peak Number	Compound Name	Code	SI	ESI	m/z	Expected RT	Real RT	Intensity	Library Name
3	Venlafaxine	p	816	837	278.2	4.90	4.83	12964	Tox_Library




Acquired Spectrum




Library Spectrum



Delta Spectrum



Library Structure



What is Mass Frontier?

- Software for small molecule structural elucidation via mass spectral interpretation
 - Predict fragmentation given a compound structure
 - Annotate spectra with fragment structures
 - Store MSⁿ spectra along with structures, peak annotations, ID numbers, pathway information, etc
 - Match unknown spectra against library entries
 - And **MUCH** more...

Tag Line:


Mass Frontier helps you to go from SPECTRA



Who should get Mass Frontier?


- Anyone who is doing small molecule structural elucidation / confirmation *via* mass spectrometry
- Examples:
 - Metabolite Identification in Drug Metabolism
 - Impurity and Degrading analysis in QC/QA
 - Endogenous Metabolite Identification in Metabolomics
 - Forensic Analyses in Federal and State Agencies
 - Doping Control in Horse Racing
 - Chemistry/Biochemistry/Pharmacy Departments in Universities doing small molecule research
 - Service labs for synthetic chemists

General Unknown Screening using Mass Frontier




Sheldon et al. Determination of Ion Structures in Structurally Related Compounds Using Precursor Ion Fingerprinting. *JASMS*, 2009, 20, 370-376



MSⁿ Spectral Trees—the **ONLY** Route to Unambiguous Structural Elucidation!





MS



MS²



MS³



MS⁴

Accurate mass information is powerful – provides a potential formula

However MSⁿ information still necessary to distinguish between structural isomers

Trees can automatically be generated by Data Dependant LC-MS/MS runs on our instruments

Component Detection from Mass Frontier can automatically deconvolute MSⁿ spectral trees!

This information collectively, uniquely defines the structure of the molecule

How Do You Get a Structure From MS Data?






MS

MS²

MS³

MS³


MS³



+



+



Mass Frontier: Toolbox for Structural Elucidation

HighChem Mass Frontier 6.0

File Edit View Tools Search Library Options Help Microsoft Office

Fragments & Mechanisms: 1

Chromatogram Processor

Easy Structural Editor

The screenshot shows the main interface of HighChem Mass Frontier 6.0. On the left, a list of peaks is displayed with columns for m/z, I_r, and I. The main window features a Total Ion Chromatogram (TIC) at the top and a mass spectrum below it. The mass spectrum shows a base peak at m/z 469.2886. On the right, a chemical structure editor displays a complex organic molecule with a protonated nitrogen atom (H⁺).

Chromatogram Processor

Easy Structural Editor

J-Component Detection Algorithm: Dog_60min_062509_RAW

Mass Merge Power: 95.00 % of Mass Resolution

Deconvolution More

General
Average Peak Width: 20 scans

Noise
 Baseline Correction
 Smoothing

Noise Modification:
 None
 Elimination
 Transition
Min. Peak Height: 0.50 % of max. peak of analysis

Advanced
Analyzed MS Stages:
 Top Stages Lower Stages All Stages
Beginning of Tree Branching: 1 MS Stage
 Retention Time Range
Start: 0.00 min. End: 0.00 min.

Calculate
Reset
Close
Default
Save
Load

Component Detection

Fragments & Mechanisms: 1

Select possible fragments with m/z 147.1168


Fragmentation Pathways

The screenshot shows the 'Fragments & Mechanisms' dialog box. It displays a fragmentation pathway starting from a precursor ion at m/z 177.1274. The pathway involves the addition of a proton (H⁺) to form a protonated intermediate, followed by a loss of a fragment (i) to yield a fragment at m/z 147.1168. A library search (Lib) is also shown for the m/z 147.1168 fragment. On the right, a list of m/z values is displayed, with m/z 147.1168 highlighted.

Fragmentation Pathways

Fragmentation Prediction: Three Knowledge Bases


1. General fragmentation rules



2. Mass Frontier Fragmentation Library™

Total number of	Mass Frontier 6.0
Fragmentation Schemes	30.936
Individual Reactions	129.229
Chemical Structures	151.762
Decoded Mechanisms	120.029

3. User Libraries



Fragmentation Library™ in 6.0 now covers >99% published literature

	Source	Volume	Year
1.	JASMS (Journal of the American Society for Mass Spectrometry)	1-17	1990-2006
2.	IJMSIP (International Journal of Mass Spectrometry and Ion Physics)	1-53	1968-1983
	IJMSIP (International Journal of Mass Spectrometry and Ion Processes)	54-175	1983-1998
	IJMS (International Journal of Mass Spectrometry)	176-255	1998-2006
3.	RCM (Rapid Communications in Mass Spectrometry)	1-20	1987-2006
4.	JMS (Journal of Mass Spectrometry)	30-41	1995-2006
5.	OMS (Organic Mass Spectrometry)	1-29	1968-1994
6.	JMSSJ (Journal of the Mass Spectrometry Society of Japan)	11-27 29-30 37-48 50-53	1964-1979 1981-1982 1989-2000 2002-2005
7.	MSR (Mass Spectrometry Reviews)	1-25	1981-2006
8.	EJMSBMER (European Journal of Mass Spectrometry in Biochemical, Medicine, and Environmental. Research)	1-2	1980-1982
9.	BMS (Biomedical Mass Spectrometry)	1-12	1974-1985
	BEMS (Biomedical and Environmental Mass Spectrometry)	14-19	1987-1990
	BMS (Biological Mass Spectrometry)	20-23	1991-1994
10.	JC (Journal of Chromatography)	181-536	1980-1991
11.	EJMS (European Journal of Mass Spectrometry)	4	1998



Predictive Fragmentation

HighChem Mass Frontier 6.0

File Edit View Tools Search Library Options Help Microsoft Office

Database Manager: 2

Fragments & Mechanisms: 1

Fragments Mechanisms [M+H]⁺

242.0958 248.0614 275.0485 277.0641 284.1394 291.0798 298.1550 300.0801 301.1005 305.0954 313.0879 314.0958 315.1036 m/z 316.1114 317.1192 318.1271 319.1111 328.1114 329.0954 331

Select possible fragments with m/z 316.1114

m/z 346.1220 m/z 346.1220 m/z 316.1114

118.0525 119.0604 120.0808 121.0886 122.0600 132.0318 133.0396 134.0475 135.0553 136.0757 144.0240 145.0475 147.0553 148.0757 149.0709 150.0913 151.0992 166.0321 167.0399 168.0478 172.0553 180.0478 181.0556 184.0427 195.0223 197.0379 198.0583 199.0662 200.0740 208.0301 209.0379 221.0379 225.0692 227.0849 233.0379 242.0958

Database Manager: 2

Reaction Info

m/z 150.0675 m/z 135.0441 m/z 107.0491 m/z 79.0542 m/z 77.0386

Results for

ID	Ac Mol. Mass	Formula	Title
2186	208.1094	C ₁₂ H ₁₆ O ₃ ¹⁺	Seiji Tobita, Susumu Tajima, Yasuko Ishihara
2187	192.1145	C ₁₂ H ₁₆ O ₂ ¹⁺	Seiji Tobita, Susumu Tajima, Yasuko Ishihara
2188	192.1145	C ₁₂ H ₁₆ O ₂ ¹⁺	Seiji Tobita, Susumu Tajima, Yasuko Ishihara
2189	150.0675	C ₉ H ₁₀ O ₂ ¹⁺	Osamu Sekiguchi, Tsutomu Noguchi, Kazuo
2190	150.0675	C ₉ H ₁₀ O ₂ ¹⁺	Osamu Sekiguchi, Tsutomu Noguchi, Kazuo

1. General Rules
2. Literature Library
3. User Libraries

How Do I Annotate Spectral Trees? ...Automatically

The screenshot displays the HighChem Mass Frontier 6.0 interface. The main window shows a mass spectrum plot with several peaks annotated with chemical structures. A context menu is open over the spectrum, with 'Auto Annotation' selected. Below the spectrum, a table lists identified peaks with their molecular masses and formulas. An 'Auto Fragment Annotation of Peaks' dialog box is overlaid on the bottom right, showing configuration options for the auto-annotation process.

ID	Ac Mol. Mass	Formula	Title
1	469.3213	$C_{32}H_{40}N_2O^+H^+$	
2	485.3163	$C_{32}H_{40}N_2O_2^+H^+$	
3	465.3057	$C_{31}H_{38}N_2O^+H^+$	
4	485.3163	$C_{32}H_{40}N_2O_2^+H^+$	
5	501.3112	$C_{32}H_{40}N_2O_3^+H^+$	

Auto Fragment Annotation of Peaks

Options

Threshold: % of highest peak

- Clear Existed Fragments
- Check Precursor
- Auto Annotation Layout
- Apply to All Nodes

Buttons: Restore Defaults, OK, Cancel

Database Manager: Integrated Knowledge Management

- All records of installed libraries are shown in Database Manager
- All records are accessible without querying
- Spectral and Fragmentation libraries are unified in Database Manager
- Searches are universal, independent of data type (structures, m/z values, names, CAS number, biological activity, etc)

One Record: Spectral tree with corresponding fragmentation mechanisms & more!



The left screenshot displays a mass spectrum plot with peaks at m/z 232.09, 246.08, and 261.06. The chemical structure shown is ACEDIASULFONE. The right screenshot displays a fragmentation tree diagram for ACEDIASULFONE, illustrating the loss of H+ and subsequent MS2 and MS3 fragmentation steps.

Match	ID Num.	Mol Mass	Formula	Title	
> 26	1000.0	1	306.0674	C ₁₄ H ₁₄ N ₂ O ₄ S	ACEDIASULFONE
27	94.6	439	307.0838	C ₁₀ H ₁₇ N ₃ O ₆ S	L-Glutathione reduced
28	79.2	63	304.9104	C ₉ H ₅ ClINO	CLIOQUINOL
29	69.7	345	152.0473	C ₈ H ₈ O ₃	
30	69.2	520	152.0473	C ₈ H ₈ O ₃	

HighChem Spectral Tree Libraries—Free with the software!


Library: HighChem ESI **Neg** 2008
Tree Count: 524
Spectra Count: 3805
Fragmentation Schemes: 263

Library: HighChem ESI **Pos** 2008
Tree Count: 1251
Spectra Count: 10180
Fragmentation Schemes: 702



- Common pharmaceutical compounds and human metabolites
- Peaks manually annotated and fragmentation mechanism elucidated

Compound Discoverer




Unknown Analysis

Identification


Statistics

Flexible Workflow




- Use common, predefined workflows or create your own
- Integrate your own nodes
- New in version 2.0: Nodes for Unknown Detection, Identification, Statistics

Predicted Composition



Conclusion



- The Orbitrap Mass Analyzer is a new type of mass analyzers with its own unique combination of analytical parameters
- Orbitraps are still evolving...
 - Higher speed
 - Higher resolving power and mass accuracy
 - Higher sensitivity
 - More routine applications
- Exciting new applications continue to emerge

Summary

- **High resolution** is a key characteristics of MS data enabling
 - Mass accuracy
 - Confident identification
 - Reliable quantitation
- **Data dependent acquisition** offers an elegant simplicity and has proven highly useful for discovery-driven proteomics
- **Mass spectrometry technology** enables comprehensive analysis of proteomics samples
 - Multiple fragmentation techniques
 - MSⁿ capability
- **Quan&Qual** experiments done on a single platform