

Determination of Derivatized Carbamate Insecticides by GC-MS/MS

Thermo Fisher Scientific Inc., Austin, TX, USA

Key Words

- GC-MS/MS
- Carbamates
- Ion Trap
- Method 632
- Pesticides
- Wastewater

Introduction

The determination of carbamate pesticides, insecticides, and herbicides has long been relegated to liquid chromatography. Liquid chromatography was chosen due to the thermally labile nature of the various carbamates. Of the three EPA LC methods referenced, two of these methods (531.1 and 8318A) utilize post column derivatization fluorescence detection while Method 632 monitors the native components. Two of these methods stipulate "...analyte identifications should be confirmed by at least one additional qualitative technique.^{1,2,3} One proposition is to use LC/MS for detection but electrospray ionization generally only produces a single protonated ion and very little fragmentation. While this is good for sensitivity, it does not address the need for confirmation.

This report details a methodology utilizing a GC-MS/MS system that provides both the sensitivity and confirmation for this application in a single injection without the need for extensive sample preparation.

The Thermo Scientific Polaris^Q external source GC-ITMS operated in EI-MS/MS was used to develop a sensitive and robust methodology for the analysis of EPA Method 632 carbamates. This was accomplished by coupling the split injection technique with flash methylation of the analytes in the injection port. By utilizing MS/MS for the N-aryl-carbamates and MS/MS/MS for the O-aryl-carbamates both sensitivity and confirmation can be achieved in standards as well as matrix.

Instrument Conditions

Polaris^Q Ion Trap

Ion source temperature: 250°C
Ionization mode: +EI; 70 eV
AGC: 50
Injection waveforms: ON (default)
Buffer gas flow: ~1.7/mL/min
MS/MS parameters: See Table 1



TRACE GC Ultra

Column: SGE BPX-50 0.25 mm ID x 60 meter, 0.25 micron film thickness
Oven: 70 °C (1 min); 10°C / min 300°C (6 min) : 30 min run time
Split/Splitless injector: Split mode
Injector temperature: 250 °C
Column flow: 1.2 mL/min
Split flow: 12 mL/min
Injection port liner: 5 mm straight; packed with 1 cm glass wool

Autosampler

Injection volume: 2 microliters

Objective

To develop a robust, simple, and sensitive method for the determination of the 19 carbamates by GC-MS/MS. After some literature searches it was determined that many experiments were performed to develop gas chromatographic methodology for the analysis of the carbamates.⁴ Of all the techniques that were attempted for this study, methylation by flash alkylation in the injector was best suited for these 19 compounds. With this technique, all 19 compounds could be chromatographed and confirmed in a single injection by using the conditions described here.

Sample Preparation

100 mL of pond water was extracted with 30 mL of methylene chloride. The methylene chloride was dried with sodium sulfate and the resulting extract was evaporated to dryness. The “blanks” were reconstituted with 100 μ L 50:50 MethElute (Pierce Chemical Company) : methanol. The “spiked” samples were prepared with 100 μ L of 500 pg/ μ L solution of method 632 standards (Accustandard) and evaporated to dryness. The “spiked” samples were then also reconstituted with 50:50 MethElute and methanol. All standards were evaporated to dryness and reconstituted with 50:50 MethElute and methanol.

Results and Discussion

Each of the compounds listed in Table 1 was injected as a single component with the methylating agent to determine the resulting spectrum and retention time. The O-aryl carbamates yielded a methyl aryl ether (Figure 1), whereas the N-aryl carbamates were methylated at the amine (Figure 2). Each individual standard was then injected to optimize collision energies to determine the correct value (the MS/MS experiment was performed as a single injection by alternating scans with varying collision energy). For most of the methyl-aryl-carbamates the MS/MS experiments produced a loss of a methyl group (15 m/z). This loss is not very diagnostic for confirmation. For each of the compounds that exhibited loss of a methyl group, another stage of MS was performed.

After each compound was optimized, the entire mixture was analyzed for interferences. A 6-point (10, 5, 2, 1, 0.5, 0.1 ng/ μ L) calibration curve was injected five times. All extracts were quantified by the external standard technique. See Table 2 for linearity data and Figure 3 and 4 for calibration curves of methylated Carbofuran and methylated Neburon.

A pond water sample was extracted and spiked at 500 pg/ μ L (500 ppt in the water) and injected 20 times to determine analytical reproducibility. These data can be seen in Table 3. The ability to detect and confirm these components can be seen in Figures 5 and 6 for methylated Carbofuran and methylated Neburon respectively. From the spectra obtained, one can easily identify and confirm the presence of these compounds.




Figure 1: Flash methylation reaction of Carbofuran




Figure 2: Flash methylation reaction of Neburon




Figure 3: Linearity of MS/MS/MS of methylated Carbofuran from 0.1-10 ng/ μ L




Figure 4: Linearity of MS/MS of methylated Neburon from 0.1-10 ng/ μ L

Carbofurans in Pond Water Extract




Figure 5: Showing the spectra and mass chromatograms for methylated Carbofuran at the 500 ppt level in both the standard and the pond water extract.




Figure 6: Showing the spectra and mass chromatograms for methylated Neburon at the 500 ppt level in both the standard and the pond water extract. Note the agreement of the chlorine isotope pattern.

| COMPOUND | RETENTION-TIME (MIN) | PRECURSOR ION | MS/MS/MS | QUAN ION |
|--------------|----------------------|---------------|----------|---------------|
| Methomyl | 8.90 | 119 | | 88 |
| Baygon | 11.95 | 166 | 124 | 109 |
| Aminocarb | 13.21 | 165 | 150 | 122 |
| Mexacarbate | 13.68 | 179 | | 148 |
| Carbofuran | 14.19 | 178 | 163 | 131 |
| Propham | 14.47 | 193 | | 151 |
| Fluometuron | 14.57 | 246 | 231 | 174 |
| Oxamyl | 15.98 | 145 | | 61 + 88 + 111 |
| Fenuron | 16.03 | 178 | 163 | 106 |
| Methiocarb | 16.07 | 182 | 167 | 152 |
| Carbaryl | 16.30 | 158 | 143 | 115 |
| Chlorpropham | 16.63 | 227 | | 185 |
| Monuron | 18.45 | 212 | | 146 |
| Swep | 18.80 | 233 | | 176 |
| Barban | 19.49 | 222 | 193 | 166 |
| Linuron | 20.23 | 231 | | 196 |
| Karmex | 20.51 | 246 | | 174 |
| Siduron | 21.66 | 260 | | 163 |
| Neburon | 22.16 | 288 | | 245 |

Table 1: Ions used for the MSⁿ experiments.

| COMPOUND | R ² |
|--------------|----------------|
| Methomyl | 0.9824 |
| Baygon | 0.9962 |
| Aminocarb | 0.9948 |
| Mexacarbate | 0.9948 |
| Carbofuran | 0.9976 |
| Propham | 0.9963 |
| Fluometuron | 0.9957 |
| Oxamyl | 0.9966 |
| Fenuron | 0.9895 |
| Methiocarb | 0.9948 |
| Carbaryl | 0.9964 |
| Chlorpropham | 0.9953 |
| Monuron | 0.9917 |
| Swep | 0.9939 |
| Barban | 0.9922 |
| Linuron | 0.9890 |
| Karmex | 0.9939 |
| Siduron | 0.9917 |
| Neburon | 0.9943 |

Table 2: Showing the R² values obtained for the linear fit for 100 pg/μL to 10 ng/mL.

| COMPOUND | CALC. AMOUNT (PPT) | % RSD (N = 20) |
|--------------|--------------------|----------------|
| Methomyl | 583 | 2.7 |
| Baygon | 470 | 2.9 |
| Aminocarb | 457 | 3.2 |
| Mexacarbate | 475 | 2.1 |
| Carbofuran | 446 | 6.0 |
| Propham | 536 | 2.3 |
| Fluometuron | 400 | 3.7 |
| Oxamyl | 452 | 3.5 |
| Fenuron | 517 | 6.8 |
| Methiocarb | 450 | 3.0 |
| Carbaryl | 412 | 4.9 |
| Chlorpropham | 490 | 2.3 |
| Monuron | 337 | 5.1 |
| Swep | 483 | 3.5 |
| Barban | 223 | 13 |
| Linuron | 554 | 3.6 |
| Karmex | 507 | 3.0 |
| Siduron | 414 | 17 |
| Neburon | 521 | 3.0 |

Table 3: The recovery and precision data for 20 replicate injections of a pond water extract spiked at 500 ppt

Conclusion

By utilizing flash alkylation in the injection port and GC-MS/MS/MS the Polaris_Q is a viable alternative to LC and LC/MS for those laboratories who wish to have a sensitive, robust and simple methodology for the determination of carbamates. This methodology allows detection and confirmation of carbamates with comparable detection limits, reproducibility, and linearity as conventional EPA methods in a single injection.

References

1. Method 632, The Determination of Carbamate and Urea Pesticides in Municipal and Industrial Wastewater US EPA, November 1994.
2. Method 531.1 Rev. 3.1, Measurement of N-Methylcarbamoyloximes and N-Methylcarbamates by Direct Aqueous Injection HPLC with Post Column Derivatization, US EPA, 1995.
3. Method 8318A, N-Methylcarbamates by High Performance Liquid Chromatography (HPLC) Rev. 1, US EPA, 2000.
4. Knapp, D. R., Handbook of Analytical Derivatization Reactions; John Wiley and Sons: New York, 1979, 364

Acknowledgement

Authors: John D Ragsdale III and Meredith Conoley

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

Australia

+61 2 8844 9500

Austria

+43 1 333 50340

Belgium

+32 2 482 30 30

Canada

+1 800 532 4752

China

+86 10 5850 3588

Denmark

+45 70 23 62 60

France

+33 1 60 92 48 00

Germany

+49 6103 408 1014

India

+91 22 6742 9434

Italy

+39 02 950 591

Japan

+81 45 453 9100

Latin America

+1 608 276 5659

Netherlands

+31 76 587 98 88

South Africa

+27 11 570 1840

Spain

+34 91 657 4930

Sweden/Norway/Finland

+46 8 556 468 00

Switzerland

+41 61 48784 00

UK

+44 1442 233555

USA

+1 800 532 4752

www.thermo.com



Thermo Fisher Scientific, Austin, TX USA is ISO Certified.

©2007 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

AN10039_E 08/07C