

USP <467>

Residual solvent determination
in pharmaceutical products.

PRESENTED BY

Vichanee Sornsa

GCMS Product Specialist

- Residual solvents in pharmaceuticals are trace level impurities of Volatile Organic Compounds in **final products** or **excipients**.
- Residual solvents originate from manufacturing **processes** related to the preparation of drug products. They can also form during product **packaging** and **storage**.
- **US Pharmacopeia (USP) method <467>** delineates the procedure to identify and quantify any residual solvent present in pharmaceuticals and excipients as a result of the production process.

The USP 467 method sets the concentration limit for each solvent and describes the analysis of 53 solvents grouped according to their health hazards:

Class 1 : Solvents to be Avoided

- Known human carcinogens
- Strongly suspected human carcinogens
- Environmental hazards

Class 1 Residual Solvent		
Solvent	Concentration Limit (ppm)	Concern
Benzene	2	Carcinogen
Carbon tetrachloride	4	Toxic and environment hazard
1,2-Dichloroethane	5	Toxic
1,1-Dichloroethane	8	Toxic
1,1,1-Trichloroethane	1500	Environment hazard



Class 2 : Solvents to be Limited

- Nongenotoxic animal carcinogens or possible causative agents of other irreversible toxicity, such as neurotoxicity or teratogenicity.
- Solvents suspected of other significant but reversible toxicities.

Class 2 Residual Solvent		
Solvent	PDE (mg/Day)	Concentration Limit (ppm)
Acetonitrile	4.1	410
Chlorobenzene	3.6	360
Chloroform	0.6	60
Cyclohexane	38.8	3880
1,2-Dichloroethene	18.7	1870
Dichloromethane	6.0	600
1,2-Dimethoxyethane	1.0	100
N,N-Dimethylacetamide	10.9	1090
N,N-Dimethylformamide	8.8	880
1,4-Dioxane	3.8	380
2-Ethoxyethanol	1.6	160
Ethyleneglycol	6.2	620
Formamide	2.2	220
Hexane	2.9	290

Class 2 Residual Solvent		
Solvent	PDE (mg/Day)	Concentration Limit (ppm)
Methanol	30.0	3000
2-Methoxyethanol	0.5	50
Methylbutyl ketone	0.5	50
Methylcyclohexane	11.8	1180
N-Methylpyrrolidone	15.3	530
Nitromethane	0.5	50
Pyridine	2.0	200
Sulfolane	1.6	160
Tetrahydrofuran	27.2	720
Tetralin	1.0	100
Toluene	8.9	890
1,1,2-Trichloroethene	0.8	80
Xylenes	21.7	2170


Class 3 Solvents with low toxic potential

- Solvents with low toxic potential to humans; no health-based exposure limit is needed. [NOTE—Class 3 residual solvents may have PDEs of up to 50 mg or more per day.]

Class 3 Residual Solvent	
Solvent	Solvent
Acetic acid	Heptane
Acetone	Isobutyl acetate
Anisole	Isopropyl acetate
1-Butanol	Methyl acetate
2-Butanol	3-Methyl-1-butanol
Butyl acetate	Methylethyl ketone
tert-Butylmethyl ether	Methylisobutyl ketone
Cumene	2-Methyl-1-propanol
Dimethyl sulfoxide	Pentane
Ethanol	1-Pentanol
Ethyl acetate	1-Propanol
Ethyl ether	2-Propanol
Ethyl formate	Propyl acetate

Instrument


- Gas chromatography (GC)
- Headspace (HS) Sampler
- Flame ionization detection (FID)
- Mass Spectrometry (MS) detector (*Optionally*)



Class 1 and Class 2 Residual Solvents Procedure

- Procedure A, Screens for and confirms the solvents' presence using a G43 (volatiles) column.
- Procedure B, Confirms the solvents' identity using a G16 (wax) column.
- Procedure C, Which also uses a G43 column, is required to quantify the amount of residual solvents present.

Diagram relating to the identification limit tests.




GC Parameter	
Column	G43 (Thermo Scientific™ TraceGOLD™ TG 624) 30 m × 0.32 mm × 1.8 μm (P/N 26085-3390) 30 m × 0.53 mm × 3.0 μm (P/N26085-3960)
Split ratio	1:5
Oven Temp.	40 °C hold 20 min. ramp to 240 °C with 10 °C/min. hold 20 min.
Injector Temp.	140
FID Temp.	250

Headspace Operating Parameter Sets			
	1	2	3
Equilibration Temp. (°C)	80	105	80
Equilibration Time (min)	60	45	45
Transfer-line Temp. (°C)	85	110	105
Carrier gas: nitrogen or helium at an appropriate pressure			
Pressureization time (s)	30	30	30
Injection volume (mL)	1	1	1


- **Class 1 Standard Solution**
 - Signal-to-noise ratio of 1,1,1-trichloroethane is not less than 5
- **Class 1 System Suitability Solution**
 - Signal-to- noise ratio of each peak is not less than 3
- **Class 2 Mixture A Standard Solution**
 - Resolution (R) between acetonitrile and methylene chloride (Dichloromethane) is not less than 1.0

Procedure A: System Suitability result

Class 1 Standard Solution



Class 1 System Suitability Solution




- 1 = 1,1-Dichloroethene S/N 126
- 2 = 1,1,1-Trichloroethane S/N 82
- 3 = Carbon Tetrachloride S/N 8
- 4 = Benzene S/N 98
- 5 = 1,2-Dichloroethane S/N 41

- 1 = 1,1-Dichloroethene S/N 769
- 2 = 1,1,1-Trichloroethane S/N 104
- 3 = Carbon Tetrachloride S/N 11
- 4 = Benzene S/N 79
- 5 = 1,2-Dichloroethane S/N 25

Class 1 system suitability solution peak-to-peak signal-to-noise (S/N) ratios for water-soluble (a) and water insoluble (b) products.

Class 2 Mixture A Standard Solution



Chromatographic resolution (R_s) between acetonitrile and dichloromethane for water-soluble (a) and water-insoluble products (b).



GC Parameter	
Column	G16 (Thermo Scientific™ TraceGOLD™ TG WaxMS) 30 m × 0.32 mm × 0.25 μm (P/N 26088-1430) 30 m × 0.53 mm × 0.25 μm (P/N 26088-1440)
Split ratio	1:5
Oven Temp.	50 °C hold 20 min. ramp to 165 °C with 6 °C/min. hold 20 min.
Injector Temp.	140
FID Temp.	250

Headspace Operating Parameter Sets			
	1	2	3
Equilibration Temp. (°C)	80	105	80
Equilibration Time (min)	60	45	45
Transfer-line Temp. (°C)	85	110	105
Carrier gas: nitrogen or helium at an appropriate pressure			
Pressureization time (s)	30	30	30
Injection volume (mL)	1	1	1


- Class 1 Standard Solution
 - Signal-to-noise ratio of benzene is not less than 5
- Class 1 System Suitability Solution
 - Signal-to-noise ratio of each peak in the is not less than 3
- Class 2 System Suitability Solution
 - Resolution (R) between acetonitrile and trichloroethylene is not less than 1.0

Procedure B: System Suitability result

Class 1 Standard Solution



Class 1 System Suitability Solution




- 1 = 1,1-Dichloroethene S/N 114
- 2 = 1,1,1-Trichloroethane S/N 376
- 3 = Benzene S/N 25
- 4 = 1,2-Dichloroethane S/N 40

- 1 = 1,1-Dichloroethene S/N 75
- 2 = 1,1,1-Trichloroethane S/N 214
- 3 = Benzene S/N 112
- 4 = 1,2-Dichloroethane S/N 23


Peak-to-peak signal-to-noise (S/N) ratios for Class 1 system suitability solutions for water-soluble (a) and water-insoluble (b) products.

Class 2 System Suitability Solution



Chromatographic resolution (R_s) between cis 1,2-dichloroethene and acetonitrile for watersoluble (a) and water-insoluble products (b).

Class 2 System Suitability Solution



Chromatographic resolution (R_s) between cis 1,2-dichloroethene and acetonitrile for watersoluble (a) and water-insoluble products (b).

Procedure C Requirements:

Perform the analysis following the instructions for **Procedure A** with the same instrument configuration on the standard solution, the test solution, and the spiked test solution.

Calculate the amount in ppm for each residual solvent found using the formula:


1. $\text{Result} = 5(C/W)[r_u / (r_{st} - r_u)]$ QUANTIFICATION FOR WATER-SOLUBLE

2. $\text{Result} = 10(C/W)[r_u / (r_{st} - r_u)]$ QUANTIFICATION FOR WATER-INSOLUBLE

C is the concentration in μg per mL of the reference standard in the standard stock solution;

W is the weight in grams of the article under test used to prepare the test stock solution;

r_u and r_{st} are the peak responses of each residual solvent obtained from the test solution and the spiked test solution, respectively.



Calibration curves

(at 1 2 5 25 50 and 100 of the concentration limit)

for some of the Class 1 residual solvents

Compound name	%RSD (n=18)		Correlation coefficient (R ²)	Residuals standard deviation (%RSD)
	Water	DMSO		
1,1-Dichloroethene	1.5	0.7	1.000	2.0
1,1,1-Trichloroethane	1.0	0.8	0.999	2.9
Carbon Tetrachloride	4.9	2.9	0.997	6.9
Benzene	0.8	0.9	0.999	3.0
1,2-Dichloroethane	1.6	1.0	0.999	3.7
Methanol	0.7	1.4	1.000	1.4
Acetonitrile	0.8	1.6	1.000	1.7
Dichloromethane	3.1	0.7	0.998	4.2
<i>trans</i> 1,2-Dichloroethene	4.0	1.2	0.999	2.9
<i>cis</i> 1,2-Dichloroethene	3.4	0.8	0.998	5.0
Tetrahydrofuran	0.9	1.4	1.000	2.2
Cyclohexane	3.6	2.8	0.999	3.0
Methylcyclohexane	3.0	2.4	1.000	2.5

Compound name	%RSD (n=18)		Correlation coefficient (R ²)	Residuals standard deviation (%RSD)
	Water	DMSO		
1,4-Dioxane	1.3	1.9	1.000	1.5
Toluene	3.6	0.8	0.997	5.6
Chlorobenzene	3.3	0.7	0.999	2.8
Ethylbenzene	3.4	0.9	0.997	5.3
<i>m</i> -Xylene	3.3	0.9	0.996	6.0
<i>p</i> -Xylene	3.3	0.9	0.996	6.0
<i>o</i> -Xylene	3.1	0.8	0.997	5.6
Hexane	1.2	0.8	0.998	5.8
Nitromethane	2.9	1.5	0.998	4.8
Chloroform	0.9	1.0	0.999	3.0
1,2-Dimethoxyethane	1.4	0.9	0.997	8.4
Trichloroethene	1.9	0.7	0.999	2.9
2-Hexanone	0.6	0.4	1.000	1.3
Tetralin	0.9	0.6	0.999	3.0

Peak area RSDs obtained from n= 18 consecutive injections using water and DMSO as diluents for the concentrated standard solutions, correlation coefficients (R²) and relative standard deviation of residuals RSD) obtained over four calibration levels at 1 2 5 25 50 and 100 of the concentration limits

- All the USP <467> method acceptance criteria were met, and the chromatograms show good peak shapes and resolution.
- These data demonstrate that headspace-gas chromatography is a quantitative technique that easily meets the requirements of the USP <467> method.
- The system configuration of the TRACE 1310 GC, TriPlus 500 Headspace, and Chromeleon CDS is a solid, fully integrated, and reliable platform to perform the analysis according to the USP <467> method with full auditing capability. Operating this system configuration and method in a regulated environment is highly possible obtaining accurate, dependable results.



ติดตามกิจกรรมของทางบริษัทได้ที่



www.scispec.co.th



[/scispec](https://www.facebook.com/scispec)



[@scispec](https://www.line.me/@scispec)



crm@scispec.co.th



ThermoFisher
SCIENTIFIC

MARKES
international



GAS

CTC Analytics



908devices

YOUR SCIENTIFIC SPECIALIST