

Material Emission Testing by TD-GC/MS

PRESENTED BY

Ratimarh Bunlorm

Ratimarh@scispec.co.th

- What is Thermal Desorption (TD)?
- Thermal desorption Sampling.
- Target analytes range.
- Example: Material emission applications.

What is Thermal desorption (TD) ?

Thermal desorption is also

- An alternative to solvent extraction for measuring VOC and SVOC compounds in many different sample matrices – solids, liquids or gases




What is Thermal desorption (TD) ?

- Thermal desorption is arguably the world's most versatile, readily-automated injector mechanism for gas chromatography.
- It is powerful in its own right and is also used as the basis for several other GC front-end technologies:
 - Head space
 - Purge-and-trap,
 - SPME,

The analytical thermal desorption process

Stage 1

- Tubes containing the sample materials or sampled sorbents are heated in a reverse flow of inert (carrier) gas, releasing the trapped compounds and sweeping them into an electrically-cooled, low thermal mass sorbent focusing trap, typically held at -30 to +30°C.




Optional 'inlet' split with re-collection

The analytical thermal desorption process

Stage 1

- Alternatively, whole air/gas samples, from canisters, bags or manifolds, are introduced directly to the focusing trap at controlled flows.



The analytical thermal desorption process

Stage 2

- Once the sample has been focussed onto the trap it is rapidly heated, at rates up to $100^{\circ}\text{C s}^{-1}$, in a reverse flow of carrier gas.
- Retained compounds are released and injected into the GC in a narrow band of vapour - delivering high sensitivity capillary GC performance.




Optional 'outlet' split with re-collection

Thermal desorption Sampling




These methods use two stages of sorbent trapping, and so are ideal for trace-level analysis.



These methods use a single stage of sorbent trapping, and so tend to be used when analyte concentrations are higher.

Image credits: 1 SP Technical Research Institute of Sweden. 2 Owlstone Medical. 3 Equipco.

Extending TD Target Analytes range





Air monitoring



Research applications



Material emission testing



**Flavour and fragrance
profiling**



**Chemical warfare and
forensic**

- Ambient air:
 - ‘Air toxics’ – US: TO-15/17 and CN: HJ 644,
 - Ozone precursors – PAMS,
 - VOC / OVOC,
 - SVOCs.
- Industrial fence line and stack emissions:
 - US Method 325, CN: HJ,
 - CEN TS 13649.
- Indoor and in-vehicle air quality:
 - ISO 16000 and 12219-series methods,
 - Ventilation tests.
- Soil gas and vapor intrusion into buildings.
- Workplace air.
- Odour monitoring – round landfill, waste water treatment and other industrial/urban centres.
- Atmospheric research.



Flavour and fragrance profiling – Food and Drink

- Product aroma profiling:
 - Shelf-life and ripeness studies,
 - Quality assurance,
 - Brand comparison,
 - Identifying taint or adulteration,
 - Confirmation of authenticity,
 - Detecting migration from packaging.
- Allergens in personal care products.
- Cigarettes and e-cigarettes.
- Essential oils in creams/ointments.
- Fragrances in soaps, shampoo and other consumer products.




- Chemical agent destruction.
- Monitoring the safety of agent stockpiles.
- Battlefield protection equipment.
- First responder / mobile labs.
- Counter terrorism / civil defence.
- Decontamination.
- Accelerants in arson debris.
- Forensicating of inks and trace contaminants.
- Detection of traces of proscribed drugs.
- Explosives and propellants.



Material emission testing

- Car trim:
 - ISO 12219, HJ 400 & VDA 278.
- Construction products:
 - ISO 16000 series & EN 16516,
 - Flooring, wall coverings, insulation (etc).
- Decorative products:
 - Adhesives, paints, sealants, etc.
- Spray applied polyurethane foam (SPF) and other PUF products:
 - ASTM D8142.
- Medical devices in contact with patient breathing space:
 - ISO 18562-3 & ISO 18562-4.
- Cleaning products.
- Semi-conductor fabrication & electronics.
- Furniture:
 - BIFMA M7.1.
- Air fresheners & scented candles:
 - EN 16738.






Example:
Material emission applications

Quick screening of small samples with direct desorption

Simplified sampling of solids and liquids

- Sample measured directly into empty TD tube or disposable tube liner and desorbed.
 - Typically at moderate temperatures.
 - Glass tubes with a central restriction hold the sample in the correct position
- Direct thermal desorption of materials can be used:
 - To completion for measuring (S)VOC content
 - Or to get a representative (S)VOC profile
- Provides a 'gas extraction' or 'dynamic headspace' alternative to conventional solvent extraction
- Minimizes sample preparation and injector contamination



High (%) concentrations of analytes accommodated using double splitting.





Stipulations of VDA 278 regarding sample size and preparation

Material type	Sample size
Foam	15 ± 2 mg
Fiber-based materials*	60 ± 20 mg
Films	30 ± 5 mg
Leather	10 ± 2 mg
Paint	50 ± 5 µm film thickness
Adhesives or similar	30 ± 5 mg

* For example, fiberglass or carbon-fiber-reinforced polymer.

Sample preparation


- Cutting the sample to a width of 3 mm.
- Weigh a sample
- Load sample to sample tube.




Sample of rubber flooring material placed into a glass TD tube, ready for direct desorption.

- VOC analysis: this involves desorbing the sample at 90 °C for 30 minutes to quantify volatile compounds up to *n*-C25.


- FOG analysis: FOG compounds (those with volatility range from *n*-C14 to *n*-C32) are then determined by leaving the sample in the desorption tube and raising the temperature to 120 °C for 60 minutes.




Step 2: Trap desorption



Automotive interior VOC and FOG emissions




VOC analysis chromatogram of the leather car trim.




FOG analysis chromatogram of the leather car trim.

Long-term reference tests with small chamber/microchamber




Sample preparation

- Cutting the sample to a width of 2 cm²
- Load sample to chamber cup




Extraction


- Temperature of 65°C and a flow of helium at 50 mL/min
- Equilibration time 20 min.
- A two-bed sorbent tube containing Tenax TA and Carbograph 5TD was then attached and emissions collected for 15 min.



Extracted-ion chromatograms for butadiene and styrene from Sample





Reproducibilities of just 5–8% RSD were obtained when Markes' μ -CTE was used to test six identical samples of plastic material.




Monitoring VOC emissions from respiratory medical devices

Sample preparation and extraction


- Cutting the sample tubing (which in each case was ~2 m long)
- Connect sample to sorbent tube and supply gas.
- Supply gas setup at 100 mL/min for 1 h at room temperature (~21°C), to transfer released VOCs onto the sorbent tube.






2 Acetone	38 Hexanal	53 2-Butoxyethanol	82 Acetophenone	92 n-Tetradecane
13 2-Methylpropan-1-ol	39 n-Octane	54 Cumene	84 n-Undecane	96 Diethyl phthalate
19 Cyclohexane	42 n-Butyl acetate	70 n-Decane	86 1,2,4,5-Tetramethylbenzene	103 Bis(2-ethylhexyl) phthalate
33 Toluene	46 Ethylbenzene	77 o-Cymene	89 n-Dodecane	
34 Pentan-1-ol	50 Cyclohexanone	79 2-Ethylhexan-1-ol	90 n-Tridecane	

Time-course study




Decay profiles for (A) cyclohexanone and (B) toluene over a period of 7 days (samples were not taken on days 5 and 6). Power-law trendlines have been added



Microplastics by TD-GC-MS

Workflow: Sample Preparation

- Simplified using TD-GC-MS



Filtering
0.2µm quartz filter



Washing
with reagents



Drying
30min at 100°C



Quick Preparation
filter into TD-tube



Direct Desorption at
320°C

1 h
Sample
preparation
time


30 mins
analytical cycle

- Hundreds of mL to liters of water are filtered
- The filter is then washed with multiple reagents to remove organic matter
- Dried filter is then placed directly into an empty thermal desorption tube for analysis

Polyethylene terephthalate (PET) in beverages

Identify marker compounds for PET by TD-GC-MS analysis of standard pellets.

- 2,4-di-tert-butylphenol used as quantitation marker
- Tetrahydrofuran is used to confirm presence of PET





ติดตามกิจกรรมของทางบริษัทได้ที่



www.scispec.co.th



[/scispec](https://www.facebook.com/scispec)



[@scispec](https://www.line.me/@scispec)



crm@scispec.co.th



ThermoFisher
SCIENTIFIC

MARKES
international



GAS

CTC Analytics



908devices

YOUR SCIENTIFIC SPECIALIST