

Per- and PolyFluorinated Alkyl Substances (PFAS)

PRESENTED BY

Niramol Jitsommai


Chromatography Product Specialist

niramol@scispec.co.th





Background – What are PFAS compounds?

aka - “Forever Chemicals”



- PFASs are Per- and PolyFluorinated Alkyl Substances. Exclusively anthropogenic.
- Structures contain a hydrophobic perfluoroalkyl backbone and a hydrophilic end group
- Include a diverse range of compounds with a variety of chain lengths and end groups



Persistent in the environment

C-F bonds are very strong and do not break down

Health Concerns* <https://www.atsdr.cdc.gov/pfas/health-effects.html>

- Affect growth, learning, and behavior of infants and children
- Endocrine disruption
- Increase cholesterol levels
- Affect the immune system
- Increase risk of cancer
- Infertility

Background – What are PFAS compounds?




Table 2-1. Discovery and manufacturing history of select PFAS

PFAS ¹	Development Time Period								https://pfas-1.itrcweb.org/
	1930s	1940s	1950s	1960s	1970s	1980s	1990s	2000s	
PTFE	Invented	Non-Stick Coatings			Waterproof Fabrics				
PFOS		Initial Production	Stain & Water Resistant Products	Firefighting foam					U.S. Reduction of PFOS, PFOA, PFNA (and other select PFAS ²)
PFOA		Initial Production	Protective Coatings						
PFNA					Initial Production	Architectural Resins			
Fluoro-telomers					Initial Production	Firefighting Foams		Predominant form of firefighting foam	
Dominant Process ³		Electrochemical Fluorination (ECF)							Fluoro-telomerization (shorter chain ECF)
Pre-Invention of Chemistry /			Initial Chemical Synthesis / Production			Commercial Products Introduced and Used			



- Drinking water Bioaccumulation => PFAS compounds build up in biological systems over time
- Food chain – found in food, feed
- Most people have been exposed => found in blood

PFAS: Source and Exposure



It's Literally Raining 'Forever Chemicals', And The Storm Could Last For Decades

ENVIRONMENT 03 August 2022 By CARLY CASSELLA




<https://www.sciencealert.com/>

The screenshot shows the website for Toxic-Free Future. The header includes the organization's logo and navigation links: About Us, Blog, Pressroom, Contact, and a search icon. Below the header are links for Our Work, Research, Get The Facts, and Take Action, along with a blue Donate button. The main content area features a purple background with a 'PRESS ROOM' label and a headline: '100% of breast milk samples tested positive for toxic "forever chemicals"'. The date 'May 13, 2021' is listed below the headline. To the right of the text is a photograph of a baby being fed. Below the main content, there is a 'Press Contact' section with the name Stephanie Stohler and the email address sstohler@toxicfreefuture.org.

<https://toxicfreefuture.org/>

Challenge to laboratory work



- Need related standards from supplier to quantify

PFAS Master List of PFAS Substances

Search for chemical by systematic name, synonym, CAS number, DTXSID or InChIKey

Identifier substring search

[CompTox Chemicals Dashboard \(epa.gov\)](https://comptox.epa.gov/dashboard/)

List Details

Description: PFASMASTER is a consolidated list of PFAS substances spanning and bounded by the below lists of current interest to researchers and regulators worldwide. For all available lists on the dashboard view [these search results](#).
 Per- and polyfluorinated alkyl substances (PFAS) represent a growing, increasingly diverse inventory of chemicals of interest to the general public, scientific researchers, and regulatory agencies world-wide. Accompanying data-gathering, testing, and environmental monitoring exercises, in turn, have led to the publication and sharing of various lists of PFAS chemicals, some exceeding several thousand substances. A major effort was undertaken by EPA researchers within the National Center for Computational Toxicology to curate and structure-annotate several public lists in DSSTox. The below list of registered PFAS lists, from within and outside EPA, encompass PFAS of potential interest based on environmental occurrence (through literature reports and analytical detection) and manufacturing process data, as well as lists of PFAS chemicals procured for testing within EPA research programs. The consolidated list contains a number of PFAS CAS-name substances, with a subset represented with defined chemical structures. There is no precisely clear definition of what constitutes a PFAS substance given the inclusion of partially fluorinated substances, polymers, and ill-defined reaction products on these various lists. Hence, PFASMASTER serves as a consolidated list of substances spanning and bounded by worldwide. This PFAS Master List will continue to expand as components are added.

Number of Chemicals: 12039

https://comptox.epa.gov/dashboard/chemical_lists/EPAPFASRL is a

https://comptox.epa.gov/dashboard/chemical_lists/EPAPFASINV is a

https://comptox.epa.gov/dashboard/chemical_lists/EPAPFAS75S1 list is a prioritized subset of this larger chemical inventory.

https://comptox.epa.gov/dashboard/chemical_lists/EPAPFASINSOL is a list of chemicals procured, but found to be insoluble in DMSO above 5mM.

https://comptox.epa.gov/dashboard/chemical_lists/PFASOECD is a list of PFAS chemicals in the OECD New Comprehensive Global Database.

https://comptox.epa.gov/dashboard/chemical_lists/PFASKEMI is a list of PFAS chemicals from a KEMI Swedish Chemicals Agency Report (provided by Stellan Fischer).

https://comptox.epa.gov/dashboard/chemical_lists/PFASRIER is a list of PFAS compiled by a community effort in 2015.

https://comptox.epa.gov/dashboard/chemical_lists/EPAPFASCAT is a list of structure-based Markush PFAS categories (capabilities under development).

https://comptox.epa.gov/dashboard/chemical_lists/PFASSTRUCT is a list of all PFAS structures containing a specific defined substructures.

https://comptox.epa.gov/dashboard/chemical_lists/PFASDEV1 is a list of PFAS chemicals without explicit structures - polymers and other UVCB chemicals.

Number of Chemicals: 12039

LC-MS/MS Methods for Quantitation of PFAS

Method Name	PFAS Targets	Sample Preparation	Reporting Limits
EPA Method 537.1	18	SPE (250 mL drinking water); reconstitute 1 mL 96% MeOH	0.53-6.3 ng/L
EPA Method 533	25	SPE (100-250 mL drinking water); reconstitute 1 mL 80% MeOH	1.4-16 ng/L
EPA 8327	24	5 mL water sample + 5 mL MeOH; filter; add 10 uL HOAc	10-50 ng/L
ASTM D7979-16	21	5 mL water sample + 5 mL MeOH; filter	1.5-107 ng/L
ASTM D7968-14	21	2 g Soil + 10 mL 50% MeOH + 0.1% HOAc; filter	2.4-258 ng/kg

PFAS workflow and analytical instrument



Sample

Sample Preparation

Analytical method

Soil sample



X-TRACTION® PFAS (PFE*) or MIX-TRACTION
For Extraction



D-EVA Concentration
For the sensor controlled evaporation of extracts




FREESTYLE™ SPE PFAS
For clean-up up to 50 mL (e.g. Dual SPE)



EluCLEAR® PFAS
SPE columns



D-EVA Concentration
For the sensor controlled evaporation to a few µL



LC-MS/MS Triple quad mass spectrometer

Aqueous sample



FREESTYLE™ XANA PFAS
For clean-up of up to 4 L



FREESTYLE™ XANA PFAS TableTop
For clean-up of up to 250 mL



EluCLEAR® PFAS
SPE columns



D-EVA Concentration
For the sensor controlled evaporation to a few µL




LC-HRAM High resolution accurate mass Q-orbitrap mass spectrometer

PFAS Analysis Kit



Thermo Scientific™ PFC-free kit (P/N 80100-62142)



- Perfluorooctanoic acid
- PFOA
 - Teflon®

Replacing the wetted Teflon™ surfaces with comparable PEEK components and installing a PFOA trapping column.


solvent bottle cap




HPLC vials to do PFAS analysis



PEEK tubing assemblies



PEEK tubing assemblies and the trap column



PEEK Tubing (from Mobile Phase Reservoirs to Vacuum Degasser)

Drinking water

thermoscientific

APPLICATION NOTE 73346

Determination of per- and polyfluorinated alkyl substances (PFAS) in drinking water using automated solid-phase extraction and LC-MS/MS

Authors: Changling Qiu, Xin Zhang, Rahmat Ullah, Wei Chen, and Yan Liu
Thermo Fisher Scientific, Sunnyvale, CA

Keywords: Perfluorinated alkyl substances, AutoTrace 280, U.S. EPA Method 537.1, Accucore RP-MS column, Vanquish Flex Duo UHPLC, TSQ Fortis triple quadrupole mass spectrometer

1. Goal
To demonstrate an efficient and reliable solid-phase extraction method with the Thermo Scientific™ Dionex™ AutoTrace™ 280 instrument for the determination of per- and poly-fluorinated compounds in drinking water per U.S. EPA Method 537.1

2. Introduction
Per- and polyfluorinated alkyl substances (PFAS) are a group of man-made chemicals including perfluorooctanoic acid (PFOA), perfluorooctanesulfonic acid (PFOS), and GenX chemicals that have been manufactured and used in a variety of industries globally.^{1,2} These compounds have a wide range of commercial product applications including industrial polymers, stain repellents, surfactants, waterproofing products, packaging, and aqueous film forming foams used for firefighting. PFAS are highly soluble in water, chemically stable, persistent in the environment, and can accumulate in the human body over time, leading to adverse human health effects.³ PFOA and PFOS are no longer manufactured in the U.S. due to their persistence and potential human health risks.



In November 2019, the United States Environmental Protection Agency (U.S. EPA) published Method 537.1 "Determination of selected per- and polyfluorinated alkyl substances in drinking water by solid phase extraction and LC/MS/MS".⁴ The method uses an offline solid-phase extraction (SPE) with liquid chromatography tandem mass spectrometry (LC-MS/MS) to extract, enrich, and determine 18 PFAS in drinking water. Currently most testing laboratories perform the sample extraction manually using a vacuum manifold, which is labor-intensive, time-consuming, and the flow rate through the cartridge is difficult to control. There is a high demand for automation of the SPE procedure.

ThermoFisher SCIENTIFIC

thermoscientific

APPLICATION NOTE

73883

Determination of per- and polyfluorinated alkyl substances (PFAS) in drinking water

Using automated solid phase extraction and LC-MS/MS for U.S. EPA Method 533

Authors: Xin Zhang, Changling Qiu, Rahmat Ullah, and Yan Liu
Thermo Fisher Scientific, Sunnyvale, CA

Keywords: Perfluorinated and polyfluorinated alkyl substances, PFAS, AutoTrace 280 PFAS, U.S. EPA Method 533, Acclaim 120 C18 column, Vanquish Flex UHPLC, TSQ Fortis triple quadrupole mass spectrometer

Goal
To demonstrate an efficient and reliable solid phase extraction method with a Thermo Scientific™ Dionex™ AutoTrace™ PFAS, an automated solid-phase extraction (SPE) system, for the determination of per- and poly-fluorinated compounds in drinking water per U.S. EPA Method 533

1. Introduction
Drinking water perfluorinated and polyfluorinated alkyl substances (collectively referred to as PFAS) occurrence studies have typically targeted perfluorooctanoic acid (PFOA) and perfluorooctanesulfonate (PFOS), and as a result, these two are the most commonly detected compounds. This discussion focuses primarily on PFOS and PFOA. However, other compounds including PFBA,



PFHxA, PFHpA, PFNA, PFUnDA, and PFHxS have also been detected in drinking water. In December 2019, the United States Environmental Protection Agency (U.S. EPA) took a key step in implementing the PFAS Action Plan by announcing a new validated method for testing PFAS in drinking water. This new validated test method complements other actions the agency is taking under the Action Plan to help communities address PFAS nationwide. It focuses on "short chain" PFAS, those PFAS with carbon chain lengths of four to twelve, and covers PFOS, PFOA, and other common PFAS that have been detected and reported. Currently most testing laboratories are

ThermoFisher SCIENTIFIC

Non-drinking water

thermoscientific

APPLICATION NOTE 65997

Direct analysis of selected per- and polyfluorinated alkyl substances (PFAS) in ground, surface, and waste water by LC-MS/MS

Authors: Cristina C. Jacob, Claudia P.B. Martins, Alan R. Atkins, Richard F. Jack
Thermo Fisher Scientific, San Jose, CA, USA

Keywords: Perfluorinated organic compounds, PFOA, PFOS, PFOA, GenX, PFOS, environmental contaminants, emerging contaminants, EPA 8327, EPA 537, EPA 637.1

Goal
To demonstrate method performance for the PFAS analysis at low levels (ng/L) in a wide variety of non-drinking water matrices by direct analysis and submit data package for EPA 8327 interlaboratory method validation.

Introduction
Per- and polyfluoroalkyl substances (PFAS) are a group of man-made chemicals that includes perfluorooctanoic (PFOA), perfluorooctyl sulfonic acid (PFOS), and hexafluoropropylene oxide dimer acid (HFPO-DA, which is part of GenX process). PFAS compounds have been manufactured since the 1940s. The most well-known PFAS compounds, PFOA and PFOS, have been the most extensively produced and studied for chemical properties and toxicological effects. Both chemicals are very persistent in the environment and accumulate in the human body over time. It is well documented that exposure to PFAS can lead to adverse human health effects¹⁻³ and are found in food packaging material as well as food processing equipment. Plants can accumulate PFAS when grown in PFAS-containing soil and/or water. These compounds are also found in a wide variety of consumer products such as

ThermoFisher SCIENTIFIC

Soil

thermoscientific

CUSTOMER APPLICATION NOTE

73937

Extraction and analysis of poly- and perfluoroalkyl substances (PFAS) from soil

Authors: Matthew S. MacLennan^{1,2}, Daniel Ng¹, David Hope^{1,3}
¹Pacific Rim Laboratories, #103, 19575-55A Avenue, Surrey, B. C. Canada V3S 8P8
²Matthew@pacificrimlabs.com
³Dave@pacificrimlabs.com






Introduction
Recent studies suggest that toxic and highly persistent poly- and perfluorinated alkyl substances (PFAS) are much more prevalent in tissue and soil than in water. The increasing length of perfluoroalkyl chain in PFAS is correlated strongly to lower water solubility/higher adsorption behavior of a particular PFAS molecule in the environment (i.e., migration of PFAS at soil/water/air interfaces) and in remediation/ filtration (i.e., choice of filtration media or sorbents). There are over 6,000 PFAS commercially available, many of which have high environmental persistence and have been found in water and soils globally. This poses a significant challenge to developing analytical methods, especially for the extraction of a variety of PFAS from soil matrices such as soil. Previously, we reported unsatisfactory (0-50%) recovery of long-chain PFAS from soil using vortex/ionization.¹

The soil samples were extracted using the Thermo Scientific™ Dionex™ ASE™ 350 Accelerated Solvent Extractor, which produced 70-130% recovery of all PFAS target compounds. Accelerated solvent extraction has outperformed commonly used, manual "shaking" extraction methods under the same conditions.

After ASE extraction, the solution from the ASE sample collection vials underwent clean-up using solid-phase extraction (SPE) and were analyzed on an LC/MS/MS in a 15-min run. Blanks contained no significant amounts of PFAS. Accelerated solvent extraction is demonstrated to be acceptable for the extraction of short- and long-chain PFAS, with a variety of polarities and head-groups, from soil in the range of 1 ng/g to 400 ng/g.

In the present study, soil was spiked with 24 PFAS (C4-C14 acids, C4-C10 sulfonates, 4:2, 6:2 and 8:2 fluorotelomers, C8 sulfonamide) at 1 ng/g, which were allowed to absorb overnight into the soil samples.

ThermoFisher SCIENTIFIC



LC-MS/MS
Triple quadrupole mass spectrometer

PFAS: Source and Exposure






Figure 2. PFAS occurrence data released by EPA for UCMR 3, using EPA Method 537 and monitoring six PFAS compounds. Data visualization by Moiz Syed. Sources: EPA and Department of Defense. <https://theintercept.com/2015/12/16/toxic-firefighting-foam-has-contaminated-u-s-drinking-water-with-pfcs/>

- the United States Environmental Protection Agency (EPA) developed EPA Method 5374 for the Unregulated Contaminant Monitoring Rule (UCMR 3) program, which collects data for contaminants suspected to be present in drinking water but that do not currently have health-based standards set under the Safe Drinking Water Act (SDWA).⁵
- In 2012, six PFASs were added to the UCMR 3 list to be monitored, including PFOS and PFOA using EPA Method 537.
- In October 2015, occurrence data from the study was released. It is important to note that this is only a small fraction of the hundreds of compounds that can potentially exist in the environment, such as the multiple branched and polyfluorinated PFASs breakdown products that have been known to be in environmental waters.

Screening for untargeted PFAS



Thermo Scientific™ Orbitrap Exploris™ 120 High Resolution, Accurate-Mass (HRAM) Quadrupole Orbitrap MS



Full scan quantitation of target PFASs


+ MS/MS confirmation

+ Screening for other contaminants

Routine quantitative workflows and non-targeted analysis can be performed in a single analysis.

thermoscientific

APPLICATION NOTE 667



A comparison between HRAM Orbitrap technology and MS/MS for the analysis of polyfluoroalkyl substances by EPA Method 537

Authors
 Ali Haghani,¹ Andy Eaton,¹
 Richard F. Jack,² Ed George,²
 Dipankar Ghosh²

¹Eurofins Eaton Analytical, Inc. Monrovia, CA; ²Thermo Fisher Scientific, San Jose, CA

Goal
 To demonstrate a liquid chromatography – high-resolution, accurate-mass (LC-HRAM) methodology using Orbitrap™ technology as a sensitive, accurate, and reliable quantitative alternative to the use of triple quadrupole mass spectrometers while simultaneously determining unknown perfluorinated compounds in the same drinking water extracts.

Introduction
 The unique water-, oil-, grease-, stain- and heat-resistant properties of perfluoroalkyl substances (PFASs) have led to their widespread use in diverse industrial applications and multiple consumer products for over fifty years. Perfluoroalkyl substances are compounds for which all hydrogens on all carbons (except for carbons associated with functional groups) have been replaced by fluorines, e.g., perfluoroalkyl acids (e.g., PFOA, PFOS). Polyfluoroalkyl substances are compounds for which all hydrogens on at least one (but not all) carbons have been replaced by fluorines, e.g., fluorotelomer-based compounds.¹ The carbon-hydrogen linkages allow for biotic and abiotic degradation in the environment. However, the C-F bond

Keywords
 Contaminants of emerging concern, CEC, perfluorinated compound, perfluoroalkyl acid, PFOA, PFOS, perfluorinated alkyl substances, PFAS, perfluoroalkyl acids, PFAA, persistent organic pollutants, POPs, fire-fighting foam, Orbitrap

ThermoFisher SCIENTIFIC

U.S. EPA Method 537 target list

PFAS compound	Critical Level (ng/L)	DL (ng/L)	LCMRL (ng/L)
PFBS	0.15	0.2	<0.5
PFDA	0.15	0.26	<0.5
PFDoA		0.47	0.73
PFHpA	0.09	0.15	<0.5
PFHxA	0.13	0.19	<0.5
PFHxS		1.7	2.4
PFNA	0.11	0.17	<0.5
PFOA		0.22	0.5
PFOS		0.26	0.5
PFTA	0.15	0.2	<0.5
PFTTrDA		0.31	0.55
PFuNA		0.38	1




Figure 8. PFOS branch ratio comparison in a sample vs. a calibration standard. These ratios are represented by overlay of the SRM transitions 499→80 (blue trace) and 499→99 (red trace).


HRAM Orbitrap technology provides equal or better quantitation in full scan as compared to traditional triple quadrupole techniques, with the additional capability to screen for unknown PFASs.

Targeted quantitation & untargeted screening of PFASs

Using a full scan approach, required detection limits or MRLs can be achieved while interrogating for other untargeted PFAS compounds.

The compounds highlighted in blue are additional analytes that are not part of the original U.S. EPA Method 537 list but were found in processed drinking water from the same UCMR3 water extracts.

PFAS compound	Critical Level (ng/L)	DL (ng/L)	LCMRL (ng/L)
PFBA		0.19	0.64
PFODA		0.55	1
PFDS	0.13	0.19	<0.5
PFHxDA		0.12	0.5
PFPA	0.18	0.19	<0.5



The screenshot displays the SciSpec software interface. On the left, a chromatogram shows a peak at 5.451 minutes. The top right shows mass spectra for peaks at 5.451, 5.467, and 5.480 minutes. The bottom left shows a table of compounds with the following data:

#	Checked	Structure	Name	Formula	Molecular Weight	Comments
1	<input type="checkbox"/>		C9 H3 F16 O3 P		493.95643	NO MATCH MS2
2	<input type="checkbox"/>		C9 H2 F16 O2		445.97993	MS2 NOT MATCHES
3	<input checked="" type="checkbox"/>		C6HF11O4S	C6 H F11 O4 S	377.94199	
4	<input type="checkbox"/>		C6HF11O3S	C6 H F11 O3 S	361.94707	C4F7O=280.98346, 9
5	<input type="checkbox"/>		C4HF7O4S	C4 H F7 O4 S	277.94838	

The bottom right shows a table with the following columns: (min), Best SFR [%], Max. # ME, # Adducts, Area, Study File ID, FISH Coverage. The data rows are:

(min)	Best SFR [%]	Max. # ME	# Adducts	Area	Study File ID	FISH Coverage
38	2	1	1	191959	F2	0.00
54	2	1	1	108988	F2	0.96
52	4	1	1	331907	F10	11.54
48	3	1	1	395253	F10	15.38
54	3	1	1	49103	F10	

The Custom Explanation Editor window is open, showing the chemical structure of C6HF11O4S and fields for Molecular weight, Formula, Name, and Comments.



PFAS, a global concern

- The vast uses of PFAS compounds has resulted in a great number of sources for these substances to enter the environment.
- Common modes of transport are runoff, especially with firefighting foams into water, or leachate from landfills, discharge from manufacturing processes or improper disposal.
- Once in the ecosystem, plants, animals and humans are all exposed to PFAS.
- The C-F bonds in the PFAS molecule are very strong, which make metabolism of these compounds very difficult.
- PFAS bioaccumulates over time and causes various health concerns
- The detection and quantification of known PFAS and the discovery of unknown PFAS substances are important



ติดตามกิจกรรมของทางบริษัทได้ที่



www.scispec.co.th



[/scispec](https://www.facebook.com/scispec)



[@scispec](https://line.me/tv/@scispec)



crm@scispec.co.th



ThermoFisher
SCIENTIFIC

MARKES
international



GAS

CTC Analytics



908devices

YOUR SCIENTIFIC SPECIALIST