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

MRL > 0.01 mg/kg

Instrument for analysis pesticide

Comprehensive monitoring of pesticide residues

GC good LC not (GC-G LC-N)	GC good LC bad (GC-G LC-B)	GC bad LC not (GC-B LC-N)	GC bad LC bad (GC-B LC-B)	GC good LC good (GC-G LC-G)	GC bad LC good (LC-G GC-B)	GC not LC good (LC-G GC-N)
Endosulfan, sum	Procymidone	Captan/Folpet	Iprodione	Cyprodinil	Dimethoate, Omethoate	Carbendazim, sum
Brompropylate	Vinclozolin	Chlorothalonil	Dicloran	Fludioxonil	Triadimefon, Triadimenol	Imidacloprid
Chlorfenapyr	Fenitrothion	Dicofol	Oxyfluorfen	Chlorpyrifos	Imazalil	Spinosad
Tetradifon	Tolclofos-methyl			Tolyfluanid	Difenoconazole	Thiacloprid
Chlorthal-dimethyl	Acrinathrin			Myclobutanil	Dimethomorph	Fenbutatin oxide
...
18 pesticides	9 pesticides	3 pesticides	3 pesticides	126 pesticides	46 pesticides	71 pesticides

Performance of the Instrument.

MRL \geq 0.01 mg/kg

Quantitative of Pesticide Residues in Fruits by LC-MS/MS

APPLICATION NOTE 73021

- 160 Compound

Liquid chromatography method

Instrumentation: Vanquish Flex UHPLC

Column: Thermo Scientific™ Hypersil GOLD™ (100 mm × 2.1 mm × 1.9 μm) (P/N 25002-102130)

Sample compartment temp.: 10 °C

Column oven temp.: 25 °C

Injection volume: 5 μL

Mobile phase: A: 2 mM ammonium formate + 0.1% formic acid in water/ acetonitrile (90:10, v/v)
B: 2 mM ammonium formate + 0.1% formic acid in water/ acetonitrile (10:90, v/v)

Total run time: 18.0 min

Gradient program:	Time	Flow Rate	%B	Curve
	0.000	0.400	1	5
	1.500	0.400	1	5
	5.000	0.400	50	5
	8.500	0.400	95	5
	13.500	0.400	95	5
	14.000	0.400	1	5
	18.000	0.400	1	5

Sr. No	Compound	RT (min)	Polarity	Precursor (m/z)	Product Ion (m/z)	Collision Energy (V)
1	3-Hydroxycarbofuran	4.41	Positive	238.1	163	19
		4.41	Positive	238.1	181	15
2	Acephate	4.70	Positive	184.1	49	35
3	Acetamidipid	4.70	Positive	184.1	143	11
		4.85	Positive	223	99	53
		4.85	Positive	223	126	29
4	Aldicarb sulfone**	1.90	Positive	240.1	86.2	28
		1.90	Positive	240.1	148.2	19
5	Ametryn	5.75	Positive	228.1	96	35
		5.75	Positive	228.1	186.1	25
6	Aminocarb	0.74	Positive	209.1	137.1	33
		0.74	Positive	209.1	152	19
7	Amitraz*	8.15	Positive	294.2	91.2	57
		8.15	Positive	294.2	148.3	22
8	Avermectin B1a	9.51	Positive	890.5	305	35
		9.51	Positive	890.5	567.5	17
9	Azoxystrobin	7.38	Positive	404.1	344.1	33
		7.38	Positive	404.1	372.1	19
10	Benalaxyl*	8.12	Positive	326.2	148.1	29
		8.12	Positive	326.2	294.1	15
11	Bendiocarb	6.81	Positive	224.1	109	33
		6.81	Positive	224.1	167.1	15
12	Benzoximate	8.71	Positive	364	105	31
		8.71	Positive	364	199	11
13	Bifenazate*	7.64	Positive	301.1	170.1	27
		7.64	Positive	301.1	198.1	13
14	Bitertanol	7.75	Positive	338.2	70	29
		7.75	Positive	338.2	269.2	13


Figure 1. Total ion chromatograms with overlay extraction ions (160 compounds) in a single window





Figure 3. (A) Extracted ion chromatogram (XIC) for quantifier ion of dimethoate (B), identification based on selection confirmed with ion ratio, and (C) calibration curve

Confident analysis of ultra-trace pesticides residues in baby food using triple quadrupole GC-MS

TRACE 1610 GC parameters	
TRACE 1610 GC parameters	
iConnect PTV parameters	
Injection temperature (°C)	70
Liner	PTV 6 baffle Siltek™ liner (P/N 453T2120)
Injection mode	Splitless
Injection time (min)	0.1
Transfer rate (°C/s)	5.0
Transfer temperature (°C)	300
Transfer time (min)	2.00
Cleaning rate (°C/s)	14.5
Cleaning temperature (°C)	320
Cleaning time (time)	5.00
Cleaning split flow (mL/min)	75
Post cycle temperature	Maintain
Split flow (mL/min)	50
Septum purge flow (mL/min)	5, constant
Carrier gas, flow (mL/min)	He, 1.2

TRACE 1610 GC parameters	
Oven temperature program	
Temperature (°C)	80
Hold time (min)	0.5
Rate (°C/min)	50
Temperature 2 (°C)	190
Rate (°C/min)	6
Temperature 3 (°C)	300
Hold time (min)	5
GC run time (min)	37.70
Column	
Trace GOLD TG-5SilMS	30 m, 0.25 mm, 0.25 µm (P/N 26096-1420)
TSQ 9610 Mass Spectrometer parameters	
Transfer line temperature (°C)	250
Ion source type and temperature (°C)	NeverVent AEI, 320
Ionization type	EI
Emission current (µA)	50
Electron energy (eV)	50
Aquisition mode	timed-SRM
Tuning parameters	AEI SmartTune
Collision gas and pressure (psi)	Argon at 70


- Linearity Range 0.05 to 500 µg/kg

Number	Pesticide	Calibration range (µg/kg)	AvCF %RSD	R ²	IDL (fg OC)	LOD (µg/kg)
1	Allidochlor	0.5–500	7.5	0.9987	112	0.112
2	Dichlobenil	0.05–500	8.0	0.9985	6	0.006
3	Biphenyl	0.05–500	14.7	0.9949	16	0.016
4	Mevinphos	0.05–500	9.6	0.9979	11	0.011
5	3,4-Dichloroaniline	0.05–500	6.7	0.9990	13	0.013
6	Etridiazole (Terrazole)	0.05–500	11.7	0.9970	34	0.034
7	Pebulate	0.5–500	7.4	0.9987	90	0.09
8	N-(2,4-Dimethylphenyl)formamide	0.05–500	10.6	0.9966	15	0.015
9	Methacrifos	0.05–500	5.9	0.9991	23	0.023

Pesticide analysis


Preparation Time : 1.0 hr. -----> Analysis Time : 20 min.
Total : 1 hr. 20 min./sample

Sample Preparation


QuEChERS
Extraction and Cleanup

Analysis

Analysis of 122 pesticide residues by GC-MS/MS

Triple Quadrupole Mass Spectrometer (GC-MS/MS)


Analyzer Mode

- Full scan
- SIM
- SRM/MRM

SRM

Selectivity ↑

Sensitivity (S/N) ↑




- ✓ Unknown screening
- ✓ Low level Quantitative
- ✓ Pesticide
- ✓ PAHs
- ✓ Dioxin (Screening)


GC-MS Component

- 1) Carrier Gas
- 2) Injector
- 3) GC Oven
- 4) Detector (MS)
- 5) Data Readout


Triple Quadrupole Mass Spectrometer (GC-MS/MS)

Q1 selects the *precursor ion*


Argon Collision Gas

Q3 selects the *product ion*



Confident Result from SRM


- Fullscan


- SRM


Result from SRM

283.8 -> 248.8@cid18.00


283.8 -> 213.8@cid30.00

285.8 -> 250.8@cid18.00

Peak	Passed	Target Ratio	Ratio
T1	●	100.00	N/A
T1C1	●	81.31	86.04
T1C2	●	63.54	66.67


Linearity Range: 5 – 500 ng/kg

Analysis of 122 pesticide residues by GC-MS/MS

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

Sample Preparation

QuEChERS (pronounced “catchers”), an acronym for Quick, Easy, Cheap, Effective, Rugged and Safe, covers a variety of sample preparation and clean-up techniques for the analysis of multiple pesticide residues in agricultural matrices.


Materials required:

- Thermo Scientific HyperSep Dispersive SPE Pre-packed Extraction Tubes
 - Thermo Scientific HyperSep Dispersive SPE Clean-up Tubes
 - Acetonitrile, 15 mL/sample
 - Surrogate and internal standard
- Centrifuge and rotor for the tubes used, minimum 3700 RPM

QuEChERS solutions

QuEChERS workflow

The QuEChERS workflow is typically a four-step process: sample homogenization, extraction with salts, followed by SPE clean-up, and analysis by GC-MS or LC-MS. While the basic procedure is more or less standardized, there are an increasing number of chemistry variations to meet different application requirements.

QuEChERS solutions

Whichever QuEChERS method you're using, we've got you covered

Original

(Anastassiades et al, 2003):
developed for non-base-sensitive compounds, using sodium chloride to reduce polar interferences and enhance extraction

AOAC

(AOAC 2007.01):
uses sodium acetate as a buffer instead of sodium chloride and is compatible with base-sensitive compounds

European

(EN 15662):
similar to the AOAC method, but uses sodium chloride to minimize polar interferences, and sodium citrate dihydrate and disodium citrate sesquihydrate instead of sodium acetate

Thermo Scientific QuEChERS Extraction Pouch EN Method

4 g Magnesium Sulfate 1 g Sodium Chloride
1 g Trisodium Citrate Dihydrate
0.5 g Disodium Hydrogencitrate Sesquihydrate

EN method clean-up kits

clean-up		Description
		General fruits and vegetables
		EN clean-up kit, prefilled 2 mL tubes with 150 mg MgSO ₄ , 25 mg PSA
		EN clean-up kit, prefilled 15 mL tubes with 900 mg MgSO ₄ , 150 mg PSA
		Pigmented fruits and vegetables
		EN clean-up kit, prefilled 2 mL tubes with 150 mg MgSO ₄ , 25 mg PSA, 2.5 mg GCB
		EN clean-up kit, prefilled 15 mL tubes with 900 mg MgSO ₄ , 150 mg PSA, 15 mg GCB
		Highly pigmented fruits and vegetables
		EN clean-up kit, prefilled 2 mL tubes with 150 mg MgSO ₄ , 25 mg PSA, 7.5 mg GCB
		EN clean-up kit, prefilled 15 mL tubes with 900 mg MgSO ₄ , 150 mg PSA, 45 mg GCB
		Fruits and vegetables with fats and waxes
		EN clean-up kit, prefilled 2 mL tubes with 150 mg MgSO ₄ , 25 mg PSA, 25 mg C18
		EN clean-up kit, prefilled 15 mL tubes with 900 mg MgSO ₄ , 150 mg PSA, 150 mg C18

What is the purpose of the chemicals?

- **MgSO₄**: removes residual water, and induces phase separation between water content in sample and acetonitrile layer
- **NaCl**: removes residual water, and induces phase separation between water content in sample and acetonitrile layer
- **NaOAc**: buffers the sample to stabilize pH
- **Disodium hydrogencitrate sesquihydrate**: buffers the sample to stabilize pH
- **Trisodium citrate dihydrate**: buffers the sample to stabilize pH
- **PSA**: removes free fatty acids and other acidic co-extractives
- **C18**: removes fats, sterols, and other non-polar interferences from sample
- **GCB**: removes pigment (not recommended for use with planar pesticides)


Analysis of ethylene oxide and 2-chloroethanol Residues in food using GC-MS/MS

Ethylene oxide residue definition in the EU: sum of ethylene oxide and 2-chloro-ethanol expressed as ethylene oxide Reg. (EU) 2015/868

Ethylene oxide

EO and its degradation products

Analysis of ethylene oxide and 2-chloroethanol Residues in food using GC-MS/MS

- Low detection limits required: The maximum residue level (MRL) depends on the commodity, and it ranges from 0.02 to 0.1 mg/kg
- EO is highly volatile:
 - Poorly retained chromatographically
 - Unstable and can be lost in vial via evaporation
- EO has a low molecular weight:
 - Difficult to determine selective SRMs
- Dry matrices:
 - Dirty extracts
 - Interferences can cause over estimation of results

Ethylene oxide

GC MS/MS method

- **TRACE 1610** gas chromatograph,
- **TSQ 9610** GC triple quadrupole and a Thermo Scientific™
- **TriPlus™ RSH** autosampler
- **TG-624SiIMS** (30m x 0.25mm x 1.40um)
- **GuardGOLD™** Capillary Columns (5m x 0.25mm)

Injector

Use this inlet

Temperature Settings

Enable temperature control

Temperature: 90 [0..450 °C]

Inlet Parameters

Operating mode: Split

Split flow control

Split flow: 5.0 [5.0..1250.0 ml/min]

Split ratio: 5

Spillover time: [0.00..300.00 min]

Purge flow control

Purge flow: 5.000 [0.500..50.000 ml/min]

Constant aspirum purge

Stop purge for: [0.00..300.00 min]

Surge pressure: [5.00..1000.00 kPa]

Surge duration: [0.00..300.00 min]

Vacuum compensation

Enable gas saver mode

Gas Saver Flow: 50.0 [5.0..500.0 ml/min]

Gas saver time: 1.50 [0.00..999.99 min]

Enable backflush

Backflush Start Time: [0.00..300.00 min]

Enable Custom Duration

Custom Duration: [0.00..300.00 min]

Set default values

PTV Ramp Settings

Pressure [kPa]	Rate [ml/min]	Temp [°C]	Time [min]	Flow [ml/min]	Back Flush
			0.80	50.0	
	12.0	250	10.00		

Injection: [] [] [] [] [] []

Evap: [] [] [] [] [] []

Transfer: [] [] [] [] [] []

Cleaning: [] [] [] [] [] []

Enable evaporation phase

Enable clean phase

Enable pressure ramp

Transfer delay time: 0.00 [0.00..300.00 min]

Post cycle temperature: Maintain

Display phase program plot

Utilities

Vapour volume calculator

Column flow calculator

Retention time alignment tool

Oven program

No	Retention time [min]	Rate [°C/min]	Target value [°C]	Hold time [min]
1	0.000	Run		
2	2.000	0.00	45.0	2.00
3	4.100	50.00	150.0	0.00
4	16.000	100.00	300.0	10.40
5		New Row		
6	16.000	StopRun		

Transfer line 250°C

Ion source temperature 270 °C


Mass transitions

Compound	CE	Transition
EtO	20	44 -> 14
EtO	5	44 -> 29
2CE	5	80 -> 31
2CE	5	80 -> 43
2CE-d4	5	84 -> 33
2CE-d4	5	86 -> 33


Analysis of ethylene oxide and 2-chloroethanol Residues in food using GC-MS/MS

Sensitivity of analysis


Analysis of ethylene oxide and 2-chloroethanol Residues in food using GC-MS/MS

Linearity for the analysis

Calibration curve for EO between 0.007 mg/kg to 16.5 mg/kg

Calibration curve for 2CE between 0.007 mg/kg to 16.5 mg/kg

Analysis of ethylene oxide and 2-chloroethanol Residues in food using GC-MS/MS

Conclusion

- Chromatography: the chromatographic method provided a very good retention of the analytes and separation from the matrix interferences
- The quantitation at MRL was easily achieved, even with 1 μ L injection volume showing excellent sensitivity
- The XLXR detector facilitates quantitation in a broad range on concentration showing good linearity between 0.007 – 16.5 mg/kg in the samples

Low-level consistent analysis of PBDEs in environmental and food matrices using triple quadrupole GC-MS/MS

Authors

Kjell Hope¹, Giulia Riccardino², Adam Ladak³, and Paul Slicock⁴

¹Pacific Rim Laboratories Inc., Canada

²Thermo Fisher Scientific, Milan, IT

³Thermo Fisher Scientific, Hemel Hempstead, UK

Keywords

Polybrominated diphenyl ethers (PBDEs), food, environment, gas chromatography-mass spectrometry, GC-MS, triple quadrupole, TSQ 9610 mass spectrometer, NeverVent Advanced Ionization ion source (AEI), TRACE 1610 GC, programmable temperature vaporizing injector, PTV, AI/AS 1610

Goal

The aim of this application note is to demonstrate the performance of the Thermo Scientific™ TSQ™ 9610 triple quadrupole mass spectrometer coupled to the Thermo Scientific™ TRACE™ 1610 GC for the determination of polybrominated diphenyl ethers (PBDEs) in environmental and food samples.

Introduction

PBDEs are classes of polybrominated hydrocarbons existing as mixtures of congeners with similar molecular structures but different chemical and physical properties (e.g., congeners with lower numbers of bromine atoms tend to be more volatile and to bioaccumulate more than higher brominated congeners).¹ Historically these compounds were widely used as flame retardants in a variety of products, such as plastics, furniture, upholstery, electrical equipment, electronic devices, textiles, and other household products, because of their capability to release bromine radicals that reduce both the rate of combustion and dispersion of fire when exposed at high temperatures.¹ These compounds enter the environment through emissions from manufacturing processes, volatilization from various products that contain PBDEs, recycling wastes, and leachate from waste disposal sites. They are considered ubiquitous persistent pollutants as they have been detected in the airborne particulate matter, bonded to sediments, surface water, fish, and other marine animals, and therefore represent a risk to human health. As a consequence, the use of certain toxic PBDEs with links to cancer and endocrine disruption (including penta-, tetra-, and deca-PBDE) have been prohibited, and are currently listed in the Stockholm Convention inventory of persistent organic pollutants.²

Authors

Richard Law¹, Alexander Schaechtele², Amit Gujar¹, Jiangtao Xing¹, and Cristian Cojocariu¹

¹Thermo Fisher Scientific, Runcorn, UK

²European Union Reference Laboratory (EURL) for Halogenated POPs in Feed and Food, Freiburg, Germany

³Thermo Fisher Scientific, Austin, Texas

⁴Thermo Fisher Scientific, Beijing, China

Keywords

Triple quadrupole GC-MS/MS, persistent organic pollutants, POPs, polychlorinated dibenzo-*p*-dioxins, polychlorinated dibenzo-*p*-furans, PCDD/Fs, dioxins, polychlorinated biphenyls, PCBs, confirmatory analysis, TSQ 9000, advanced electron ionization, AEI

Goal

To demonstrate the utility of the Thermo Scientific™ TSQ™ 9000 triple quadrupole GC-MS/MS system with Thermo Scientific™ Chromeleon™ Chromatography Data System (CDS) software for the routine and regulatory compliant analysis of polychlorinated dibenzo-*p*-dioxins and polychlorinated dibenzo-*p*-furans (PCDD/Fs), dioxin-like polychlorinated biphenyls (PCBs), and indicator PCBs in food and feed samples.

Introduction

Dioxins and dioxin-like compounds are highly toxic substances classed as persistent organic pollutants (POPs). Due to their high fat-solubility, dioxins accumulate in the fatty tissues of animals. As a result, more than 90% of human exposure to dioxins is through food, especially meat, dairy, fish, etc. Therefore, accurate monitoring of food and feed is essential to control dioxin uptake from the food chain.¹

In 2014 a change in European Commission regulations^{2,3} permitted gas chromatography-triple quadrupole mass spectrometry (GC-MS/MS) to be used as an alternative to gas chromatography-high resolution mass spectrometry (GC-HRMS) for confirmatory analysis and for the control of

Thank You

Q & A

 The world leader in serving science