


# Elemental Analysis in Soil and Fertilizer

PRESENTED BY

Kantima Sitlaothaworn



# What to analyze in Soil and Fertilizer?



## Fertilizer

In the production process, the elemental composition of fertilizers is periodically monitored for their characterization.

- Raw materials
- Finish products
- N, P, K, C/N ratio

## Major nutrients

- Nitrogen (N)
- Phosphorus (P)
- Potassium (K)

## Secondary nutrients

- Calcium (Ca)
- Magnesium (Mg)
- Sulfur (S)

## Micronutrients

- |                 |               |
|-----------------|---------------|
| Boron (B)       | Cobalt (Co)   |
| Copper (Cu)     | Chromium (Cr) |
| Iron (Fe)       | Vanadium (V)  |
| Molybdenum (Mo) | Sodium (Na)   |
| Zinc (Zn)       | Silica (Si)   |
| Manganese (Mn)  |               |

## Toxic


- Arsenic (As)
- Cadmium (Cd)
- Chromium (Cr)
- Mercury (Hg)
- Lead (Pb)
- Copper (Cu)
- Manganese (Mn)
- Nickel (Ni)
- Selenium (Se)



โลหะหนักที่อยู่ในรูปประจุบวก Pb As Cd Hg



ธาตุที่ตรึงบนผิวเม็ดดินที่อยู่ในรูปประจุลบ เช่น P N O C



# Soil and Fertilizer Regulation

## Soil: Regulation in Thailand

ประกาศคณะกรรมการสิ่งแวดล้อมแห่งชาติ เรื่องกำหนดมาตรฐานคุณภาพดิน 2564

| Element | mg/kg                       |                                    |
|---------|-----------------------------|------------------------------------|
|         | ดินประเภทที่ 1<br>อยู่อาศัย | ดินประเภทที่ 2<br>เพื่อการเพาะปลูก |
| As      | 6                           | 25                                 |
| Cd      | 67                          | 762                                |
| Cr      | 17.5                        | 212                                |
| Cu      | 2920                        | 35040                              |
| Pb      | 400                         | 800                                |
| Hg      | 22                          | 263                                |
| Ni      | 436.5                       | 5205                               |
| Mn      | 1710                        | 19640                              |
| Se      | 365                         | 4380                               |



Animal feed



Human food

## Fertilizers: Regulation in Thailand


ประกาศกรมวิชาการเกษตร เรื่อง กำหนดเกณฑ์ปุ๋ยอินทรีย์ / ปุ๋ยหมัก 2557

| Element | mg/kg |
|---------|-------|
| As      | 50    |
| Cd      | 5     |
| Cr      | 300   |
| Cu      | 500   |
| Pb      | 500   |
| Hg      | 2     |


| Element   | % โดยน้ำหนัก |
|-----------|--------------|
| N         | >1           |
| P         | >0.5         |
| K         | >0.5         |
| C/N ratio | 20:1         |
| Na        | <1           |




# Heavy metal



| Toxic elements | Major uses and sources of soil contamination  |
|----------------|---|
| Arsenic (As)   | Pesticides, plant desiccants, animal feed additives, coal and petroleum, mine tailings and detergents |
| Cadmium (Cd)   | Electroplating, pigments for plastics and paints, plastic stabilizers and batteries, fertilizers      |
| Chromium (Cr)  | Stainless steel, Chrome-plated metal, pigments and refractory brick manufacture                       |
| Lead (Pb)      | Combustion of oil, gasoline and coal: Iron and steel production                                       |
| Mercury (Hg)   | Pesticides, catalysts for synthetic polymers, metallurgy, thermometers                                |
| Nickel (Ni)    | Combustion of coal, gasoline, and oil; alloy manufacture, electroplating, batteries                   |



- พิษต่อตับและไต
- ตับอักเสบ
- ทำลายสมอง




- โรคไต
- กระดูกผุ ปวดบริเวณเอวและหลัง
- เป็นสารก่อมะเร็ง โดยเฉพาะมะเร็งปอด มะเร็งต่อมลูกหมากและไต
- โรคอัมพาต อัมพาต



# Techniques routinely used for soil / Fertilizer analysis

|                      | Organic element analyzer | Discrete analyzer | IC | Accelerated Solvent Extraction | GC & GC-MS | X-ray Fluorescence | ICP-OES & ICP-MS |
|----------------------|--------------------------|-------------------|----|--------------------------------|------------|--------------------|------------------|
| Nutrient Analysis    | ●                        | ●                 | —  | —                              | —          | ●                  | ●                |
| Metal Contaminants   | —                        | —                 | —  | —                              | —          | ●                  | ●                |
| Inorganic Anions     | —                        | —                 | ●  | —                              | —          | —                  | —                |
| Organic Contaminants | —                        | —                 | —  | ●                              | ●          | —                  | —                |




trace (ppm to ppb)  
ultra-trace (ppb to ppt)

100 ppm to 100 %



CHNS



Oxygen




1800°C  
Flash Combustion

Organic Elemental Analysis (Combustion)


- The Organic Elemental Analyzer is used for Carbon, Hydrogen, Nitrogen, Sulphur, and Oxygen analysis
- The ICP-OES and ICP-MS are used to provide information of major plant nutrients (N, P, K), secondary plant nutrients (Ca, S, Mg), micronutrients such as B, Mn, Fe, Cu, Zn, Mo and Se, also the toxic elements (As, Cd, Pb and Hg)

# What is Inductively Coupled Plasma (ICP)?

## PERIODIC TABLE CHART




- The Inductively Coupled Plasma (ICP) is an ionization source that fully decomposes a sample into its constituent elements and transforms those elements into ions.
- It is typically composed of argon gas, and energy is "coupled" to it using an induction coil to form the plasma.



# What is ICP-OES?


Sample introduction

Delivers all the sample to the plasma which provide a *non-effect to plasma stability*



Inductively Coupled Plasma

Convert all forms of sample into *free Atom/Ions* and Excite free Atom/Ions turn to *Excited state*



Spectrometer



Wavelength Separator

Distinguish emission light from each elemental Atom/Ions and Detect the intensity of each wavelength

Detector

Readout

Provide all important information when required but in a simple way.



# What is ICP-OES?









**Sample introduction**  
Delivers all the sample to the plasma which provide a *non-effect to plasma stability*

**Inductively Coupled Plasma**  
Convert all forms of sample into *free Atom/Ions* and Excite free Atom/Ions turn to *Excited state*


**Spectrometer**  
Distinguish emission light from each elemental Atom/Ions and Detect the intensity of each wavelength

**Readout**  
Provide all important information when required but in a simple way.






# What is ICP-MS?




Sample preparation


1. Sample introduction system to form a fine aerosol mist from the liquid sample
2. Plasma (ICP) to convert the elements in the sample aerosol to ions
3. Interface to extract the ions into the vacuum system
4. Ion lens to focus the ions and separate them from background signals
5. Collision/reaction cell (CRC) to resolve the analyte ions from interfering ions
6. Mass spectrometer (MS) to filter the analyte ions by mass
7. The electron multiplier detector & Data processing



# Plasma



# How to choose Techniques?



## Analysis Speed

|   |  |
|---|--|
| ICP-OES<br>2-3 min per sample             | ICP-MS<br>2-3 mins per sample            |
| Flame AA<br>15 sec per element per sample | GF AAs<br>2-3 min per element per sample |



PERFORMANCE



SENSITIVITY AND  
DETECTION LIMIT



ANALYSIS  
SPEED



EASE OF USE  
AND  
MAINTENANCE



COST PER  
SAMPLE

### ICP-OES

- Simultaneous
- Multi-elements
- Sub ppb to %
- Highest matrix tolerance
- High sample throughput

### ICP-MS / ICP-MSMS

- Multi-elements
- High dynamic range and lowest limits of detection
- ppb to ppt range
- High sensitivity
- High sample throughput

- ✓ ICP-MS offers high dynamic range and lowest limits of detection
- ✓ ICP-OES has highest matrix tolerance
- ✓ GFAAs offers sensitivity, best for few elements
- ✓ Flame AA offers fast analysis, and economical sol.

# How to choose Techniques?




AAs



ICP-OES




ICP-MS



# Elemental analysis workflow



Standard preparation




Sample preparation



Sample dilution



Sample introduction



Sample analysis



Data



Be aware of contamination sources.



Minimize handling and transfer steps.



Use high-purity reagents.



Use ultrapure water.



Measure weights and volumes with accuracy.



Apply proper skill, technique, and attention to detail.



# Sample Preparation for Soil and Fertilizers



Test Method 3010A, 3051A

## Hazardous Waste Test Methods / SW-846

- EPA Method 200.7: Determination of Metals and Trace Elements in Water and Wastes by [Inductively Coupled Plasma-Atomic Emission Spectrometry](#)
- EPA Method 200.8 Determination of Trace Elements in Waters and Wastes by [Inductively Coupled Plasma-Mass Spectrometry](#)
- EPA Method 6010D (SW-846) [Inductively Coupled Plasma - Atomic Emission Spectrometry](#)
- EPA 6020B ((SW-846) ) [Inductively Coupled Plasma Mass Spectrometry](#)

- **Fertilizer** : Sampling of the fertilizer material should follow accepted guidelines e.g. the International Organization for Standardization (ISO), the Association of Official Agricultural Chemists (AOAC 2017.02), EN method or other applicable protocols

**2.6.35**  
**AOAC Official Method 2006.03**  
**Arsenic, Cadmium, Cobalt, Chromium, Lead, Molybdenum, Nickel, and Selenium in Fertilizers**  
**Microwave Digestion and Inductively Coupled Plasma-Optical Emission Spectrometry**  
 First Action 2006  
 Final Action 2009

INTERNATIONAL  
STANDARD



ISO  
22862

First edition  
2021-11

**Fertilizers and soil conditioners —  
Compound fertilizer — General  
requirements**

*Engrais et amendements — Engrais composé — Exigences générales*

### Method 200.8, Revision 5.4

- Used for drinking water and wastewater compliance monitoring programs under the SDWA and CWA.
- Applies to the analysis of dissolved elements in groundwaters, surface waters, and drinking water.
- Applies to the analysis of total recoverable elements in groundwaters, surface waters, wastewaters, sludges, and soils.



Wet digestion  
Hot plate



Block Acid  
digestion



Microwave  
digestion

### Method 6020B (SW-846)

- In support of RCRA, SW-846 Method 6020B was developed for the ICP-MS analysis of various environmental samples: soils, sediments, sludges, groundwater, surface water, aqueous samples with suspended solids, and industrial wastes requiring a measure of the total leachable elements.
- This is a performance-based method, used for guidance, not for compliance testing.

Sample preparation methods are given in the SW-846 Compendium, Chapter Three, "Inorganic Analytes." The following methods are specified for the digestion of samples prior to ICP-MS analysis: Methods 3005A, 3010A, 3015A, 3020A, 3050B, 3051A, and 3052.


Method 3050B – Acid Digestion of Sediments, Sludges and Soils

Method 3051A – Microwave Assisted Acid Digestion of Sediments, Sludges, Soils and Oils

# Sample preparation techniques

## Open vessel acid digestion





Hot plate



Hot block



- ✓ Simple and inexpensive set-up involving the use of commonplace laboratory apparatus and a hot plate
- ✓ Procedures are standardized and uncomplicated
- ✓ Higher sample sizes (e.g., > 1 gram)





- ✓ Reduced sample handling and transfers
- ✓ Exposure to contamination is reduced
- ✓ Elimination of issues associated with glassware

(adsorption of analytes to container walls or leaching of elemental impurities to the sample solution, are eliminated)

X However, the duration of the digestion process is still extensive, reagent consumption is high, and exposure to contamination from the atmosphere may be possible since it is an open system

## Closed vessel acid digestion

Microwave



- ✓ Speed of digestion
- ✓ Quality digestion
  - ✓ Complete decomposition of the matrix
  - ✓ Cleaner digestate
  - ✓ Better analyte recoveries
- ✓ Reduced exposure to contamination
- ✓ Reduced reagent consumption
- ✓ Retention of analyte

# Quality Control

Quality Control Tests

| Name                             | Description                         |
|----------------------------------|-------------------------------------|
| <b>Blank Tests</b>               |                                     |
| CCB                              | Continuing Calibration Blank        |
| ICB                              | Initial Calibration Blank           |
| MTB                              | Memory Test Blank                   |
| PRB                              | Preparation Blank                   |
| <b>Calibration Tests</b>         |                                     |
| CCV                              | Continuing Calibration Verification |
| ICV                              | Initial Calibration Verification    |
| LCS                              | Laboratory Control Standard         |
| QCS                              | Quality Control Standard            |
| <b>Paired Sample Tests</b>       |                                     |
| DUP                              | Duplicate                           |
| SER                              | Serial Dilution                     |
| <b>Paired Sample Tests (EPA)</b> |                                     |
| DUP EPA                          | Duplicate (EPA)                     |
| SER EPA                          | Serial Dilution (EPA)               |
| <b>Spike Tests</b>               |                                     |
| LFB                              | Laboratory Fortified Blank          |
| MXS                              | Matrix Spike                        |
| PDS                              | Post Digestion Spike                |
| <b>Spike Tests (ARC)</b>         |                                     |
| MXS ARC                          | Matrix Spike (ARC)                  |
| <b>Continuous Tests</b>          |                                     |
| RCV                              | Regression Coefficient Verification |
| RSV                              | Relative Stability Verification     |
| <b>Interference</b>              |                                     |
| SIC                              | Spectral Interference Check         |
| <b>Internal Standard Test</b>    |                                     |
| IST                              | Internal Standard Test              |

Test details for CCB

Number of analyte failures to generate a QC failure: 1

Number of analyte warnings to generate a QC failure: 1

If this QC fails: Ignore and continue from the next sample

If this QC fails again: Ignore and continue from the next sample

If this QC fails a final time: Ignore and continue from the next sample

Test Parameters

| Enabled                             | Analyte           | Warning Limit | Failure Limit |
|-------------------------------------|-------------------|---------------|---------------|
| <input checked="" type="checkbox"/> | Al 396.152 (Aque) | 1             | 2             |
| <input checked="" type="checkbox"/> | Ca 393.366 (Aque) | 1             | 2             |
| <input checked="" type="checkbox"/> | Ca 422.673        |               |               |
| <input checked="" type="checkbox"/> | Co 228.616        |               |               |

Sample List estimated runtime: 1 hours 28 minutes 44 seconds

| No | Date / Time        | Sample Type | Label          | Status | Repeats | SemiQuant | Sample Type |
|----|--------------------|-------------|----------------|--------|---------|-----------|-------------|
| 1  | 2/22/2024 11:23:36 | BLK         |                |        |         |           | BLK         |
| 2  | 2/22/2024 11:25:20 | BLK         |                |        |         |           | BLK         |
| 3  | 2/22/2024 11:26:56 | STD         |                |        |         |           | STD         |
| 4  | 2/22/2024 11:26:56 | STD         | STD 1          |        | 3       |           | STD         |
| 5  | 2/22/2024 11:26:56 | STD         | STD 2          |        | 3       |           | STD         |
| 6  | 2/22/2024 11:30:09 | STD         | STD 3          |        | 3       |           | STD         |
| 7  | 2/22/2024 11:31:45 | STD         | STD 4          |        | 3       |           | STD         |
| 8  | 2/22/2024 11:33:22 | STD         | STD 5          |        | 3       |           | STD         |
| 9  | 2/22/2024 11:34:58 | QC - ICB    | ICB            |        | 3       |           | QC          |
| 10 | 2/22/2024 11:36:34 | QC - QCS    | QCS_L          |        | 3       |           | UNKNOWN     |
| 11 | 2/22/2024 11:38:11 | QC - QCS    | QCS_H          |        | 3       |           | UNKNOWN     |
| 12 | 2/22/2024 11:44:37 | UNKNOWN     | Blank Dup 2    |        | 3       |           | UNKNOWN     |
| 13 | 2/22/2024 11:46:14 | UNKNOWN     | Blank Spk 1    |        | 3       |           | UNKNOWN     |
| 14 | 2/22/2024 11:47:51 | UNKNOWN     | Blank Spk 2    |        | 3       |           | UNKNOWN     |
| 15 | 2/22/2024 11:49:28 | UNKNOWN     | Sample A Dup 1 |        | 3       |           | UNKNOWN     |
| 16 | 2/22/2024 11:51:05 | UNKNOWN     | Sample A Dup 2 |        | 3       |           | UNKNOWN     |
| 17 | 2/22/2024 11:52:43 | UNKNOWN     | Sample A Spk 1 |        | 3       |           | UNKNOWN     |
| 18 | 2/22/2024 11:54:21 | UNKNOWN     | Sample A Spk 2 |        | 3       |           | UNKNOWN     |
| 19 | 2/22/2024 12:02:32 | QC - ICB    | ICB            |        | 3       |           | QC          |
| 20 | 2/22/2024 12:05:45 | QC - QCS    | QCS_L          |        | 3       |           | QC          |

- Calibration
- Requirements for initial and continuing calibration verification are given in Sections 7.24 and 7.25 of Method 6020B Quality control
  - ICB, ICV, CCB, CCV

- Spiked method blank
- Spiked sample (Matrix spiked)
- Spiked Internal standard

Concentrations

| No | Date / Time        | Sample Type | Label | Sc 424.683 (A) | Al 396.152 (Aq) | Cr 206.149 (Aq) |
|----|--------------------|-------------|-------|----------------|-----------------|-----------------|
| 1  | 2/22/2024 11:23:36 | BLK         |       | 100.0%         | -0.002          | -0.001          |
| 2  | 2/22/2024 11:25:20 | BLK         |       | 100.0%         | 0.000           | 0.000           |
| 3  | 2/22/2024 11:26:56 | STD         |       |                |                 |                 |
| 4  | 2/22/2024 11:26:56 | STD         | STD 1 | 99.2%          | 0.110 (0.100)   | 0.117 (0.100)   |
| 5  | 2/22/2024 11:26:56 | STD         | STD 2 | 99.0%          | 0.507 (0.500)   | 0.539 (0.500)   |
| 6  | 2/22/2024 11:30:09 | STD         | STD 3 | 100.1%         | 5.065 (5.000)   | 5.315 (5.000)   |
| 7  | 2/22/2024 11:31:45 | STD         | STD 4 | 99.7%          | 25.493 (25.000) | 26.043 (25.000) |
| 8  | 2/22/2024 11:33:22 | STD         | STD 5 | 96.9%          | 49.747 (50.000) | 49.447 (50.000) |


Calibrations

| No | Date / Time        | Sample Type | Label          | Sc 424.683 (A) | Al 396.152 (Aq) | Cr 206.149 (Aq) |
|----|--------------------|-------------|----------------|----------------|-----------------|-----------------|
| 8  | 2/22/2024 11:34:58 | QC - ICB    | ICB            | 100.5%         | 0.000           | 0.003           |
| 9  | 2/22/2024 11:36:34 | QC - QCS    | QCS_L          | 100.0%         | 1.002 (100.2%)  | 1.058 (105.8%)  |
| 10 | 2/22/2024 11:38:11 | QC - QCS    | QCS_H          | 97.4%          | 29.925 (99.8%)  | 30.559 (101.9%) |
| 11 | 2/22/2024 11:44:37 | UNKNOWN     | Blank Dup 2    | 101.0%         | 0.057           | 0.001           |
| 12 | 2/22/2024 11:46:14 | UNKNOWN     | Blank Spk 1    | 100.4%         | 19.444          | 20.074          |
| 13 | 2/22/2024 11:47:51 | UNKNOWN     | Blank Spk 2    | 94.7%          | 19.239          | 19.911          |
| 14 | 2/22/2024 11:49:28 | UNKNOWN     | Sample A Dup 1 | 100.1%         | 0.082           | 0.006           |
| 15 | 2/22/2024 11:51:05 | UNKNOWN     | Sample A Dup 2 | 99.1%          | 0.032           | 0.000           |
| 16 | 2/22/2024 11:52:43 | UNKNOWN     | Sample A Spk 1 | 88.1%          | 0.993           | 0.982           |
| 17 | 2/22/2024 11:54:21 | UNKNOWN     | Sample A Spk 2 | 104.1%         | 1.036           | 1.090           |
| 18 | 2/22/2024 12:02:32 | QC - ICB    | ICB            | 98.5%          | -0.001          | 0.000           |
| 19 | 2/22/2024 12:04:08 | QC - QCS    | QCS_L          | 97.7%          | 0.978 (97.8%)   | 1.027 (102.7%)  |
| 20 | 2/22/2024 12:05:45 | QC - QCS    | QCS_H          | 96.5%          | 29.570 (98.6%)  | 30.205 (100.7%) |





## Soil, Sediment and Solid Waste EPA Method 6010D (SW-846)



## Fertilizer


AOAC 2017.02: Arsenic, Cadmium, Calcium, Chromium, Cobalt, Copper, Iron, Lead, Magnesium, Manganese, Molybdenum, Nickel, Selenium, and Zinc in Fertilizers: Microwave Acid Digestion and ICP-OES Detection

### Elemental Categories in Fertilizer

| Class                  | Elements                             |
|------------------------|--------------------------------------|
| Primary Nutrients      | N, P, K                              |
| Secondary Nutrients    | Ca, Mg, Fe, Mn, Na, Cu, Zn, Mo, B, S |
| Micronutrients         | Al, Co, V, Se, Ni                    |
| Harmful trace elements | As, Cd, Cr, Pb                       |




# Herbal Big four: As Cd Pb Hg



|  | Type                      | Pb       | As      | Cd        | Hg        |
|--|---------------------------|----------|---------|-----------|-----------|
|  | Herbal medicines          | 10 mg/kg | -       | 0.3 mg/kg | -         |
|  | Raw herbal materials      | 10 mg/kg | 5 mg/kg | 0.3 mg/kg | 0.2 mg/kg |
|  | Raw herbal materials      | 10 mg/kg | 5 mg/kg | 1.0 mg/kg | 0.5 mg/kg |
|  | Finished herbal materials | 10 mg/kg | 5 mg/kg | -         | 0.5 mg/kg |
|  | Finished herbal materials | 20 mg/kg | 5 mg/kg | -         | 0.5 mg/kg |
|  | Raw herbal materials      | 10 mg/kg | 5 mg/kg | 0.3 mg/kg | 0.5 mg/kg |



# Speciation analysis of As in Rice




หน่วยงาน: FAO/WHO  
Guideline: Maximum levels for inorganic As in rice  
ความเข้มข้น: 0.2 mg kg<sup>-1</sup>

หน่วยงาน: EU  
Guideline: Maximum levels of inorganic As in food for infants and young children  
ความเข้มข้น: 0.1 mg kg<sup>-1</sup>

Biswas *et al.* (2020)  
<https://mx.nimt.or.th/?p=13683>

# Herbal Big four: As Cd Pb Hg



|  | Type                      | Pb       | As      | Cd        | Hg        |
|--|---------------------------|----------|---------|-----------|-----------|
|  | Herbal medicines          | 10 mg/kg | -       | 0.3 mg/kg | -         |
|  | Raw herbal materials      | 10 mg/kg | 5 mg/kg | 0.3 mg/kg | 0.2 mg/kg |
|  | Raw herbal materials      | 10 mg/kg | 5 mg/kg | 1.0 mg/kg | 0.5 mg/kg |
|  | Finished herbal materials | 10 mg/kg | 5 mg/kg | -         | 0.5 mg/kg |
|  | Finished herbal materials | 20 mg/kg | 5 mg/kg | -         | 0.5 mg/kg |
|  | Raw herbal materials      | 10 mg/kg | 5 mg/kg | 0.3 mg/kg | 0.5 mg/kg |



# speciation analysis of As in Rice






Table 1. Ion Chromatography operating conditions.

| Parameter        | Value  |
|------------------|--|
| Columns          | Dionex IonPac AG7 (4x50 mm) and AS7 (4x250 mm)   |
| Injection Volume | 100 µL   |
| Eluents          | A: Ultrapure water<br>B: 100 mM (NH <sub>4</sub> ) <sub>2</sub> CO <sub>3</sub> + 3% MeOH, pH 10.3 |
| Gradient         | Time min    Flow mL·min <sup>-1</sup> % A    % B   |
|                  | 0-5            0.7            95    5  |
|                  | 5.1-10        1.5            70    30  |
|                  | 10.1-13      1.5            40    60   |
|                  | 13.1-14      0.7            95    5  |

# Techniques routinely used for soil / Fertilizer analysis

|                      | Organic element analyzer | Discrete analyzer | IC | Accelerated Solvent Extraction | GC & GC-MS | X-ray Fluorescence | ICP-OES & ICP-MS |
|----------------------|--------------------------|-------------------|----|--------------------------------|------------|--------------------|------------------|
| Nutrient Analysis    | ●                        | ●                 | —  | —                              | —          | ●                  | ●                |
| Metal Contaminants   | —                        | —                 | —  | —                              | —          | ●                  | ●                |
| Inorganic Anions     | —                        | —                 | ●  | —                              | —          | —                  | —                |
| Organic Contaminants | —                        | —                 | —  | ●                              | ●          | —                  | —                |

100 ppm to 100 %



CHNS



Oxygen



1800°C  
Flash Combustion


Organic Elemental Analysis (Combustion)

- The Organic Elemental Analyzer is used for Carbon, Hydrogen, Nitrogen, Sulphur

# What is Organic Elemental Analyzer?

Elemental characterization for agricultural purposes gives information useful for determining agronomy management plans. The determination of Nitrogen and carbon, Total Organic Carbon (TOC), and sulfur enables the characterization of

- Soils
- Leaves, Plants, Crops, and other materials
- Fertilizer




Based on the combustion of the sample. Upon combustion, the sample generates uniform compound gases of the elements C, H, N, and S. These combustion products are measured using gas chromatography. Thus, the ratio of the elements in the original sample is determined. C, H, N, and S can all be determined simultaneously, whereas Oxygen by pyrolysis.



|                                      |                |
|--------------------------------------|----------------|
| Quantification of the sample         | Weighting      |
| Quantitative oxidation of the sample | Combustion     |
| Reduction of combustion gases        | Reduction      |
| Separation of the oxidation gases    | Chromatography |
| Generation of signal                 | Detection      |



# CHNS/O Analyzer principle



# Agronomy: Soil / Fertilizer



**Authors**  
Dr. Libera Fazio, Dr. Francesco Lorus and Dr. Guido Cicciocioppo  
Thermo Fisher Scientific, Milan, Italy


**Keywords**  
Argon, CHN Determination, Elemental Analysis, Inductively Coupled Plasma, Soils

**Introduction**  
Carbon, hydrogen and nitrogen are regularly characterized in soils and plants to determine agricultural and environmental practices. As the demand for soils and plants testing has grown in the last years, the classical analytical methods showed to be no longer suitable for their time-consuming sample preparation and for their use of hazardous reagents. For this reason a simple and automated technique is the requirement for modern laboratories dealing with routine analysis.

**Goal**  
This application note shows nitrogen, carbon and hydrogen determination for soils and plants with the FlashSmart EA using argon as carrier gas.

**This note presents data on CHN determination in soils and plants reference materials with different concentration to show the performance of the system using argon as carrier gas and to show the reproducibility of the results obtained.**

ThermoFisher SCIENTIFIC






# Biomass, Biochar and Carbon credits

- Biomass is renewable organic material that comes from plants and animals.
- Biomass is used for heating and electricity generation and as a transportation fuel. Biomass is an important fuel in many countries, especially in developing countries for cooking and heating.
  - Combustion
  - Gasification
  - Fermentation
  - Hydrothermal liquefaction

## Properties of Biomass Sources

- Moisture Content <50%
- Calorific value
- Proportion of Fixed Carbon and Volatile
- Ash/ Residue Content
- Alkali metal: Na K Mg Ca P
- Cellulose/ Lignin Ratio
- Size and Bulk density



Source: U.S. Energy Information Administration (public domain)


- ✓ ทรัพยากรหมุนเวียน : เป็นแหล่งชีวมวลที่ยั่งยืน
- ✓ เป็นกลางทางคาร์บอน : ลดการเพิ่มของก๊าซเรือนกระจกสุทธิสู่ชั้นบรรยากาศ
- ✓ การลดของเสีย
- ✓ ความมั่นคงด้านพลังงาน
- ✓ ประโยชน์ทางเศรษฐกิจ
- ✓ ความคล่องตัว ชีวมวลสามารถแปลงเป็นพลังงานรูปแบบต่างๆ ได้ เช่น ความร้อน ไฟฟ้า และเชื้อเพลิงชีวภาพ (เช่น เอทานอลและไบโอดีเซล)
- ✓ การปล่อยมลพิษต่ำ : มี S น้อยไม่สร้างปัญหาฝนกรด
- ✓ การปรับปรุงดิน เช่น ไร่ข้าว
- ✓ การผลิตพลังงานจากพื้นที่ที่เพาะปลูกยาก



# Biochar and Carbon credits



- Thermal conversion of Biomass
- Biochar is produced during pyrolysis or gasification (350-700 °C)



Source: Patel and Panwar (2023), Castillo *et al.* (2022)

# Agricultural residue and Biochar analysis


## Agro-Residue

- Organic carbon (%)
- C/N ratio
- Total P, K, N

## Biochar

- Organic carbon (35-95%)
- C/H < 0.7 (0.4 for Premium)
- C/O < 0.4
- C/N
- Total P, K, Mg, Ca, Fe, S
- Heavy metal

| รายการ | โลหะปนเปื้อน  | เกณฑ์การปนเปื้อนโลหะสูงสุด (mg/kg น้ำหนักแห้ง) |                 |                 |
|--------|---------------|--|-----------------|-----------------|
|        |               | เกษตรกรกรม                                     | อุตสาหกรรม      |                 |
|        |               |  | ชั้นคุณภาพที่ 1 | ชั้นคุณภาพพิเศษ |
| 1      | สารหนู (As)   | 20   | ไม่กำหนด        | 13              |
| 2      | แคดเมียม (Cd) | 5  | ไม่กำหนด        | 1.5             |
| 3      | โครเมียม (Cr) | 200  | ไม่กำหนด        | 100             |
| 4      | ปรอท (Hg)     | 2  | ไม่กำหนด        | 1               |
| 5      | นิกเกิล (Ni)  | 100  | ไม่กำหนด        | 50              |
| 6      | ตะกั่ว (Pb)   | 300  | ไม่กำหนด        | 120             |
| 7      | ทองแดง (Cu)   | 200  | ไม่กำหนด        | 140             |
| 8      | สังกะสี (Zn)  | 1 000  | ไม่กำหนด        | 420             |





Follow us



[WWW.SCISPEC.CO.TH](https://www.scispec.co.th)



SCISPEC



@SCISPEC



CRM@SCISPEC.CO.TH



LC-MS/MS

GC-MS

ICP-MS

Sample Prep

thermo scientific

Authorized Distributor

MARKES international

SepSolve Analytical



GAS

CTC Analytics

LC Tech 908 devices

VISIT US

[www.scispec.co.th](http://www.scispec.co.th)

scispec

02-454-8533