

Microplastic Analysis by Gas Chromatography-Mass Spectrometry (GC-MS)

PRESENTED BY

Ratimarth Bunlorm

GC GCMS Product Specialist


- Microplastics (MPs)
- Gas Chromatography – Mass Spectrometry (GC-MS)
- Microplastic Analysis by TD-GC/MS
- Microplastic Analysis by PY-GC/MS
- Q&A

Microplastics (MPs)


Plastics are found throughout the environment



Microplastics are any type of plastic polymer less than 5 mm in length*



Due to their small size they can be easily ingested and accumulate in the human body.



- **UN:** Commissioned a report - Legal Limits on Single-Use Plastics and Microplastics: A Global Review of National Laws and Regulations
- **Europe:** ECHA's Committee for Risk Assessment (RAC) currently in a consultation period. Expected to amend REACH Annex XV to include some MPs.
- **WHO:** Released report in 2019, covers risk to human health, identifies knowledge gaps, recommendations for management actions.
- **ISO/CD 24187:** "Principles for the development of standards for investigation procedures of plastics in environmental matrices and related materials"
- **ASTM WK67788:** New Test Method for Identification of Microplastic Particles and fibres in Municipal Wastewater using Pyrolysis-GC/MS.



Gas Chromatography – Mass Spectrometry

GC-MS



Gas Chromatography


- Gas Chromatography (GC) : Chromatography technique which gas is used as mobile phase
- Sample will be injected into the system, Injection port where all components are vaporized and swept into the column
- Sample components will then be separated according to the interaction with stationary phase and eluted to detector.



- The production of ions that are subsequently separated or filtered according to their mass-to-charge (m/z) ratio and detected.
- The resulting mass spectrum is a plot of the (relative) abundance of the produced ions as a function of the m/z ratio.”




Characteristics of a Mass Spectrum

- Graph of Ion Intensity vs. m/z
- Ion Fragments detail structure and molecular weight of compound
- Example : Mass spectrum of Carbon Tetrachloride



Total Ion Chromatogram (TIC)

Extracted Ion Chromatogram (EIC), and Mass Spectrum



Spectrum peak
at RT 12.7 min




Library Search

GC/MS (EI) Library Search



#	Lib.	Match	Prob. (%)	Name
1	M	942	31.8	Caryophyllene
2	M	920	12.6	Bicyclo[5.2.0]nonane, 2-methylene-4,8,8-trimethyl-4-vi...
3	M	907	8.13	Bicyclo[7.2.0]undec-4-ene, 4,11,11-trimethyl-8-methyl...
4	M	889	4.19	Bicyclo[7.2.0]undec-4-ene, 4,11,11-trimethyl-8-methyl...
5	M	888	4.03	1R,3Z,9S-4,11,11-Trimethyl-8-methylenebicyclo[7.2.0...
6	M	876	2.68	Longifolene-(V4)
7	M	875	2.58	(1R,9R,E)-4,11,11-Trimethyl-8-methylenebicyclo[7.2.0...
8	M	873	2.38	(4aS,9aR)-3,5,5,9-Tetramethyl-2,4a,5,6,7,9a-hexahyd...
9	M	863	1.68	β -Longipinene
10	M	863	1.68	Bicyclo[5.3.0]decane, 2-methylene-5-(1-methylvinyl)-8...
11	M	861	1.55	Spiro[5.5]undec-2-ene, 3,7,7-trimethyl-11-methylene-, ...
12	M	861	1.55	Bicyclo[4.3.0]nonane, 7-methylene-2,4,4-trimethyl-2-vi...
13	M	860	1.49	Alloaromadendrene
14	M	858	1.37	Longifolene
15	M	853	1.10	1H-Cycloprop[e]azulene, decahydro-1,1,7-trimethyl-4-...
16	M	851	1.02	Cycloheptane, 4-methylene-1-methyl-2-(2-methyl-1-pr...
17	M	848	0.90	4,11,11-trimethyl-8-methylenebicyclo[7.2.0]undec-3-ene
18	M	846	0.83	(3R,4aS,5R)-4a,5-Dimethyl-3-(prop-1-en-2-yl)-1,2,3,4,...
19	M	846	0.83	1H-Benzocycloheptene, 2,4a,5,6,7,8,9,9a-octahydro-3-...
20	M	844	0.76	Aciphyllene
21	M	844	0.76	(S,1Z,6Z)-8-Isopropyl-1-methyl-5-methylenecyclodec...
22	M	843	0.73	(-)-Tricyclo[6.2.1.0(4,11)]undec-5-ene, 1,5,9,9-tetramet...
23	M	843	0.73	Naphthalene, 1,2,3,5,6,7,8,8a-octahydro-1,8a-dimethy...
24	M	841	0.68	1,4-Methano-1H-indene, octahydro-4-methyl-8-methyl-...
25	M	840	0.65	...



▲ aaaaaa Plant_sample#2917 RT: 1 | Head to Tail MF=942 RMF=942 | Caryophyllene | 942 31.8P


Difference Head to Tail Side by Side Subtraction /

Name: Caryophyllene
Formula: C₁₅H₂₄
MW: 204 **Exact Mass:** 204.1878 **CAS#:** 87-44-5 **NIST#:** 291486 **ID#:** 82944 **DB:** ma
Other DBs: Fine, TSCA, RTECS, HODOC, NIH, EINECS
Contributor: NIST Mass Spectrometry Data Center, 1998.
InChIKey: NPNUFJAVOONJE-WDZFZDKYSA-N **Non-stereo**
Related CAS#: 8007-38-3, 1407-53-0
120 m/z Values and Intensities:

15	2	26	3	27	103	29	94	30	2
37	1	38	7	39	288	40	71	41	769
42	46	43	85	44	3	50	17	51	82
52	63	53	299	54	40	55	356	56	76

(mainlib) Caryophyllene | Plot/Text of Hit / Plot of Hit

Thermal Desorption (TD)-GC/MS



Thermal desorption (TD)

- Thermal desorption is arguably the world's most versatile, readily-automated injector mechanism for gas chromatography.
- An alternative to solvent extraction for measuring VOC and SVOC compounds in many different sample matrices – solids, liquids or gases



Thermal desorption Sampling





Image credits: 1 SP Technical Research Institute of Sweden. 2 Owlstone Medical. 3 Equipco.

The analytical thermal desorption process

Stage 1

- Tubes containing the sample materials or sampled sorbents are heated in a reverse flow of inert (carrier) gas, releasing the trapped compounds and sweeping them into an electrically-cooled, low thermal mass sorbent focusing trap, typically held at -30 to +30°C.



Optional 'inlet' split with re-collection

Stage 2

- Once the sample has been focussed onto the trap it is rapidly heated, at rates up to 100°C/s, in a reverse flow of carrier gas.
- Retained compounds are released and injected into the GC in a narrow band of vapour - delivering high sensitivity capillary GC performance.



Optional 'outlet' split with re-collection


A background image showing a close-up of water splashing into a glass, with many bubbles and droplets. The image is semi-transparent and serves as a backdrop for the text.

Microplastic Analysis by TD-GC/MS

Sample Preparation for analysis MPs (PET) in drinker water



DiffLok caps



Filtering
0.2µm quartz filter



Washing
with reagents



Drying
30min at 100°C



Quick Preparation
filter into TD-tube



Direct Desorption at
320°C

1 h
Sample
preparation
time

30 mins
analytical cycle

- Hundreds of mL to liters of water are filtered
- The filter is then washed with multiple reagents to remove organic matter
- Dried filter is then placed directly into an empty thermal desorption tube for analysis

Workflow: TD-GC/MS analysis

- After trap desorption the compounds separate within the GC column and are detected and identified with the MS.
- This produces a chromatogram which can be used as a chemical fingerprint
- Marker compounds are identified using standards and peak areas can be used to create a calibration curve



TD-GC/MS system



Polyethylene terephthalate (PET) in water

Identify marker compounds for PET by TD-GC-MS analysis of standard pellets.

- 2,4-di-tert-butylphenol (2,4-DTBP) used as quantitation marker
- Tetrahydrofuran is used to confirm presence of PET

Create calibration curve to quantify PET in samples



- Mass of PET vs DTBP peak area




Microplastics in beverages: Case study

Bottled water sample analysis:

- Both the quantification and confirmation markers were identified
- Concentration of PET in a still water sample was quantified as 46 $\mu\text{g/L}$




EIC
2,4-DTBP




EIC
Tetrahydrofuran

Polyethylene terephthalate (PET) in water


- In addition to markers for PET, a number of compounds used in the process of manufacturing plastics have been tentatively found including dimethyl ether, acrolein and cyclopentene and could assist with source profiling.
- Bisphenol A (BPA) is an additive found in plastics to help with hardening. Research suggests this may be an endocrine disruptor so it is a compound of interest in assessing toxicity.



Pyrolysis (Py)-GC/MS




Pyrolysis and pyrolyzates of polymeric materials




Pyrograms of typical polymers


Polyethylene




Polypropylene



Polystyrene




Higher methacrylate copolymer




Pyrogram
fingerprint

Configuration of pyrolysis-GC/MS system




Polymeric material and pyrolysis

Pattern diagram of typical polymeric material




- : Additives
- : Polymer
- : Inorganics



- X: Additives
- Y: Polymer
- Z: Inorganic




Pyrolysis (Single-Shot)



Mixed Information of polymer and volatiles


Sometimes difficult to interpret results !!

Polymer search using pyrolyzates of pyrogram



Pyrogram of unknown polymer (A) and (B) below were identified by library search on NIST library; however, (C) and (E) could not be identified.

Search result for peak (E)



Pyrolyzate candidates generated upon library search.

Polymer candidates which generate 2-Phenethyl-4-phenylpent-4-enitrile by pyrolysis.

F-Search Software


Specifications for F-Search System (Japanese patent 3801355, US patent 6444979)

Product name (P/N)	F-Search "All-In-One" (PY-1110E-181)	Optional libraries (search software F-Search (Ver. 3.6) (PY-1111E-181) required)			
		EGA-MS18B (PY-1112E-181)	PyGC-MS18B (PY-1113E-181)	Pyrolyzate-MS18B (PY-1115E-181)	ADD-MS16B (PY-1114E-161)
Analytical technique	Package contains F-Search (Ver. 3.6) and all four libraries	Evolved gas analysis (EGA-MS)	Pyrolysis-GC/MS (Py-GC/MS) and Thermally assisted hydrolysis and methylation-GC/MS (THM-GC/MS)	Pyrolysis-GC/MS (Py-GC/MS) and Thermally assisted hydrolysis and methylation-GC/MS (THM-GC/MS)	Pyrolysis-GC/MS (Py-GC/MS) and Thermal desorption-GC/MS (TD-GC/MS)
Number of polymers/additives		1,000 polymers	1,000 polymers (THM data in 33 polymers)	268 polymers (THM data in 33 polymers)	494 additives (Py and TD data in 110 additives)
Stored chromatogram		Thermogram	Pyrogram/chromatogram		
Number of mass spectra		c.a. 1,900	c.a. 2,800	c.a. 5,500	c.a. 4,800
Other		Contains all polymers listed in "Pyrolysis - GC/MS Data Book of Synthetic Polymers -Pyrograms, Thermograms and MS of Pyrolyzates-" S. Tsuge , H.Ohtani and C. Watanabe, 2011, Elsevier Inc.			

A background image showing a close-up of water splashing into a glass, with many bubbles and droplets. The image is semi-transparent and serves as a backdrop for the text.

Microplastic Analysis by PY-GC/MS

PY-Sample preparation



Sample Preparation Protocol for analysis MPs in environment sample

MPs in Sediment

Mass reduction by Sieving



Density separation

(Na_2WO_4 , NaCl, NaI, ZnCl_2 aq soln.)



Remove of nature debris from the MPs surface (Chemical; H_2O_2 , enzymatic)

MPs in Water

Sampling with a plankton net (ca. 0.3 mm)



Mass reduction by Sieving




Remove of nature debris from the MPs surface (Chemical; H_2O_2 , enzymatic)

Visual analysis (Shape, color, size, *etc.*) by Microscope


Identification of MPs by FTIR, Raman, PY-GC/MS, *etc.*

Pyrograms


Polyethylene




Polypropylene




Polystyrene



Higher methacrylate copolymer







MPs (Mix of 12) Chromatogram



- How to obtain such as clean chromatogram for separation of trace level MP?
- One click solution for identification and quantification?

GC coupled to high-resolution, accurate-mass spectrometry

Thermo Scientific™ Exactive™ GC Orbitrap™



TraceFinder deconvolution

TraceFinder deconvolution browser showing α -methylstyrene (RT=5.1 min)

Tentative identification based on library (NIST) match (reverse search index, SI 901), fragment rationalization with a confidence score > 97% and mass accuracies of measured fragments

The screenshot displays the TraceFinder deconvolution browser interface. It is divided into several sections:

- Sample List (A):** Shows the sample name 'X3_spikedFish'.
- Peak List (B):** A table of detected peaks. The peak at RT 5.102 is identified as α -Methylstyrene.
- Peak Identifications (C):** A table of identified chemicals. The top entry is α -Methylstyrene with a score of 97.8 and a reverse search index (SI) of 901.
- Extracted Ions (E):** A mass spectrum plot showing relative abundance versus m/z. The base peak is at m/z 119.081056.
- Spectra (D):** A comparison of the acquired spectrum (top) and the NIST library spectrum (bottom) for α -methylstyrene. The x-axis is m/z (49.000000 to 121.000000) and the y-axis is relative abundance (-100 to 100).

Component Name	RT	Reference m/z	BP Area	BP Height
Ethaneperoxy acid, 1-cya...	5.077	105.070160	9702267	490601
α -Methylstyrene	5.102	117.069916	202814689	10573991
1-Propoxy-3,3-diethyltriazen...	5.138	103.041763	10006027	458871
Peak@5.14458	5.145	58.041393	2848492	14194
Cyclooctane, 1,4-dimethyl-...	5.169	69.069977	20773177	18594
3,4-Dimethyl-3-pyrrolin-2-one	5.187	111.068008	541962	3563
2,5-Furandione, dihydro-3-(1...	5.192	107.073051	1136095	5882
Benzene, 1-(1-ethenyl)-1,3...	5.216	117.069916	4686007	273001

Score	Matched Compound	Formula	CAS	SI	HRF Score	RS
97.8	α -Methylstyrene	C9H10	98-83-9	901	99.3559	
97.7	Benzene, cyclopropyl-	C9H10	873-49-4	895	99.3559	
97.7	Benzene, 1-propenyl-	C9H10	637-50-3	895	99.3559	
97.6	α -Methylstyrene	C9H10	98-83-9	892	99.3559	
97.5	Tetracyclo[3.3.1.0(2.8).0(4.6)]n...	C9H10		886	99.3559	
97.5	Benzene, 1-ethenyl-3-methyl-	C9H10	100-80-1	889	99.3559	
97.5	Benzene, 1-propenyl-	C9H10	637-50-3	887	99.3559	

Active	Measured m/z	Area	Height	Fragment ID	Theo m/z
<input checked="" type="checkbox"/>	119.081139	18412105	9609631	C(12)8 C(13)1 H10	119.081056
<input checked="" type="checkbox"/>	118.077736	185529321	96502021	C(12)9 H10	118.077701
<input checked="" type="checkbox"/>	118.073318	19201674	9986735	C(12)8 C(13)1 H9	118.073231
<input checked="" type="checkbox"/>	117.069916	202814689	105739996	C(12)9 H9	117.069876
<input checked="" type="checkbox"/>	116.062309	4878402	2575034	C(12)9 H8	116.062051
<input checked="" type="checkbox"/>	116.057854	9331874	4720205	C(12)8 C(13)1 H7	116.057581
<input checked="" type="checkbox"/>	115.054214	96483220	50379743	C(12)9 H7	115.054226

Samples processed (A), peaks detected (B), identified chemicals (C), acquired versus library spectra (D), and deconvoluted mass spectra for α -methylstyrene (E) are indicated.

Easily identify and quantify unknown microplastic (MPs)

- Quick identification of polymer types for unknown MPs by a patented search algorithm*
- Automatic creation of calibration curve and quick quantification
- Library for 12 commonly used polymers

*Japan Patent #6683335



Microplastic Calibration Standard (MPCS)

MPCS for calibration curve on MPs analysis

- Mixture of 12 polymers

PE / PP / PS / ABS / SBR / PMMA



PC / PVC / PU / PET / N6 / N66

- Micro spatula contained




Automate qualification and quantification work

1. Search marker peaks for each polymer
2. Quantify each polymer based on calibration curve



Pyrograms of the reference polymer mixture and the real marine microplastic sample

*K.Matsui et al. J. Anal. Appl. Pyrol. 149 (2020) 104834



Pyrograms of the reference polymer mixture (a) and of the real microplastic sample (b).

Results of Qual/Quant

Library Match Quality

(%) Plastic Name

Quant Result (µg)


Relative Amount (%)

Peak Area

Polymer	Prob. [%]	Qty [ug]	Ratio [%]	Area	RT [m...]	LOQ [ug]
PE	99.5	11.20	42.5	31420	16.36	7.60
PVC	92.5	9.355	35.5	146285	10.57	2.70
PET	7.8	2.562	9.73	21353	14.10	1.20
SBR	18.8	0.917	3.48	7107	11.50	1.30
PP	89.9	0.691	2.62	4116	6.46	3.90
PS	98.2	0.601	2.28	75144	21.33	0.51
PMMA	99.2	0.375	1.42	39050	4.82	0.69
PU	96.1	0.276	1.05	81556	18.01	0.69
ABS	57.6	0.150	0.57	2697	18.02	0.76
N66	94.1	0.138	0.52	6349	6.23	0.55
N6	61.6	0.058	0.22	3745	11.50	0.23
PC	69.5	0.018	0.07	5027	11.24	0.67
			(100)			

Bar graph of Relative Amount


(%)




Confirmation of searching result

Peak shapes and Mass Spectrum

Marker peak of Sample




Marker peak of MPs Calibration standard



Summation mass spectrum of sample peak



Summation mass spectrum of Reference polymer in library



Each polymer has its own marker peaks and ions

Check peak shape, RT, Mass Spectrum, Integrated peak area

- GC/MS technique can analysis MPs in all particle size is only limited by the size of filter chosen
- TD-GC/MS and PY-GC/MS Simple sample preparation workflows that can be applied to wide ranging sample types. And large sample sizes for enhanced repeatability and sensitivity.
- TD-GC/MS technique can identification and quantification of MPs and provides simultaneous information on targets and non-targets compounds. Additional toxicity information and source profiling.
- Py-GC-MS technique provides detailed information about polymers, additives, and even contaminants.
- Py-GC/MS technique And Software F-Search MPs can identification and quantification of MPs is rapid and automated.



ติดตามกิจกรรมของทางบริษัทได้ที่



www.scispec.co.th



[/scispec](https://www.facebook.com/scispec)



[@scispec](https://www.line.me/@scispec)



crm@scispec.co.th



ThermoFisher
SCIENTIFIC

MARKES
international



GAS

CTC Analytics



908devices

YOUR SCIENTIFIC SPECIALIST