The RadEye SX is a modern and compact multi-purpose meter that fits to most any external scintillation detector

RadEye SX

Survey Meter for external scintillation probes



Key Features

- The RadEye SX can operate with dual phosphor α / β , Nal(Tl) and plastic probes
- Up to 16 different probe configurations are selectable in a submenu
- Easy to replace probes in-the-field, with a simple button-push
- Weighs only 160 g (5.6 oz) with rubber protection, without cable, 110 x 67 x 62 mm (4.3" x 2.6" x 2.4")
- Traditional probes can be used for one and two hand operation
- Versatile operation modes



RadEye BTcom cover for Bluetooth® communication



The RadEye SX is a modern compact multi-purpose survey meter for external scintillator counter tubes. General count rate and surface contamination measurements can be performed as well as dose rate measurements.

As part of the growing RadEye product family of high-end stand-alone meters, the RadEye SX is designed to exceed the most demanding user expectations.

Due to the clear and large display all essential functions and software parameters can be easily accessed. The display and the alarm-LED can be seen while the instrument is worn in the transparent case.

All settings and the data analysis can be done by an optional WindowsTM-

based PC-software and a reader device. The last 1500 mean and maximum values of the count rate or dose rate are recorded internally and can be read out via a serial interface. Additionally the RadEye SX logs the last 250 alarms, error messages and changes of the configuration.

Operation modes:

- •Scaler / Timer with preset countand preset time for sample measurements
- Continuous ratemeter mode for frisker operation
- \cdot Simultaneous α / β and gross / ROI operation (distinct audible indication for each measuring channel)



RadEye SX Specifications

nuclye on opeenioutions			
Measured quantities	Count rate (cps, cpm), surface contamination (Bq, dps, dpm, Bq/cm²), dose rate (R/h, Sv/h, rem/h).		
Background subtraction	In count rate and contamination mode.		
Measuring range	Default: 100,000 cps - extendable by individual probe calibration.		
Probe cables	RG 58, max. 1.5 m (59") – MHV connector.		
High voltage range	100 V1400 V with output impedance 2 MΩ.		
Probe library	16 different detectors with corresponding high voltage, calibration factor, dead time correction, overload threshold, detector area and timeout for detector failure.		
Alarm threshold	Two alarm thresholds for count rate, activity, dose and dose rate each.		
Audible alarm intensity	80 dB at a distance of 30 cm (11.8").		
Working temperature	-20°C + 50°C (-4°F 122°F).		
Relative humidity	10 90 % at 35°C (95°F) not condensing.		
Operating voltage	1,8 4 V, battery low voltage starting from 2.3 V.		
Scaler/Timer	Preset count, preset time.		
EMC	Disturbance emission: EN 61000-6-3, Immunity: EN 61000-6-2.		
Size	110 mm x 67 mm x 62 mm (4.3" x 2.6" x 2.4"), with rubber protection, without cable.		
Weight	Around 160 g (5.6 oz) including AAA cells and protection sleeve.		
Internal memory	The last 1500 measured values are saved and can be read out via PC program. Max and mean value of count rate an dose rate. The time interval is factory preset to 120 s by default. As well scaler measurements and momentaneous readings can be stored manually. Logbook with 250 entries for changes of configuration, occurring alarms and errors		
Averaging filters	Ratemeter filter type: Advanced Digital Filter (ADF), Digital RC-Filter with time constant, 1s180 s, depending on coun rate and count rate change.		
Battery life time	Depends on external detector (required high voltage and dynode chain impedance), > 500 h for instrument alone.		

RadEye SX Accessories

Transparent case with lanyard # 425067044	The RadEye BTcom cover with manual and RadEye.exe update installation file # 425067087	
RadEye SX adapter for DP6 probe # 425069360	Holder with vacuum cup and goose neck # 425522005	
Probe cables MHV to MHV and MHV to PET, various lengths available	Right-angle-plug converter # KT 162245107	
Battery lid # 425067034 for inductive chargers	Desktop holder and inductive charger # 425067083	

© 2010 Thermo Fisher Scientific Inc. All rights reserved. Windows is a registered trademark of Microsoft Corporation in United states and other countries. Bluetooth is a trademark of Bluetooth SIG, Inc., Bellevue, Washington, U.S.A. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidaries. Results may vary under different operating conditions. Thermo Fisher Scientific makes no warranties, expressed or implied, in this product summary. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representatives for details. LIT100816_DB_RadEye SX-e-V1.1

27 Forge Parkway

Franklin, MA 02038 USA

Europe, Africa, Middle East & Countries Not Listed

 Frauenauracher Strasse 96
 +49 (0) 9131 998-226

 D 91056 Erlangen, Germany
 +49 (0) 9131 998-172 fax

 customerservice.eid.erlangen@thermofisher.com
 +49 (0) 9131 998-172 fax

China 7th Floor, Tower West, Yonghe Plaza +86 10 8419 3588 No. 28 Andingem E. Street, Beijing, 100007 China +86 10 8419 3581 fax info.eid.china@thermofisher.com +86 10 8419 3581 fax

Singapore

11 Biopolis Way, Helios, Units #12-07/08 Singapore 138667 info.eid.singapore@thermofisher.com

United Kingdom

Bath Road, Beenham, Reading RG7 5PR United Kingdom customerservice.eid.beenham@thermofisher.com

+44 (0) 118 971 5042 +44 (0) 118 971 2835 fax

+65 6478 9728

+65 6478 9505 fax

USA, Canada, Mexico, Central & South America

info.eid@thermofisher.com India Plot No. C -327, T.T.C. Industrial Area, Pawne

Plot No. C -327, T.T.C. Industrial Area, Pawn Navi Mumbai 400 705, India info.eid.india@thermofisher.com

www.thermoscientific.com/rmp

+1 (508) 553 1700 +1 (800) 274 4212 US toll-free

+1 (508) 520 2815 fax

+91-22-41578800 +91-22-41578801 fax

